

Smart Diagnostics with Universal Array Technology



PaxView[®] COVID-19 real-time RT-PCR Kit

*Innovative and Chimeric Technologies
Leading to a New World of Healthcare*



PaxGenBio

PaxView® COVID-19 real-time RT-PCR Kit

Intended Use

- In vitro diagnostics(IVD) Kits
- Designed for the qualitative detection of COVID-19 viral RNA
- Oropharyngeal swab, nasopharyngeal swab, sputum and bronchoalveolar lavage (BAL)

Introduction

- COVID-19 (Corona Virus disease 2019) is the respiratory infectious disease caused by a novel coronavirus that has not been previously identified in humans.
- The outbreak began in December 2019.
- Reported illness have ranged from mild respiratory symptoms with fever, cough, and shortness of breath to severe illness and death.
- Most estimates of the incubation period for COVID-19 range from 1-14 days.
- The case fatality rate is around 2%.
- The novel coronavirus (nCoV) named as SARS-CoV-2
- beta-coronavirus like MERS-CoV and SARS-CoV.
- All three of these viruses have their origins in bats.
- Sequencing analysis of multiple viral strains from COVID-19 patients suggested that wild animal acted as the natural reservoir for the virus.
- The virus can spread from person to person. The virus is mainly transmitted through contact with respiratory droplets.

Kit Components



Components	4X RT-PCR PREMIX	PRIMER/PROBE MIX	POSITIVE CONTROL
Image			
Volume	240 µl/tube	240 µl/tube	50 µl/tube
Tubes	2 tubes	2 tubes	1 tube



PaxView® COVID-19 real-time RT-PCR Kit

Product description

- Detects the novel coronavirus.
- COVID-19 using one-step real-time RT-PCR in a single tube.
- The target genes for the detection of the virus are Orf1ab and N gene.
- The sequence of primers and probes for the amplification were adapted from sources published by China CDC and WHO.
- The kit includes primers and probe for the amplification of the human RNase P gene that serves as an internal positive control for the real-time RT-PCR.

[References]

1. National institute for viral disease control and prevention (China); http://ivdc.chinacdc.cn/kyjz/202001/t20200121_211337.html
2. Laboratory testing for 2019 novel coronavirus (2019-nCoV) in suspected human cases Interim guidance (WHO)
3. Real-Time RT-PCR Panel for Detection 2019-Novel Coronavirus (CDC)
4. CDC protocol of real time RT-PCR for influenza A(H1N1) (WHO)

- ▶ **Compatible real-time PCR machines:** CFX96(Bio-Rad), ABI 7500/7500Fast(Thermo Fisher Scientific), SLAN-96S(Hongshi), Rotor-Gene Q(QIAGEN)

Performance Characteristics

Analytical Sensitivity (Limit of Detection)

- Serial dilution of in vitro transcription RNA from 10^5 to 0.1 copies/ μ l was tested.
- The limit of detection (LoD) is 10 copies/ μ l (100 copies/reaction).

Analytical Specificity (Cross-reactivity)

- Genes of 12 respiratory viruses and 14 bacterial strains were tested.
- The microorganisms tested were not detected.

Key Features and Benefits

- **User friendly & Simple method**
- **One tube & One step RT-PCR**
- **Use internal control with Rnase P gene from RNA extraction**
→ monitor the entire process including the sample collection, RNA extraction, and RT-PCR


Company	Sample Collection	Nucleic acid extraction	Nucleic acid amplification	Remarks
PaxGenBio	O	O	O	Rnase P gene
A Company	X	O ¹⁾	O	
B Company	X	X	O ²⁾	

1) Internal control is added to the extraction reagent


2) Internal control is added to the PCR reaction

Procedure

PCR protocol



Result analysis



Results Interpretation

1. Parameter setting

Fluorescence	Parameter	Level		
		7500/7500Fast	CFX96-DX	SLAN96S
FAM	Threshold	30,000	350	100
	Baseline	Auto	Auto	Auto
HEX	Threshold	5,000	150	50
	Baseline	Auto	Auto	Auto
Cy5	Threshold	25,000	200	50
	Baseline	Auto	Auto	Auto

2. Ct. value for controls

Control	FAM (Orf1ab)	HEX (N)	Cy5 (IC)
Positive control	≤35	≤35	≤35
Negative control	No Ct.	No Ct.	No Ct.

3. Interpretation

Case	FAM (Orf1ab)	HEX (N)	Cy5 (IC)	Determination
1	+	+	+	Positive
2	+	+	-	Positive*
3	+	-	+/-	Inconclusive**
4	-	+	+/-	
5	-	-	+	Negative
6	-	-	-	Invalid***

* When the amount of PCR products (FAM or HEX) is substantial on positive result data, the signal of internal control (Cy5) may or may not be observed.

** 1) Recommended to retest by increasing the sample concentration.
2) Recommended to proceed with sequencing.

*** Invalid case

- Sample collection fail
- RNA extraction fail



https://youtu.be/E0p_103cQn4

Clinical Evaluation

■ Malaysia - Institute For Medical Research



INSTITUT PENYELIDIKAN PERUBATAN
(Institute For Medical Research)
Jalan Pahang
50588 KUALA LUMPUR
MALAYSIA



Telefon : 03-2616 2666
Faks : 03-2693 9335
http : //www.imr.gov.my

Ruj.kami: IMR/P/15/1501/0046/03 (11)
Tarik : 24 April 2020

Performance Analysis

Test		Tested Kit Assay		Interpretation
		Detected	Not Detected	
IMR In-house Panel	SARS-CoV2 Positive (COVID-19)	15	0	Sensitivity = 100%
	SARS-CoV2 Negative (COVID-19)	0	12	Specificity = 100%

Comments

Use 10ul of extracted RNA while in house assay use 5 ul of extracted RNA

Disclaimer:

The content of this report does not imply that this product is endorsed or recommended by the Institute For Medical Research (IMR). This report may not be reviewed, abstracted, quoted, reproduced or translated, in part or in whole, without the prior written permission of IMR. IMR does not allow the use of this report or reference in any manner in the labeling of this product or kit package.

■ Kazakhstan – National Public Health Center

КАЗАХСТАН РЕСПУБЛИКАСЫ
ДЕНСАУЛЫК САҚТАУ МИНИСТРЛІГІНІҢ
«ҚОҒАМДЫҚ ДЕНСАУЛЫҚ САҚТАУ ҰЛТТЫҚ
ОРТАЛЫҒЫ» ҒАРУАНЫҢ АТЫ АЛҒАН
ҚУРМАНБАЙЛАҒЫ РЕСПУБЛИКАЛЫҚ
МЕДИКЕТТИК КОМПЛЕКСІНІҢ
«САНИТАРИЯЛЫҚ-ЭПИДЕМИОЛОГИЯЛЫҚ
САРАПТАМА ЖӘНЕ МОНИТОРИНГ ҒЫЛЫМИ-
ПРАКТИКАЛЫҚ ОРТАЛЫҒЫ» ФИЛИАЛЫ



ФИЛИАЛ «НАУЧНО-ПРАКТИЧЕСКИЙ ЦЕНТР
САНИТАРИО-ЭПИДЕМИОЛОГИЧЕСКОЙ
ЭКСПЕРТИЗЫ И МОНИТОРИНГА»
РЕСПУБЛИКАНСКОГО ГОСУДАРСТВЕННОГО
ПРЕДПРИЯТИЯ НА ПРАВЕ ХОЗЯЙСТВЕННОГО
ВЕДЕНИЯ «НАЦИОНАЛЬНЫЙ ЦЕНТР
ОБЩЕСТВЕННОГО ЗДРАВООХРАНЕНИЯ»
МИНИСТЕРСТВА ЗДРАВООХРАНЕНИЯ
РЕСПУБЛИКИ КАЗАХСТАН

050008, Алматы қаласы, Фрунзе көшесі 84
факс/тел.(727) 375 61 55
e-mail: npc@npc-ses.kz

050008, г. Алматы, ул. Ауэзова, 84
факс/тел.(727) 375 61 55
e-mail: npc@npc-ses.kz

15-37/4802

14.05.2020

Результаты верификационных исследований тест – набора
«PaxView COVID-19 real-time RT-PCR Kit»

Лаб. №	Результаты 3 серий испытаний тест – набора «PaxView COVID-19 real-time RT-PCR Kit» (Корея)									Тест-система «BGI Biotechnology», Китай		Референсная методика «TIB MOLBIOL», (Германия)	
	1 серия (08.05.2020)			2 серия (11.05.2020)			3 серия (11.05.2020)			Ct	Результат	Ct	Результат
	Ct (Green)	Ct (Yellow)	Результат	Ct (Green)	Ct (Yellow)	Результат	Ct (Green)	Ct (Yellow)	Результат				
6329	23,35	21,06	пол	25,42	23,29	пол	24,39	21,82	пол	17,92	пол	27,73	пол
6204	36,35	32,83	пол	36,57	33,94	пол	36,19	33,17	пол	20,77	пол	32,25	пол
6210	27,07	23,90	пол	27,14	24,53	пол	27,56	24,53	пол	18,82	пол	21,43	пол
2179	27,33	24,83	пол	28,75	26,73	пол	28,54	25,68	пол	23,67	пол	25,28	пол
2182	19,82	17,35	пол	21,33	19,46	пол	20,50	18,26	пол	16,40	пол	19,16	пол
6696	-	-	отр	-	-	отр	-	-	отр	-	отр	-	отр
6895	-	-	отр	-	-	отр	-	-	отр	-	отр	-	отр
7165	-	-	отр	-	-	отр	-	-	отр	-	отр	-	отр
7252	-	-	отр	-	-	отр	-	-	отр	-	отр	-	отр
7244	-	-	отр	-	-	отр	-	-	отр	-	отр	-	отр
KЭ	-	λ	отр	-	-	отр	-	-	отр	-	отр	-	отр
K-	-	-	отр	-	-	отр	-	-	отр	-	отр	-	отр
K+	31,37	30,46	пол	32,05	31,23	пол	30,95	30,19	пол	27,97	пол	29,46	пол

- Comparison kit :TIB MOLBIOL(Germany), BGI Technology(China)
- Equipment used: Rotor-Gene Q(Qiagen)
- Sample used: Positive 5, Negative 5
- Test results: 100% match

Certificate

Certificate

The Certification Body of
TÜV Rheinland LGA Products GmbH

herby certify that the organization
PaxGenBio Co., Ltd.
#803~807, 361, Simin-daero, Dongan-gu,
Anyang-si, Gyeonggi-do 14057
Republic of Korea

has established and applies a quality management system for medical devices
for the following scope:

In accordance with EN ISO 13485:2016 Medical devices
-Quality management systems-Requirements for regulatory
purposes (ISO 13485:2016)
See attachment(s) for the scope

Proof has been furnished that the requirements specified in
EN ISO 13485:2016

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date: 2020-02-05
Certificate Registration No.: SX 60146977 0001
An audit was performed. Report No.: 12031532 001
This Certificate is valid until: 2020-10-29

Certification Body
Masahiro Asami

Date: 2020-02-05

TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
Tel.: +49 201 806 1271 Fax: +49 201 808 3833 e-mail:cert.safety@tuv.com http://www.tuv.com/safety

**Certificate
of EU product notification**

Herewith we confirm that

MT Promed Consulting GmbH
Altenhofstraße 80
66386 St. Ingbert
Germany

has taken over the function of an European Authorized Representative according to the
requirements of Article 10 of the IVDD 98/79/EC for

PaxGenBio Co., Ltd.
#803~807, 361, Simin-daero,
Dongan-gu, Anyang-si
Gyeonggi-do, 14057
Republic of Korea

MT Promed Consulting GmbH has made the product notification at the relevant competent
authority according to Article 10(3).
The in vitro diagnostic medical devices of the manufacturer, covered by the notification, are listed
in Annex I of this certificate.
This certificate does not attest the conformity of the medical devices with the above mentioned
directive. The conformity is stated in the respective product-related Declarations of Conformity
signed under the sole responsibility of the manufacturer.

24. March 2020

Dr. Michael Rinck
- Managing Director -

**Enclosure
Annex I**

인증번호 (No.) : KTC-AQB-5881

**의료기기 제조 및 품질관리 기준 적합인정서
(Certificate of GMP)**

■ 업체명/허가번호 (Company name of Applicant / License No.)
(주)박스젠바이오(제 5295 호)

PaxGenBio Co., Ltd

■ 대표자 (Representative)
박영석 (PARK YOUNG SEK)

■ 업체 소재지 (Company address of Applicant)
경기도 양주시 용인구 시민대로 361 800-806호(사부실 및 연구소), 807호(직접소, 시험실, 창고)

#803-807 361, Simin-daero, Dongan-gu, Anyang-si, Gyeonggi-do, Rep. of Korea

■ 제조소명 (Name of Manufacturer)

제조사 : (주)박스젠바이오(PaxGenBio Co., Ltd)

■ 제조소 소재지 (Address of Manufacturer)

제조사 : 경기도 양주시 용인구 시민대로361, 800-806호(사부실및연구소), 807호(직접소, 시험실, 창고)
#803-807 361, Simin-daero, Dongan-gu, Anyang-si, Gyeonggi-do, Rep. of Korea

■ 품목군 (Category)
의료진단 목적기름 시약류(Reagent for In-Vitro Diagnostic Device)

의료기기 제조 및 품질관리기준에 적합함을 인정합니다.
(We hereby certify that the above manufacturer complies with Korea
Good Manufacturing Practices of Medical Devices for the product group listed above)

발행일자(Date of Issue) : 2019. 07. 03
유효기간(Date of Expiration) : 2020. 12. 28

경인지방식품의약품안전청
GYEONGGI REGIONAL FOOD AND DRUG ADMINISTRATION

한국기계전기전자시험연구원
Korea Testing Certification

Document Number : SFCW-HFM-CGSP-MSK

Seoul Health Technology Administration Office,
147 Gokseonggyeong-ro, Gyeongseong, Neungdeok-gu,
Chongju-si, Chungcheongnam-do, Korea, 28186
Tel: 182-48-719-2356, Fax: 182-48-718-5588

No. of Certificate : 2020052883 Date : 2020/04/02

Certificate of Free Sales

Exporting(certifying) country : Republic of Korea
Importing-requesting country : United States of America

The Ministry of Food and Drug Safety, certifies that the following firm is authorized to
manufacture medical devices under the Medical Device Act and the following item(s) is/are
permitted to be freely sold in overseas markets.

Manufacturer (Registered No. : 5295)

PaxGenBio Co., Ltd.

#803, 361, Simin-daero, Dongan-gu, Anyang-si, Gyeonggi-do, Korea Rep.

Product-License No.	Classification
20-035	IMP Reagent for detection of nucleic acids of SARS-CoV-2 (real-time polymerase chain reaction) (in vitro diagnostic)

Manufacturer : List of Product Classification and Model

Director of High-Tech Medical Device Division
Department of Medical Device Evaluation
National Institute of Food and Drug Safety
Evaluation
Ministry of Food and Drug Safety

This certificate is issued on the internet, you can check whether to forge or modify in homepage(www.mfds.go.kr).
Furthermore, You can also check it by barcode existing document check program for scanner.

Kit Package layout



Shipping Packing Box

■ Thermal Insulation Container

- EPS (Expanded Polystyrene, **Neopor®**)
- * **Neopor**: A new material with high thermal insulation performance by mixing graphite, an inorganic material with existing EPS



■ Cold Chain System: PCM(Phase Change Material)

- * PCM is a phase change material (liquid->solid, solid->gas, etc.), A substance that absorbs heat when the surrounding temperature rises, and when the temperature decreases, it crystallizes and releases heat.

■ Carton Dimension

- Measurement: 45.5(H) x 45.5(D) x 50.0(H) Cm
- Total Weight: 21 Kg
- 40 Kits per Carton
- Shipping Stability for 7 days

Ordering Information

Product Name	Cat. No.	Package	Storage	Expiration date
PaxView® COVID-19 real-time RT-PCR Kit	R0501N	96tests/Kit	-25 °C ~ -15 °C	12 months

*Innovative and Chimeric Technologies
Leading to a New World of Healthcare*



EC	REP
----	-----

 MT Promedt Consulting GmbH
Altenhofstrasse 80, 66386 St. Ingbert, Germany

PaxGenBio Co., Ltd.

#803~807, 361, Simin-daero, Dongan-gu, Anyang-si,
Gyeonggi-do, 14057 Rep. of Korea
Customer services: (+82) 31-244-5036, info@paxgenbio.com