



TS-205

# Black Eye™

## Endoscopic **Marker**

### Indication for use

Black Eye is indicated for endoscopically marking lesions in the GI tract permanently.

Black Eye's biocompatible, non-toxic formulation contains no harmful ingredient. Black Eye consists of medical grade carbon black which lasts for more than 10 years in the body!

### Permanent Ink

### Stable Suspension

Poor suspension can make the solution separated, which leads to unstable marking!  
"No shake" Necessary = Easier, Faster, Safer!

Poor emulsion can choke the tube of the syringe or body fluid!

### Stable Emulsion

Proven in over 50 countries





# Black Eye™

## Endoscopic Marker

### STORAGE

Store the product at room temperature.

### SHELF LIFE

24 months.

### PHYSICAL PROPERTIES

- 1) Color: Black
- 2) Density: 1.05g/mL
- 3) Boiling Point: 97°C (207°F)
- 4) Freezing Point: -0.5°C (33°F)

### STERILIZATION

- 1) Steam sterilization
- 2) Sterility Assurance Level (SAL):  $10^{-6}$

### RAW MATERIALS

Water for injection  
Glycerin  
Carbon Black

### DIRECTIONS FOR USE

A 23-25 gauge sclerotherapy needle with a needle length of 4mm or less is recommended for this procedure.

### IMPORTANT

When the needle is properly positioned, insert the tip 30°-45° tangentially into the submucosal space of the colon.

### CAUTION

Inject 0.5~0.75ml of the marker into each injection site. This volume typically produces a mark on the colon measuring 1~2 cm in diameter. The maximum recommended dose is 8ml per patient.



### 5ml syringe, box of 10

Rectangular canister for secure and compact storage. Blister packaging acts as a sterile barrier and protects against cross-contamination.



### THE STANDARD CO., LTD.

120, Gunpocheomdansaneop 2-ro, Gunpo-si, Gyeonggi-do, 15880, Korea  
Tel : +82-2-838-5533 Fax : +82-2-838-5523 E-mail : sales1@thestd.co.kr