




**Signers Co, Ltd.**

**Address** A-208, Keumkang-Pentarium IT Tower, 215, Galmachi-ro, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea, 13217  
**Tel** +82-70-8808-7478  
**Fax** +82-31-754-7479  
**E-mail** libby@signers.kr, speedgun72@signers.co.kr  
**Website** http://www.signers.co.kr

**Technology and Product(s)**



SUPiA CR System offers such a better clinic environment with high resolution image quality and affordability. It is very light and fast. Our console software 'gram' is simple, fast, and compatible. You can connect with any DICOM printer and paper printer. SUPiA gram offers easy multilingual support and special tools.

- Compact & Lightweight : 436 x 636 x 196mm, 21.5kg (47.4lbs)
- High Throughput : Up to 72 IPs/Hour (14 x 17" / 160 μm)
- Extremely Simple Structure Design
- Authentic Rigid Type among Clinic CRs
- Able to bear up to 150kg of weight on cassette
- Single Cassette Feed
- 14 x 17" : 2.1kg, 10 x 12" : 1.3kg
- Single Phase 50 ~ 60Hz, AC 90 ~ 264V



**Stemlab, Inc.**

**Address** #402, Suntech City Building, 474, Dunchon-daero, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea  
**Tel** +82-70-7586-9615  
**Fax** +82-70-7586-9696  
**E-mail** sales@sm-eng.net  
**Website** http://www.stemlab.co.kr

**Technology and Product(s)**



**SUPRA ELITE (BMD)**

SUPRA ELITE is a pencil-beam densitometer using DEXA (Dual energy X-ray Absorptiometry) based on central bone type for a higher degree of accuracy, intuitive user interface and outstanding cost efficiency. In addition, the system allows users to inspect high-risk regions including spinal column, pelvis and wrist, and It provides reliable analytical results for both clinical and research purposes. SUPRA ELITE uses an auto-calibration algorithm to adjust counter values and maintain its best performance in addition to user-friendly menu driven and easy operation. The result sheet includes BMD data, other useful information and trend chart that may provide important assistance in clinical treatment.



**DEXINO (BMD)**

DEXINO is a central bone type equipment of the pencil beam type that uses dual-energy X-ray absorptiometry (DEXA), DEXINO enables a more effective analysis of results and forecasting of trends. With DEXINO, scanning particular parts, such as the spine and the femur region, would be so much easier. DEXINO is more accurate bone density measurement, which is equipped with various and effective features.



### Zeronitec Co., Ltd.

**Address** #202-602, 388, Songnae-daero, Bucheon-si, Gyeonggi-do, Republic of Korea.  
**Tel** +82-32-811-8890  
**Fax** +82-32-819-8890  
**E-mail** info@zeronitec.com, zeronitec@gmail.com  
**Website** www.zeronitec.com

#### Technology and Product(s)



#### Product portfolio

- PURE-50 (50L) : Floorstanding type
- PURE-120 (120L) : Floorstanding type, Automatic sliding door
- PURE-120D (120L): Double door system

#### Features

- Economical cost & Reliable quality product (CE & ISO) : We provide our global partners with reliable & high-end quality products at more competitiveness. Our products are in compliance with CE and use high-end quality components for our sterilizer system which give our partners and customers benefits of our product reliability.
- Safe cassette-type sterilization agent : PURE sterilization agents (H2O2) are absolutely sealed cassette type and multiple-cycle(7cycles/ cassette) cassette, it's safe to be transported and stored and convenient and safe to be treated by user.
- Graphic User interface : Most of operating menus are displayed both graphically and numerically. User can recognize easily the status of sterilizer such as sterilant cassette, process monitoring, process parameters with less human error.
- Front large process monitoring display  
 Large process monitoring LED device located in sterilizer door can provide user with more convenient use of our sterilizer system. So, user can know which cycle phase is now in process and can recognize remaining cycle time easily even from a distance..



## Physiotherapy & Home care Equipment

AMI GLOBAL CO., LTD

ATEX CO., LTD.

ChoiTechnology Co., Ltd.

Color Seven Co., Ltd.

CYBERMEDIC CO., LTD

DAEYANG MECHANICS

DAESUNG MAREF CO., LTD.

Hanmed Co., Ltd.

Homelec Korea Co., Ltd.

K1MED CO.,LTD

KTMED INC.

MAXSTAR INDUSTRIAL CO., LTD.

NEOFECT Co.,Ltd.

openM inc

Optowell Co.,Ltd.

P&S Mechanics Co., Ltd.

SEALTECH

TAEHA MECHATRONICS Co., Ltd.

(주)웰뷰텍

Won Industry.,Co.,Ltd.



**AMI GLOBAL CO., LTD**

**Address** 22, JANGNIMBEONYEONG-RO, 104BEON-GIL, SAHA-GU, BUSAN, 49465, SOUTH KOREA  
**Tel** +82-51-293-6472~3  
**Fax** +82-51-293-6474  
**E-mail** drmedkorea@gmail.com  
**Website** www.drmedbrace.com

**Technology and Product(s)**



We AMI Global Company, have been a global medical device company specialized and experienced in the field of comfortable orthopedics, rehabilitation, healthcare and sporting supports since June, 1965 with its medical line, Dr.MED and sport line. The operating system of factories is accredited by ISO13485 and its products are certified by FDA(USA) in 2003, CE mark in 2006 and GMO(medical device). Our own brand DR.MED has become various international partners for all those working at medical support/brace all over the world and will bring you improvement of person's life quality, comfort and pleasure. And then you can see more detail about our product line in the website as below.  
 (www.demedbrace.com)



**ATEX CO., LTD.**

**Address** 144, Jangjisan-ro, Gwangtan-myeon, Paju-si, Gyeonggi-do, 10955, Korea  
**Tel** +82-31-944-0546  
**Fax** +82-31-942-6182  
**E-mail** atexsales@atexmedical.com  
**Website** www.AtexSports.com

**Technology and Product(s)**




Kinesiology Tape is designed to mimic the human skin in both thickness and elasticity to allow a wide range of mobility and flexibility without restricting your body movements. When the tape is properly applied, it will lift the skin to create increased circulation of blood, lymph, and tissue fluids which results in faster relief of muscle and joint pain. It is becoming one of the must have essential tools in today's pain management and preventive practices.



**ChoisTechnology Co., Ltd.**


**Address** #S-1404, Songdo Technopark IT center, 32, Songdogwahak-ro, Yeonsu-gu, Incheon, Korea (21984)  
**Tel** +82-32-346-3409  
**Fax** +82-32-246-3406  
**E-mail** inquiry@choistec.com  
**Website** www.choistec.com

**Technology and Product(s)**



Smart monitoring thermometer 'Thermosafer(XST200)':

- (via dedicated app) It has real-time temperature monitoring and abnormal body temperature alarm and it can record the temperature measurement result and check history with graph
- More stable and accurate temperature measurement is possible comparing with Non-contact type thermometer (it's error range is  $\pm 0.1^{\circ}\text{C}$ )
- Superior measurement distance by applying Bluetooth technology (max 15m)
- Using low irritant silicon patch for sensitive skin.



**Color Seven Co., Ltd.**

**Address** 4F, Geoam Bldg, 62, Banpo-daero 4-gil, Seocho-gu, Seoul, Korea  
**Tel** +82-2-3462-5618  
**Fax** +82-2-3462-5619  
**E-mail** sales@color7.co.kr  
**Website** www.color7.co.kr/eng

**Technology and Product(s)**



**VPOD - New Light Therapy Patch**

VPOD is a unique healthcare light therapy device which is consumable model. It is a skin adhesive Light Therapy Device based on the Low Level Light Therapy technology. This is an effective and safe electric device for women & men's health and also for beauty purpose. The principle is called PAMS technology (Color Seven's own intellectual property). Light from the device relaxes the smooth muscles in the tissue and recovers blood flow.



**Woman Stress Free (COLOR DNA WSF)**

Woman Stress Free is the first and unique in the world for Dysmenorrhea treatment medical device. It is a skin adhesive Color Light Therapy based on the Low Level Light Therapy technology that proved scientifically the principal of treatment. This is an effective and safe medical device for women who have painful menstruation stress due to PMS (PreMenstrual Syndrome) & dysmenorrhea. The principle of treatment is based on PAMS technology (Color Seven's own intellectual property) to relax the smooth muscles in uterus tissue. Women Stress Free was medically verified through clinical evaluation with menstrual pain patients in two University hospitals and the result has been published in Springer, Vol. 286:947-952, also this evaluation presented "Issues in Gynecology, Obstetrics, Fertility and Pregnancy Research in 2013 edition"





**CYBERMEDIC CO., LTD**

**Address** 203/1, Standard factory, 20-23, Yackchon-ro 8-gil, Iksan-si, Jeollabuk-do, Republic of Korea  
**Tel** +82-63-835-5114  
**Fax** +82-63-834-0154  
**E-mail** Daniel.lee@cybermedic.co.kr  
**Website** www.cybermedic.co.kr

**Technology and Product(s)**



MotoCog is Interactive ADL & Cognitive rehabilitation training system for elderly, hypofunction of ADL, damaged of cognitive function, central nerve injury and dementia patients



An instrument for analyzing and strengthening the anterior and posterior movement of the tongue to overcome dysphagia. Improved dysphagia by increasing strength and coordination of tongue, which is important for food movement also pressure making and holding, when swallowed. An instrument for measuring and strengthening laryngeal muscle resistance for overcoming dysphagia.



A System of exercise for muscle strength, muscular hypertrophy, muscle endurance, pain modulation and muscle coordination training according to measurement of maximum muscle strength using electronic pulley(Apply center of professional sport and musculoskeletal center and physio therapy exercise).



The FES for dysphagia patient who cannot easily swallow food and liquid, has a built-in treatment pattern and stimulate the muscles necessary to swallow. Exclusive pattern helps to exercise a swallowing and to reach a normal swallow.



**DAEYANG MECHANICS**

**Address** DaeYang Mechanics 32 Muhak-ro 35 Gil , Dongdaemun-gu, Seoul 130-823 South Korea  
**Tel** +82-2-929-2222  
**Fax** +82-2-928-4411  
**E-mail** info\_hubex@naver.com  
**Website** www.hubex.com

**Technology and Product(s)**



**3D Care Exercise System for the Pelvis, Back and Abdomen**

A multi rehabilitation system that trains the core area with the extensor and refractive exercise on the hip joint and the lumbar joint while completely relaxing the body muscles. As the upper body maintains zero weight by the gravity offset program, the neck, shoulder and back muscles do not become tense but relaxed. The back and the knees are safely protected so the elderly and patients can easily respond to the musculoskeletal disorders. Hospital, public health, occupational health, women's fitness and silver care



**DAESUNG MAREF CO., LTD.**

**Address** 298-24, Gongdan-ro, Gunpo-si, Gyeonggi-do, Korea.  
**Tel** +82-70-8740-6056  
**Fax** +82-31-459-7215  
**E-mail** medical@dsmaref.com  
**Website** www.dsmaref.com

**Technology and Product(s)**



**DTS-3000 (Smart Tourniquet System)**

We are providing the highest functional pneumatic tourniquet systems in market. It provides suitable functions for an operating room under the optimized conditions. (Recording, IOP sensing, LCD Touch Screen, Safety Alarm System, Battery mode).



**DVT-2600 (Venous Assist System)**

DVT-2600 IPC system can be used individually as an alternative medicine for a stroke patient, or a patient who cannot use orally administer drugs, or any patient who cannot use anticoagulants such as heparin or warfarin, to prevent DVT and PE(Pulmonary Embolism)



**Hanmed Co., Ltd.**

**Address** 80-121, Golden Root-Ro, Juchon-Myeon, Gimhae-Si Gyeongsangnam-Do, KOREA  
**Tel** +82-55-331-0575  
**Fax** +82-55-331-0547  
**E-mail** kinetrac@naver.com  
**Website** www.kinetrac.com  
 www.hanmed.net

**Technology and Product(s)**



**Decompression System Kinetrac DAVINCI**

- Traction With Twist  
Traction With Twist directly focused on the Axis of Spine
- Traction With Target Lordosis a new-concept Kinetic Decompression technology have created the new theory of disc decompression.
- Traction With Muscle relaxation therapy  
Pelvic joint, Lumbar Muscle Flexion/Extension.



**3D Gravitation Muscle Training and evaluation System 3D NEWTON**

- 3D Newton is a new conception of exercise system for rehabilitation with 3D Gravitation Training an system. It is possible of 3D space exercise both 360° rotation and 0° ~ 60° inclination control.
- Real-time Bidirectional Monitoring training system



**Homelec Korea Co., Ltd.**

**Address** A-522, Tera Tower, Songpa-daero 167, Songpa-gu, Seoul, Republic of Korea  
**Tel** +82-2-521-2288  
**Fax** +82-2-2054-3588  
**E-mail** info@homelec.co.kr  
**Website** www.homelec.co.kr

**Technology and Product(s)**



**CM07 : The Best Premium Cordless Massager**

- Dual motor : powerful percussion motor with 3,700RPM, gentle micro-vibration motor with 11,000RPM
- Adjustable massage strength with dial button
- Patented air cushion stick for joints and bones
- Slim design : easy to maneuver over the body with an extra long handle for hard to reach area
- Rechargeable with quality Li-ion battery : 180 minutes continuous use
- 780g, ultra-light weight for easy use
- 6 Multiple massage tips



**K1MED CO.,LTD**

**Address** Rm.205, Woolim E-Biz Center 2, 12, Digital-ro 33 gil, Guro-gu, Seoul, 08377, Korea.  
**Tel** +82-2-871-0657  
**Fax** +82-2-871-0658  
**E-mail** overseas@k1med.com  
**Website** www.k1med.com

**Technology and Product(s)**



LUBISONE is new technology to recover lesion in body tissue.

It promotes revascularization and removing calcareous near tendons or joints though physical changing of cell membranes while stimulates musculoskeletal system where has lesion by shock wave.

- Name of product : LUBISONE
- Model Number : K1W-1501
- Rated Voltage : 230V~, 50/60Hz
- Power Consumption : 150VA
- Weight (Main unit only) : 10kg
- Dimension (Main unit only) : 397\*370\*200mm





**KTMED INC.**

**Address** #602, Byucksan Digital Vally 1-cha, 24, Digital-ro  
27-gil, Guro-gu, Seoul, Korea  
**Tel** +82-2-2108-7890  
**Fax** +82-2-2108-7894  
**E-mail** sales@ktmed.com  
**Website** www.ktmed.com

**Technology and Product(s)**



**NePlus(NE-SM1)**

SUPER MESH NEBULIZER “Clinical Test Completed!”  
(Doi: 10.1111/bcpt.12499)

- Most strong mesh(No metal dust & No hole block)
- MMAD 5µm particles sprays more than 55%
- Noiseless operation
- Preventing infection by boiling or autoclave for Mask, Mouth Piece
- No disposable kit(Reusable mask, medication bottle, battery)
- Very portable size with rechargeable Li-Ion Battery
- Residual volume is under 1% by world patented technology



**AnyNeb(NE-SM2)**

SUPER MESH NEBULIZER “Simple is the Best”

- Most strong mesh(No metal dust & No hole block)
- MMAD 5µm particles sprays more than 55%
- Noiseless operation
- Preventing infection by boiling or autoclave for Mask, Mouth Piece
- No disposable kit(Reusable mask, medication bottle)
- Very portable size with normal Alkaline batteries.



**MAXSTAR INDUSTRIAL CO., LTD.**

**Address** 152-12, HWANGGEUM-RO 23 BEON-GIL,  
YANGCHON-EUP, GIMPO-SI, GYEONGGI-DO,  
KOREA, 10048  
**Tel** +82-31-989-3543  
**Fax** +82-31-989-3548  
**E-mail** overseas@maxstar.co.kr  
**Website** www.maxstar.co.kr

**Technology and Product(s)**



UAM-12000: BEAUTY USE



UAM-8500 : MEDICAL USE



UAM-8100 : HOME USE



UAM-8300N : HOME USE



UAM-9306-BK: SPORTS USE



UAM-9406-BK: SPORTS USE





**NEOFECT Co.,Ltd.**

**Address** 401 West Hall, Dankook University, 152, Jukjeon-ro, Suji-gu, Yongin-si, Gyeonggi-do, Korea, Republic of  
**Tel** +82-31-889-8521  
**Fax** +82-31-8005-4121  
**E-mail** contact.eu@neofect.com  
**Website** www.neofect.com/en

**Technology and Product(s)**



**RAPAEL Smart Glove**

The wearable device-based digital solution for hand rehabilitation




**RAPAEL Smart Kids**

The wearable device-based digital rehabilitation solution for pediatric patients



**RAPAEL Smart Board**

The digital solution for upper limb rehabilitation



**RAPAEL Smart Pegboard**


Digital pegboard for functional and cognitive rehabilitation



**openM inc**

**Address** 413, Expo-ro, Yuseong-gu, Daejeon, 34051, Republic of Korea  
**Tel** +82-42-823-7301  
**Fax** +82-42 823 7302  
**E-mail** sales@openm.com  
**Website** www.openm.com

**Technology and Product(s)**



OPENCAST is a next generation, shower-friendly support frame.

**Benefits of OPENCAST**

**Shows with cast on**

OPENCAST allows daily washing and showers with cast on. Water can be easily dried off with a towel and a hair dryer.

**Odor and Itch free**

Perspiration quickly evaporates to prevent itchiness and unwanted odor.

**Visual check**

Doctors and patients can visually check the condition of the skin through openings with cast on.

**Repeated wear and removal**

OPENCAST can be repeatedly opened and be reworn as needed.

**Stylish design**

OPENCAST is made of fine materials and has a stylish design.

**Comfortable and Hygienic**

OPENCAST is Comfortable and Hygienic.



**Optowell Co.,Ltd.**

**Address** Optowell Co.,Ltd.  
**Tel** +82-63-214-1186  
**Fax** +82-63-214-1186  
**E-mail** gabriel@optowell.com  
**Website** www.comalgum.com

**Technology and Product(s)**



Comalgum is the rhinitis curer applying Low Level Laser Therapy(LLLT) which is not harm to human body. Comalgum is able to cure trouble in nasal cavity that is the root cause of rhinitis and intensified trouble in nasal inside in the meantime because it emit over 2 wave-length of laser from the laser detector. Comalgum makes easy to breath by making the breath hole in the nasal inserted area.

Comalgum is recommended for

**Children**

Children's health and fitness may be impaired if they suffer from frequent colds which can lower their immunity to further infections and weaken their constitution. Children's normal growth may be adversely affected by constant nasal infections and consequent breathing problems. Nasal congestion can inhibit concentration and lead to developmental and concentration problems.

**People who taking examinations**

Rhinitis, hay fever, runny noses and excessive sneezing are all well-documented causes of memory and concentration problems. The resulting fatigue also impairs exam performance.

**Office workers**

The office environment, particularly air-conditioned offices is a well-known source of nasal problems and infections. Acute rhinitis and allergic rhinitis frequently develop as a result of colds caught in the office



**P&S Mechanics Co., Ltd.**

**Address** Rm. 509, Woori Venture Town II, 70, Seonyu-Ro, Yeongdeungpo-Gu, Seoul, South Korea (zip) 07294  
**Tel** +82-2-2632-3084  
**Fax** +82-2-2632-3085  
**E-mail** walkbot@walkbot.co.kr  
**Website** www.walkbot.co.kr

**Technology and Product(s)**



Walkbot is robot-assisted gait rehabilitation training system for locomotion-impaired patients who are suffered from stroke, SCI, TBI, multiple sclerosis and cerebral palsy, etc. Products are lined up Walkbot\_S (for adult), Walkbot\_K (for pediatric), Walkbot\_G (for both adult and pediatric, module type) and recently launched Walkbot\_Premium (for adult) having the most advanced features of Walkbot series. Based on effectiveness on the training and the excellent quality of Walkbot, it has been installed in USA, Europe, Russia, Asia region and Korea.

- The world-first activated hip/knee/ankle joint drive motors based on kinematic study makes accurate and natural gait pattern close to the human being.
- Individually customized gait pattern maximizes its clinical effects
- More sessions and efficient therapy by automatic leg length adjustment
- Useful training mode (interactive, assisted training mode)
- Auditory(Metronome), visual(Front Camera), tactile(Game/Sport) stimulation for the neuroplasticity encouragement
- Various assessment tools (ROM, Stiffness, Force)
- Real-time biofeedback through augmented reality software
- CE, FDA, MET, CFDA, ISO 13485, ISO 9001 certified
- Clinically proven equipment





**SEALTECH**

**Address** 51-13 Myobong-ro 173beon-gil Idong-eup  
Cheoin-gu Yongin-si Gyeonggi-do  
**Tel** +82-31-431-7915  
**Fax** +82-505-009-7916  
**E-mail** evoss@esealtech.com  
**Website** www.evoss.co.kr

**Technology and Product(s)**




1. Brand : EVOSS
2. Category : Anti-bedsore cushion and mattress
3. Raw material : 100% natural rubber
4. Type : Air-cell cushion and mattress
5. Product
  - 5-1. Cushion
    - 5-1-1. Octagonal 8\*9 380\*420\*60mm
    - 5-1-2. Octagonal 7\*8 460\*410\*80mm
    - 5-1-3. Octagonal 7\*7 410\*410\*80mm
    - 5-1-4. Hexagonal 8\*9 420\*380\*60mm
    - 5-1-5. Hexagonal 8\*9 450\*410\*60mm
    - 5-1-6. Hexagonal 8\*9 510\*490\*61mm
    - 5-1-7. Hexagonal 8\*9 510\*490\*87mm
  - 5-2. Mattress : Octagonal 810\*1840\*80mm
6. Incoterms : EX-WORK, FOB, CIF, DAP, DDP
7. Payment : T/T in advance, LC
8. Price : The price depends on the quantity.
9. Country of origin : Made in Korea
10. Address : 51-13 Myobong-ro 173beon-gil Idong-eup  
Cheoin-gu Yongin-si Gyeonggi-do



**TAEHA MECHATRONICS Co., Ltd.**

**Address** 421, Bakdal-ro, Manan-gu, Anyangsi, Gyeonggi-do  
**Tel** +82-31-463-7989  
**Fax** +82-31-463-0529  
**E-mail** petergo@taeha.co.kr  
**Website** www.stexfitness.com

**Technology and Product(s)**



**Simple and User Friendly Design**

- Easy to use
- 19 cm height of step
- Guide hand rail for safety
- Track-Master protocol with RS232 connection

**Main Specification**

- 2.5HP Treadmill Duty AC Motor
- Speed: 0.8 Km/n ~ 25 Km/h
- Elevation: 0% ~ 25%
- Maximum User Weight: 180 Kg
- Warranty: 3 years
- Dimension: W 1,010mm \* L 2,100mm \* H 1,150mm
- Running Space: W540 mm \* L1,500mm
- Power: AC110 or AC220V / 50/60Hz
- Product Weight: 170 Kg



주월뷰텍

**Address** 14, Baekbeom-ro 677 beon-gil, Seo-gu, Incheon, Korea  
**Tel** +82-32-870-0934  
**Fax** +82-32-870-0939  
**E-mail** sales@welbu.com  
**Website** www.Welbu.com

Technology and Product(s)



Seven Liner ZAM ZAM will make you most beautiful body line!  
 Seven Liner ZAM ZAM is a limb cyclical pressure device, helping to prevent blood clots in blood vessels with 5 chambers. In order to help increase blood circulation, periodically, high air pressure pumping to inflate the limbs and sleeves prevents swelling and daily fatigue. Especially, Hose Plug Sensor Switch (cuff connection recognition) and Cuff Deflation function (removal of remaining air in the cuff after use) are our own technology patented in Korea.



Won Industry.,Co.,Ltd.

**Address** 7-12,7-6,Bangadari-gill, Siheung-si,Gyeonggi-Do, KOREA  
**Tel** +82-31-404-3093  
**Fax** +82-31-318-5573  
**E-mail** lead\_care@naver.com  
**Website** www.lead-care.com

Technology and Product(s)




Morning Life

- Relaxes, detoxifies, and rejuvenates; but it can also be used as a pre-workout to loosen up muscles and a post-workout to help your muscles recover by ridding them of lactic acid build-up
- Digital or analog type
- Time Range : 0~60min
- Dimension : 280x218x172mm
- Weigh : 4.2Kg




FoxyL

- Relaxes, detoxifies, and rejuvenates; but it can also be used as a pre-workout to loosen up muscles and a post-workout to help your muscles recover by ridding them of lactic acid build-up
- Dimension : 210x200x125mm
- Anti-overheating function
- Weigh : 3.5Kg
- Comfortable touch display function



Thor

- Rated voltage/Power consumption : 8.4V, 1A / 13W
- Mode : Stage 1 to 5 Dynamic Control
- Operation Time : 15 Minutes Operation Auto Timer
- Configuration : Body + Self massage toolX3
- Three types of replaceable safety bars are available in a variety of areas.
- Easy portable wireless massager



Surginet


- \*Feature and Benefit
- Unwoven strands
- Washable and Reusable
- Simple to use especially for head, breast, shoulder, waist and thigh
- Hold the gauge placed on the affected part comfortably even when you move.
- Free cutting for your purpose



## Surgical & Medical Supplies

ACMEDI KOREA Co., Ltd  
AMT Engineering Co. Ltd.  
DAEJUMEDITECH  
ENGINEERING CO.,LTD  
DAIHAN LABTECH CO., LTD.  
DAIWHA CORP., LTD.  
DENTIS CO.,LTD.  
ENDOVISION CO.,LTD.  
Everaid Co., Ltd.  
EWHA BIOMEDICS Co., Ltd.  
Firson Co., Ltd  
GENOSS  
IDS Ltd  
IMEDICOM Co., Ltd  
LIVSMED Inc.  
LOWTEM Co., Ltd.  
MEDICALPARK CO., LTD.  
MEDONICA

MEINNTECH CO., LTD.  
Nanoomtech Co., Ltd.  
Renosem Co.,Ltd.  
SAEUM MEDITEC  
SAMSUNG MEDICAL ENG. CO., LTD.  
Samyang Biopharmaceuticals Corporation  
ScienceMedic Co., Ltd.  
Sejong Medical  
SM ENG. CO., LTD.  
StraTek Co.,Ltd.  
SUNGKWANG MEDITECH CO.,LTD  
SurgiCore Co., Ltd.  
TaeWoong Medical  
VOTEM CO., LTD  
Wonbiongen Co., Ltd.  
WONTECH CO., LTD.  
YOUNG CHEMICAL CO., LTD.  
ZERONE CO., LTD.



## ACMEDI KOREA Co., Ltd

**Address** #611, Hyundai 41 Tower, 293, Mokdong dong-ro, Yangcheon-gu, Seoul, Korea

**Tel** +82-70-5158-3446

**Fax** +82-2-6971-9764

**E-mail** sales@acmedikorea.com

**Website** www.acmedikorea.com


### Technology and Product(s)



#### VegaPlax

Advanced Plasma Technology for Medical Purpose. It is the ideal wireless Plasma device made in Korea for Blepharoplasty with Plasma Arc. It will help patients to improve eye ptosis without surgery. Moreover, it can be used for skin lifting by delivering plasma energy into the skin tissue.

- Output Frequency : 70kHz
- Power Level 3 Steps (0.7W, 1.0W, 2.0W)
- Output Mode : Continuous & Pulse
- Pulse Mode Frequency 10, 20, 40, 100, 200Hz
- Running Time : When Full charged, Max 3 hours
- Weight : Body 130g
- Size : 197 x 40 x 38mm
- Screen : 1.5" Full Color OLED
- This product is guaranteed for 1 year from the date of purchase.



## AMT Engineering Co. Ltd.

**Address** #811, Byoksan Technopia, 560 Dunchon-daero, Jungwon-gu, Seongnam-si, Gyeonggi-do, 13230 KOREA

**Tel** +82-31-737-8066

**Fax** +82-31-737-8069

**E-mail** wikim@amtpro.co.kr

**Website** www.amtpro.co.kr

### Technology and Product(s)



#### Lucas plus: Q-Switched ND:YAG Laser System

LUCAS plus is the best Tattoo removal & Laser Toning device.



#### Olive plus: 810nm Diode laser device

OLIVE PLUS is the best hair removal laser device.



#### A-TONE: "PREMIUM" Multi-Waved(I2PL) + PIP Toning System

A-TONE is multi-functional aesthetic devices for toning, hair removal, pigment and vascular lesions.



#### BOTOSONIC PLUS: Aesthetic HIFU Device

Botosonic is best Face Lifting & Tightening.



## DAEJUMEDITECH ENGINEERING CO.,LTD

**Address** (02066) 487,Mangu-ro, Jungrang-gu, Seoul, Korea

**Tel** +82-2-2208-0905

**Fax** +82-2-438-0905

**E-mail** infilux@daejumedi.co.kr

**Website** <http://www.infilux.co.kr>

### Technology and Product(s)



#### Co2 Fractional Laser

It is an ULTRA PULSE CO2 FRACTIONAL LASER , which is a new generational version of the CO2 Fractional laser and its wave length is 10600nm, which was replaced by a stronger laser for severe scars, enlarged pores etc. It is Ultra Pulse CO2 Fractional Laser, which is a new generational version of the CO2 Fractional laser for severe scars, enlarged pores, and various kinds of skin disorders.



#### Diode Hair removal Laser 600W

- What is hair removal laser?  
Our Diode Hair removal Laser 600w has high energy pulsed 808nm diode laser for hair removal. Aroma machine can reduce pain while reducing skin damage by Cooling effect to reduce skin's pain and can get perfect result by quick and simple treatment of hair removal.



#### HIFU for face lifting and tightening and body

- What is HIFU?  
HIPRO is HIFU machine that has various cartridges and HIPRO delivers HIFU energy to SMAS layer which is located between subcutaneous tissues and muscle. And then, make 0.7~1.0mm size of TCP(Thermal Coagulation Points) zone without surgery or any disrupting skin surface to stimulate collagen growth.



## DAIHAN LABTECH CO., LTD.

**Address** 247-16 GAGOK-RI, HWADO-EUP, NAMYANGJU-CITY, KYONGGI-DO, KOREA

**Tel** +82-31-593-5961

**Fax** +82-31-593-9563

**E-mail** [Info2@labtech.co.kr](mailto:Info2@labtech.co.kr)

**Website** [www.labtech.co.kr](http://www.labtech.co.kr)

### Technology and Product(s)



#### Bench Top Autoclave (Class B)

With automatic dual door & LCD monitor display for whole procedure.  
Available capacity: 12/25/40 liters



#### Medium & Large autoclave

With pre & post vacuum system & strong hatch door for user safety  
Available capacity: 40/60/80/100/200/300 liters



#### Pass through door autoclave (Class B)

With full touch screen & round/square chamber selectable  
Available capacity: 200/300/500 liters



#### Plasma Sterilizer

With full touch screen & square chamber, it's suitable for sterilizing medical tools/equipments that are weak at heat/high temperature such as endoscopes, catheters, etc.

**DAIWHA**  
www.daiwha.com



### DAIWHA CORP., LTD.

**Address** Daiwha Bldg, 3, Yeoksam-ro 33 gil, Gangnam-gu, Seoul, South Korea  
**Tel** +82-2-558-1711  
**Fax** +82-2-554-0317  
**E-mail** overseas03@daiwha.com  
**Website** www.daiwha.com

#### Technology and Product(s)



##### DI-6000/DS-5000 (Open/Close System)

- Innovative Design
- High Convenience – TFT-LCD, Jog Dial
- Various Safety Functions : IV Set Loading(Infusion), Easy Menu, Auto Syringe Detection(Syringe)
- 1~1500ml/h Range of flow rate
- Smart Docking Station System using multiple pumps(~4 units).
- Battery : Operating 10 hours, Recharging 6 hours
- Weight : Approx. 1.9(±0.05)kg



##### DE-400 MEDITOM RF Generator with APC

- Modern Design : improved convenience LCD touch screen
- Various Safety Function : D.R.M(Daiwha patient Return plate monitoring system)
- Various Output mode
  - 2 modes of Bipolar Cutting(APC-General Cut)
  - 2 modes of Bipolar Coagulation(APC-Standard, Low)
  - 4 modes of Mono-polar Cutting(Pure,Blend 1/2/3)
  - 3 modes of Mono-Polar Coagulation(Forced, soft Low)
- 2 Monopolar and 2 Bipolar Plug with play instrument



##### DUB-100 & DB-137/ DU137/127/Scissor

- DUB-100 : 1 Ultrasonic Plug, 1 Bipolar Plug, 1 Monopolar Plug
- DU-137 : Bipolar Vessel Sealing Energy Device
- DU-137 : Ultrasonic Energy Device
- DU-127/Scissor : Ultrasonic Open-Type Energy Device  
- At Target Hospitals- early 2019

**LUVIS**



### DENTIS CO.,LTD.

**Address** Seong Seo Seo Ro 99, DalSeo Gu, DaeGu, Korea  
**Tel** +82-2-869-6020  
**Fax** +82-2-869-6021  
**E-mail** luvis@dentis.co.kr  
**Website** http://luvis.co.kr/eng/index.php

#### Technology and Product(s)



##### LUVIS-L200(Surgical light)

- Intensity : 160,000 Lux
- Color Temperature : 3800K / 4300K / 4800K
- CRI : 96



##### LUVIS-M200(Surgical light)

- Intensity : 120,000 Lux
- Color Temperature : 4300K
- CRI : 95



##### LUVIS-S200(Minor surgical light)

- Intensity : 100,000 Lux
- Color Temperature : 3800K / 4300K / 4800K
- CRI : 90



##### LUVIS-E100(Examination light)

- Intensity : 60,000 Lux
- Color Temperature : 3500K / 4300K / 4800K
- CRI : 90





## ENDOVISION CO.,LTD.


**Address** #214, 2nd Floor MJ Technopark, 20, Dalseo-daero 109-gil, Dalseo-gu, Daegu, 42709, Korea  
**Tel** +82-53-751-1955  
**Fax** +82-53-751-1964  
**E-mail** braelyn@endovision.co.kr  
**Website** www.endovision.co.kr

### Technology and Product(s)




#### NEXOSEAL

Hemorrhagic hemostasis product during general surgery, easy use during operation and securing visibility, fast hemostatic effect and excellent prevention of additional hemostasis



#### DUAL PATTY

Bleeding management during spinal or brain surgery enables quick and accurate management of bleeding during surgery and easy protection of surrounding organs



#### UBE Tool kit

A set of bi-directional spinal endoscopic instruments for various spine surgery.



#### Flex-car R/C

RACZ catheter used for neuroplasty (PEN) with excellent pushability and trackability




#### TOME FORCEP

Rongeur used to remove bones or remove tissue during spinal surgery



#### FREETOME PUNCH

It is a punch used for removal of bone and tissue during spine surgery. It is convenient to rotate due to shaft rotation



## Everaid Co., Ltd.

**Address** 36, Dulleom-gil, Yangsan-si, Gyeongsangnam-do, Korea  
**Tel** +82-55-363-2122  
**Fax** +82-55-363-4683  
**E-mail** hq@everaid.com  
**Website** www.everaid.com

### Technology and Product(s)



**AIDFoam is HDI type hydrophilic foam dressing.** This HDI foam dressing has a high absorption & retention capacity with its special fluid-lock structure. The foam can absorb fluid about 10 to 14 times of its own weight (EN13726-1), and minimize exudates leaking. The high MVTR top film, about 2000~2500g/24hr (EN13726-2), increase exudates control capacity.

On top of that, AIDFoam comfort, silicone adhesive foam dressing has very high tack adhesive silicone, adhering well on body contours. The silicone adhesive provides long wear time and minimizing pain and damage in changing dressing.



**AIDerm is high MVTR PU transparent film dressing.** The sterile thin and non-glossy film dressing can perfectly protect skin and wound from outside contaminants such as bacteria, gem and fluid to prevent 2nd infection.



**NeoDerm is high tack adhesive hydrocolloid dressing.** This high tack adhesive make the dressing adhere well to skin and minimize skin trouble compared to common hydrocolloid dressings. NeoDerm Changes absorbed exudate to gel promptly to create optimal moist wound environment. The high MVTR top film increases exudate control capacity as well.



## EWHA BIOMEDICS Co., Ltd.

**Address** #901, 9F, Sungwoo BLDG, 49, Mapodaero, Mapo-gu, Seoul, South Korea

**Tel** +82-2-2055-3820

**Fax** +82-2-2055-3819


**E-mail** trade@ewhainc.co.kr

**Website** [http://www.anapa.com/en\\_index.php](http://www.anapa.com/en_index.php)

### Technology and Product(s)



disposable infusion device ANAPA & ANAPLUS  
ANAPA is the world's first and unique disposable infusion device using self-produced CO<sub>2</sub> as power source. By its unique design, it is safe of the many shortcomings of mechanical, electronic and elastomeric devices. It offers the best combination of accuracy, safety and convenience in pain palliation.



- 1ml graduation for easy monitoring of infused medication volume.
- Negligible amount of medication remaining in device after infusion.
- Patented end-cap that allows for safe and time saving self-priming.



ANAPA CTx - Chemotherapy is the world's first and unique disposable infusion device using self-produced CO<sub>2</sub> as power source. By its unique design, it is safe of the many shortcomings of mechanical, electronic and elastomeric devices. It offers the best combination of accuracy, safety and convenience in chemotherapy.

- Unique filter to prevent crystallization of chemotherapy's medication.



## Firson Co., Ltd

**Address** 47, Baekseokgongdan 1-ro, Seobuk-gu, Cheonan-si, Chungcheongnam-do, Korea

**Tel** +82-1833-8200

**Fax** +82-2-2138-8280

**E-mail** global@firson.co.kr

**Website** [www.firson.co.kr](http://www.firson.co.kr)

### Technology and Product(s)



- Use life : Up to 14 days
- Shelf life : 24 months
- Storage information : 15~30 °C
- Open bottle shelf life : 75 days

- Rapid disinfection (5min)
- Proven safety & Efficacy (CE approval)
- High reusability
- Cost-saving

**GENOSS**  
For Patients & Doctors



## GENOSS

**Address** 1st floor, Gyeonggi R&DB Center, 105,  
Gwanggyo-ro, Yeongtong-gu, Suwon-si, Gyeonggi-  
do, South Korea

**Tel** +82-10-3348-9391

**Fax** +82-31-888-5136

**E-mail** yrchoi@genoss.com

**Website** www.genoss.com

### Technology and Product(s)



#### GENOSS MONALISA Filler

MONALISA is a high purity hyaluronic acid filler which keeps the excellent volume effect for a long time at the beginning of the procedure, being decomposed naturally over time. Depending on particle size, you can choose appropriate products from soft, mild to hard type.



#### GENOSS Inflator B30

- Easy fixing and releasing with the button type Ergonomic design for comfortable handling
- 20cc barrel allows rapid and easy deflation



#### GENOSS Inflator B40

- Easy fixing and releasing with the button type Ergonomic design for comfortable handling
- 30cc barrel allows rapid and easy deflation



#### GENOSS PTCA Balloon Catheter

- Great pushability and trackability with PTFE coated hypotube
- High crossability with ultra-thin profile
- Specially engineered tri-wrapped balloon for excellent re-wrapping

**ids**  
Innovation for better life



## IDS Ltd

**Address** #105, Byoksan Digital Valley 3 cha, 271, Digital-  
ro, Guro-gu, Seoul, Korea

**Tel** +82-2-2025-1151

**Fax** +82-2-2025-1154

**E-mail** info@ids.co.kr

**Website** www.ids.co.kr

### Technology and Product(s)



#### SmaXel

SmaXel is a scanning type fractional ultra-pulse CO2 laser, in which the surgeon can control freely Scan Area, Density Level and Beam Degree. SmaXel is an advanced ultra-pulse CO2 laser, which increases treatment effects while minimizing the influence to neighboring tissues as a fine laser beam is irradiated.



#### Q-Switched Nd:YAG laser (Max : 1.8J)

Application : Tattoos, Vascular lesion Removal, Nevus of OTA, Freckle, Nano peeling, Quasi long peeling, Hair removal using long pulse. Etc



#### Alexandrite

This hair removal system is a light therapy device specifically designed to eliminate undesired hair from all parts of the body.

- Easy operation
- Less painful
- Treatment: Arm, Beard, Breast, Armpit, Groin, Leg



#### Bu:tera

High-Intensity Focused Ultrasound (HIFU) is a highly precise medical procedure that applies high-intensity focused ultrasound energy to locally heat and destroy diseased or damaged tissue through ablation.





## IMEDICOM Co., Ltd

**Address** #612, 172 LS-ro, Gunpo-si, Gyeonggi-do,  
Republic of Korea

**Tel** +82-31-479-1156

**Fax** +82-31-479-1157

**E-mail** leejames@imedicom.co.kr

**Website** www.imedicom.co.kr

### Technology and Product(s)



#### Balloon Kyphoplasty System

MEDINAUT/MEDINAUT Plus Balloon Kyphoplasty System is intended to be used as conventional bone tamps for the reduction of Vertebral Compression Fractures (VCFs) and/or creation of a void in cancellous bone in the spine during Balloon Kyphoplasty.



#### Epidural Catheter

For patients who are afraid to take spine surgery against disc disease or spinal stenosis, SPINAUT-E & SPINAUT-V are good alternative surgical procedures for patients and surgeons. At lower spinal column, SPINAUT-E & SPINAUT-V treat pains with lysis and delivering pain medication.



#### Powered Acetabular Cup Remover

EZX helps you remove an acetabular cup safely, precisely, easily and efficiently. You can extract a cup with less bone, tissue loss and prevent the creation of additional acetabular fracture.



#### Surgical Power Handpiece

IMEDICOM has two kinds of power tool system, ZAGUAR and Leopard. ZAGUAR powered instruments system is battery operated system and Leopard power tool system is pneumatic system which is operated by compressed air or nitrogen. These two system is intended for Trauma and Joint replacement surgery.



## LIVSMED Inc.

**Address** 304, D-dong, Pangyo-ro 700, Bundang-gu,  
Seongnam-si, Gyeonggi-do 13516, Republic of  
Korea



**Tel** +82-70-4450-2717

**Fax** +82-31-706-3211

**E-mail** saebomi\_lee@livsmed.com

**Website** www.livsmed.com

### Technology and Product(s)



Since the existing general laparoscopic instruments have no joint structure, instruments cannot approach at certain angle such as side of tissue or lesion. In other words, motions are very limited during operation and surgeons need more time to operate.

- It is newly developed the multi-joint and articulating laparoscopic surgical instrument that can overcome the drawbacks of a straight-type device and allow intuitive manipulation. This product is a manual laparoscopic surgical instrument with intuitive multi-joint function. Product name is <ArtiSential>.
- Due to the multi-joint structure of end tool, <ArtiSential> can move 90 degrees to upward, downward, leftward and rightward. The movement of the end tool, the movement of the control part, and the movement of user's wrist and fingers are exactly synchronized. Because of these advantages, intuitive articulating movement can be achieved.
- And <ArtiSential> allows access to a narrow surgical site and allows access to a variety of angles that conventional laparoscopic surgical instruments can not.





### LOWTEM Co., Ltd.

**Address** #A-301, 302, Gunpo IT Valley, 17, Gosan-ro  
148beon-gil, Gunpo-si, Gyeonggi-do, Korea

**Tel** +82-31-427-0763

**Fax** +82-31-427-0793

**E-mail** lowtem@lowtem.com

**Website** www.lowtem.com

#### Technology and Product(s)



LOWTEM Series plasma sterilizer Crystal/Smart 30 (33Liters) / Crystal 40R (40Liters) / Crystal 50 (50Liters) / Crystal 100 (100Liters) /Crystal 120M(D) (130Liters) /Smart 150(D) (150Liters) LOWTEM plasma system uses low-temperature, hydrogen peroxide and plasma technology to sterilize a wide range of medical instrument effectively and safely.

- Rapid Warmup & Dry System : The world's first patented allows to overcome the humidity related problem, removing residual moisture and chamber's cold point.
- Convenient and Safe Sterilizing Agent : Multi-cycle "Tank type" sterilant is safe for users to store, to deliver, and to use. The automatic RFID system for using sterilant tank is convenient and reliable for users.
- Auto Interlock & Open System : One gentle push makes the chamber door locked, and the chamber door can be opened automatically when the door sensor detects finger movement.
- Plug and Play : Warming up time is very rapid (Maximum 15 min) after main power is on, so that standby power consumption is zero during off-duty or overnight.
- USB History Memory : Sterilization history can be downloaded to your USB memory stick and reprint anytime you want to.
- 6. Login Function : Providing User ID and PW, so that hospital can manage users who access the sterilizer.



### MEDICALPARK CO., LTD.

**Address** #601/623/624/901, Knowledge-Industry Center,  
Bundang Suji U-Tower, 767, Sinsu-ro, Suji-gu,  
Yongin-si, Gyeonggi-do, Korea

**Tel** +82-70-4800-5660

**Fax** +82-70-8266-1108

**E-mail** info@medicalpark.co.kr

**Website** http://bexcore.com

#### Technology and Product(s)



Bexcore is to inspect or remove a tissue without leaving any incision or scar. It is a method to make a 4~5mm incision and aspirate tissue into a needle by placing it near to a lesion and then the tissue cut by an inner blade of the needle is discharged automatically by vacuum suction. By real-time ultrasound, an accurate diagnosis and treatment are possible.



## ■ MEDONICA

**Address** 37-9, Techno 1-ro, Yuseong-Gu, Daejeon, Korea

**Tel** +82-2-808-0234

**Fax** +82-2-808-0235

**E-mail** info@medonica.com

**Website** www.medonica.com

### ○ Technology and Product(s)



#### Plasma sterilizer ZEUS-50 & 130S

Three kinds of cycles can be selected depending on the application and the load volume

Short cycle : 25minutes ( $\pm 5$ ) : ZEUS-50 & 130S

Standard cycle : 35 minutes ( $\pm 5$ ) : ZEUS-50 / 45 minutes ( $\pm 5$ ) : ZEUS-130S,

Advanced cycle : 45 minutes ( $\pm 5$ ) : ZEUS-50 / 55 minutes ( $\pm 5$ ) : ZEUS-130S,

Installable in a small area

Operates with economy and efficiency even with a small load : however, a large load can also be handled

The best low temperature sterilizer to start with a reasonable price



## ■ MEINNTech CO., LTD.

**Address** 5th Floor, A-dong, Digital Empire B/D, 387 Simindaero, Dongan-gu, Anyang-si, Gyeonggi-do, Korea 14057

**Tel** +82-31-381-7076

**Fax** +82-31-381-7053

**E-mail** overseas@meinntech.com

**Website** <http://www.meinntech.com/eng/main>

### ○ Technology and Product(s)



ANYFUSION is the product combining the syringe pump and infusion pump into one unit - Cylinder pump that solved problems of accidental bolus, free flow, inadvertent bolus, siphoning phenomenon etc. that lasted for over 60 years due to the technical limitation by using rotary motion of pistons on the contrary to direct motion; and capable of medicine injection such as contrast medium with continuous accuracy and precise volume control. Moreover, the cylinder cartridge maintains high accuracy within  $\pm 2\%$  deviation for long period without replacing the cartridge for 3days, thus reduces unnecessary in use that distress for users; and increases optimization of and convenience of application in use as the device is not affected by gravity nor its location with advantage of own mechanism.

- High accuracy within  $\pm 2\%$  of deviation
- Injection capable regardless of location
- Directly connected to the container without frequent replacement of consumables
- Fundamental prevention of free flow, accidental bolus, inadvertent bolus and siphoning phenomenon
- Minimum alarm, removing stress on over/under dose
- Integrated syringe & infusion mode use



**nanomtech Co., Ltd.**

**Address** 57, Cheomdan Venture So-ro, Buk-gu, Gwangju, 61003, South Korea  
**Tel** +82-62-955-8588  
**Fax** +82-62-955-8589  
**E-mail** altoall@altochem.com  
**Website** www.nanoomtech.co.kr

**Technology and Product(s)**



Heart+ResQNT-381 is a compact, portable device that analyzes the heart's rhythm and delivers a defibrillation shock to re-establish the normal heart rhythm. Heart+ResQNT-381 guides the user through each step of the defibrillation process by providing voice/ visual prompts. Heart+ResQNT-381 of Nanoomtech Co., Ltd provides electrode pads that can be used for both adults and children.  
 -5 year-warranty for main unit, 2year-warranty for electrode pads.



**Renosem Co.,Ltd.**

**Address** Unit 103-806, 397, Seokcheon-ro, Bucheon-si, Gyeonggi-do, Korea  
**Tel** +82-32-222-0180  
**Fax** +82-32-222-0110  
**E-mail** sale@renosem.com  
**Website** www.renosem.com

**Technology and Product(s)**



RENO PLASMA STERILIZER is the Best Solution for Infection Control.  
 Proven Strong Penetrability  
 Fast Sterilization Cycle & Gentle Sterilization  
 Non Toxic Residue  
 User-friendly Interface & Low running Cost






## SAEUM MEDITEC

**Address** 102-1406, 36, Bucheon-ro 198beon-gil, Bucheon-si, Gyeonggi-do, 14557, Republic of Korea

**Tel** +82-32-661-0804

**Fax** +82-32-661-0805

**E-mail** admin@saeum.net

**Website** WWW.SAEUMMEDITEC.COM

### Technology and Product(s)



#### EPI-Detection® is Medical pressure monitor to detect the epidural space

The EPI-Detection® is Medical pressure monitor and use it to detect the epidural space for accuracy of procedure of "Epidural block"

When the Epidural needle reaches the epidural space, it detects the pressure change of the muscle layer and epidural space. the pressure sensor recognizes it, and then displayed on LED so that the user can identify it.

It has a structural feature that allows the Racz's catheter to be inserted right away.

To put it simply, the device supports the safe treatment of patient by allowing the Epidural needle to precisely reach the epidural space.

This product was developed in collaboration with Severance Hospital, and so far clinical research has been actively carried out. According to the study 100 cases (EPI-Detection ver.1) in 2016, epi- detection was approved its effectiveness and stability. With this result, In July, 2017, we obtained a KOTRAN FDA approval for a medical pressure gauge.




## SAMSUNG MEDICAL ENG. CO., LTD.

**Address** 12, OKGUCHEONSEO-RO 237BEON-GIL, SIHEUNG-SI, GYEONGGI-DO, KOREA

**Tel** +82-31-499-0679

**Fax** +82-31-499-0675

**E-mail** s040431@chol.com

**Website** www.ssmed.net

### Technology and Product(s)



#### Wall suction

SME Wall suction is collecting a flood, spill, phlegm and flood phlegm from the patient body after surgery, there are 2 regulated vacuum range, 200mmHg and 700mmHg. It contains unbreakable polycarbonate 1200cc collection bottle, over-flow protector, hose and adaptor, Has CE marked.



#### Flowmeter with humidifier

SME Flowmeter with humidifier connects a gas outlet and supplies an oxygen or air to the patient who is needed it in the ward, ICU, CCU, OR and ER. The flow rate is 0~15LPM. It contains unbreakable polycarbonate bottle, audible alarm and adaptor. Has CE marked.



#### Gas outlet

SME Gas outlets are installed on the wall, bed head unit and pendant where needs to delivery a gas to the patient. All of gas outlet obtained the CE mark.



#### Bed head unit

SME bed head unit is a universal type, and shall be contained power socket, medical gas outlet, nurse call as buyer's requirement. It is customized products.



**Samyang Biopharmaceuticals Corporation**


**Address** 295 Pangyo-ro, Bundang-gu, Seongnam-si, Gyeonggi-do 13488, Korea  
**Tel** +82-2-2157-9817  
**Fax** +82-2-2157-9062  
**E-mail** www.samyangbiopharm.com  
**Website** sunwoo.kim@samyang.com

**Technology and Product(s)**




**Synthetic Suture Material (Bulk Thread)**

- Absorbable Braided Suture : TRISORB, NEOSORB
- Absorbable Mono Filament Suture : MONOSORB, MONOFAST



**Surgical Mesh for Hernia Repair**

- Proflex Mesh : Partially Absorbable PP+PGCL mesh



**Hemostat**

- SurgiGuard Oxidized Regenerated Cellulose (ORC) Hemostat



**Adhesion Barrier**


- InterGuard for Adhesion Barrier



**ScienceMedic Co., Ltd.**

**Address** #5F ScienceBuilding 104, 24beon-gil, Jeonpa-ro, Anyang City, Gyeonggi-do, Korea  
**Tel** +82-31-479-1194  
**Fax** +82-31-479-1194  
**E-mail** sgjun@sciencemedic.com  
**Website** www.sciencemedic.com

**Technology and Product(s)**



**Skin Closure**

- No scar for Staple or Suture
- Stronger than sutures in protecting incision from distraction forces
- No staple or suture tissue punctures that may allow bacteria into wound site
- No need to visit Hospital or clinic for removal(Remove patient home)
- Faster than sutures
- Fewer wound complications



**Battery Powered Handpiece**

Battery-typed power Handpiece is portable, so that it is suitable for delicate surgical operating.



**Saw Blade**

Saw Blades have a sharp edge which are used for cutting bone by vibrating horizontally. There are Linvatec, Stryker, Aesculap, Synthes, Desoutter and etc. and we could manufacture all spec that customers want. It is very competitive for the price and we are supplying over 60% of local market.



**Pneumatic Powered Handpiece**

Power Handpiece is used for drilling, reaming and pinning a bone and there are drill, reamer, wire, pin and saw attachments for the Multi Handpiece can be assembled & disassembled easily. It is ergonomic design, strong power, durability and also cheap in price.



### Sejong Medical

**Address** 11 Sinchon 2-ro, Paju-si, Gyeonggi-do, 10880 South Korea  
**Tel** +82-31-945-8191  
**Fax** +82-31-945-8190  
**E-mail** info@sejongmedical.com  
**Website** www.sejongmedical.com

### Technology and Product(s)



Laparoscopic Trocar



Laparoscopic Specimen Retrieval Bag



Single Port for Minimal Invasive Surgery



Laparoscopic Suction & Irrigation



### SM ENG. CO., LTD.

**Address** SM ENG. CO., LTD.  
**Tel** +82-51-305-8016  
**Fax** +82-51-305-8021  
**E-mail** sales@sm-eng.net  
**Website** www.sm-eng.net

### Technology and Product(s)

Suture needles made of SUS 304 stainless steel are more flexible and has high penetration performance. With our high technology, we also guarantee solidity of needle and thread attachment. We offer complete range of absorbable and non-absorbable suture.





**StraTek Co.,Ltd.**

**Address** RM705~707, 721~724, Yucheon Factopia, 30, Jeonpa-ro, Manan-gu, Anyang-si, Gyeonggi-do, Rep. of Korea 14086  
**Tel** +82-31-467-3137~8  
**Fax** +82-31-467-3139  
**E-mail** www.stratek.co.kr  
**Website** stratek1998@gmail.com

**Technology and Product(s)**



**Model: SPMI-330**  
: Combination Electrical Therapy (ICT+SSP)



**Model: STS-1000**  
Muscle Tap Apparatus (for massage)



**Model: CPS-2000**  
: Knee C.P.M



**Model: MC-1**  
: Arm & Leg Trainer active or passive motion for the rehabilitation



**Model: 3D-CENTOUR**  
: 3D Spinal Stabilization Device Computer System



**STL-101**  
: Digital Therapy Laser (for pain relief)



**SUNGKWANG MEDITECH CO.,LTD**

**Address** 1014 - 20, Bukhansan-ro, Jangheung-myeon, Yangju-si, Gyeonggi-do, Korea  
**Tel** +82-31-855-0627  
**Fax** +82-31-855-7549  
**E-mail** skmeditech@skmeditech.com  
**Website** www.skmeditech.com

**Technology and Product(s)**



## SurgiCore Co., Ltd.

**Address** Bldg. A, 43-10, Jingoksandan 5beon-Ro,  
Gwangsan-Gu, Gwangju 62464 Korea

**Tel** +82-62-714-3850

**Fax** +82-62-714-3861

**E-mail** www.surgicore.co.kr

**Website** www.surgicore.co.kr



### Technology and Product(s)



#### LigaForce

[Specification]

Shaft diameter: 10mm, Clip: M/L, 20 implantable clips preloaded.

[Features]

- 340mm shaft length
- The safety interlock
- Anti clip drop system
- Ergonomic grip design



#### SurgiTractor

[Available sizes]

XX-Small, X-Small, Small, Medium, Large, X-Large, XX-Large

[Features]

- Providing pure circle-Maximum opening with minimum incision
- Easy to fold
- 12mm port equipped high transparency silicone cover is available
- Easy to attach and detach the cover to the retractor



#### SurgiPort

[Available sizes]

5.5mm, 11.5mm

[Features]

- Breakage free
- Ergonomic design of the obturator head
- Funnel type cannula entry
- Patented screwing structure
- Easy to screw the device into the thoracic space

## TaeWoong Medical

**Address** Bldg.C 4F, Daebang Triplaon Business Tower,  
158 Haneulmaeul-ro, Ilsandong-Gu, Goyang-Si,  
Gyeonggi-Do, Korea 10355

**Tel** +82-31-904-6196 (ext.622)


**Fax** +82-31-904-6157

**E-mail** contact@stent.net

**Website** www.taewoongmedical.com




### Technology and Product(s)



#### TaeWoong Niti-S Biliary Stents

- Consists of the implantable metallic stent and the introducer system.
- Both Covered and Uncovered stents are available.
- The D-Type Stent is an exclusive design which maintains full luminal patency in the tortuous biliary anatomy for optimal drainage.



#### TaeWoong Niti-S Esophageal Stents

- Consists of the implantable metallic stent and the introducer system.
- Different types and wide ranges of stents provide various options to clinicians for different indications.
- The Double-Type Stent, Proximal release delivery system, and the TTS delivery system are our exclusive designs.



#### TaeWoong Niti-S Duodenal / Colonic Stent

- Consists of the implantable metallic stent and the introducer system.
- Both Covered and Uncovered stents are available.
- Palliative treatment for malignant Pyloric / Duodenal and colorectal obstructions are available.



## VOTEM CO., LTD

**Address** 27, Geodudanji 1-gil, Dongnae-myeon,  
Chuncheon-Si, Gangwon-Do, Korea 24398

**Tel** +82-33-910-0701

**Fax** +82-33-911-0701

**E-mail** votem@votem.kr

**Website** www.votem.kr

### Technology and Product(s)



#### VP-700 / VP-1000 / VP-1200

Display : 7"/10.4"/12.1 TFT LCD(800X600), 10 wave trace on/pff, Full ECG 7 wave display

Parameters : ECG, SpO2, NIBP, 4IBP, Resp., 4TEMP, EtCO2, ICO, EKG, Multi-gas, BIS

Interface : RS-232 port, VGA port, ECG output(opt.) port, SD memory card, 12 Languages

Feature : 7days trend data, 6 kinds of virtual screens

#### VM15, Modular Patient Monitor

Display : 15.1" TFT LCD

4.3" Transport Monitor, Easy S/W upgrade using USB (16GB)



## Wonbiongen Co., Ltd.

**Address** 56, 1 Gongdan-ro 5-gil, Gumi-si,  
Gyeongsangbuk-do, Korea

**Tel** +82-54-462-7965

**Fax** +82-54-471-7963

**E-mail** wonbiogen@wonbiogen.co.kr

**Website** http://www.wonbiogen.co.kr

### Technology and Product(s)




#### Therasorb algiplus Hydrophilic/ Adhesive dressing

Alginate ingredient, when mixed( interacted) with blood in the wound, forms alginate-blood gel and exchange ions with exudation in the wound, thereby producing an effect of stopping bleeding.



#### Therasorb AD


Therasorb AD is coated with silicone adhesive which is a lot of holes on top of polyurethane foam dressing, allowing for self-driven adhesion. Thus secondary dressing is not needed for fixation.



#### Therasorb Net

Porous film on wound contact layer does not leave residue on wound site during dressing replacement, and does not cause pain during replacement.

- Prous film minimizes expansion, allowing for excellent retention of exudate, leading to optimized humidity.



#### Therasorb Silver Hydrophilic Dressing

Polyurethane foam containing silver that has antibacterial effect and hemostasis effect of alginate which included in Therasorb AlgiPlus.





**WONTECH CO., LTD.**

**Address** 3F, Pangyo Global R&D Center, 22, Daewangpangyo-ro 712beon-gil, Bundang-gu, Seongnam-si, Gyeonggi-do, Republic of Korea / 13488  
**Tel** +82-31-1670-1450  
**Fax** +82-70-7882-8658  
**E-mail** neosoul7@wtlaser.com  
**Website** www.wtlaser.com

**Technology and Product(s)**



**VEINCARE**

Varicose veins laser treatment veincare(1470nm)



**VEINCARE1940**

New generation for treating varicose veins by using 1940nm wavelength.



**GI LASER**

ESD (Endoscope submucosal dissection) Laser for colorectal cancer and gastric cancer



**THY\_S**

New generation of surgical thyroid treatment.



**YOUNG CHEMICAL CO., LTD.**

**Address** 80-93, Golden root-ro, Juchon-myeon, Gimhae-si, Gyeongsangnam-do, Korea  
**Tel** +82-55-365-1361  
**Fax** +82-55-323-9523  
**E-mail** yc@youngchemical.co.kr  
**Website** www.youngchemical.co.kr

**Technology and Product(s)**





**ZERONE CO., LTD.**

**Address** (Shinil IT UTO Bldg., Dangeong-dong) #810, LS-Ro 13, Gunpo-Si, Gyeonggi-Do, Korea

**Tel** +82-31-453-3512

**Fax** +82-31-453-5444

**E-mail** eddie@01zeus.com

**Website** www.01zeus.com

**Technology and Product(s)**



**ARGON-Z**

Argon Plasma Coagulator, AC100V~AC240V, 50/60Hz, Gas flow: 0.1l/min~9.9l/min, Class 1, Type CF.



**ZEUS VISION**

400W, 50/60Hz, 99 memories, Class 1, Type CF, Endo-Cut(Intra Cut), Bipolar T.U.R.P.



**SENS-Z**

Specially designed for E.N.T. Available cut and coagulation in single unit. 30W AC100V~AC240V, 50/60Hz, 10 memories, Class 1, Type CF.



**ULTRA-Z**

Ultrasonic Surgical Unit for Aesthetic lipo-sculpture, 90W AC110V or AC 230V, 50/60Hz, Class 1, Type BF.



**INFU-Z 1000**

AC100~240V, 50/60Hz, DC12~18V, Class 1, Type CF  
 IPX2 standard, F. range: 0.1~999.9ml/h, V. range: 0.1~9999ml,  
 Accuracy: ± 5%, Macro/Micro mode, Battery: Li-ion(OT: 10 hours, RT: 8 hours)