

# PORTABLE MINE ALNU




**40~80kV**  
**3mA**  
**0.4mm**  
FOCAL SPOT

## HIGH- RESOLUTION IMAGES

MINE ALNU is safe and powerful.  
User-friendly, intuitive UI, with light weight and operator safety.  
Compact size but good quality, MINE ALNU is powerful.

**m-SENSOR**  
**17 x 17 DETECTOR**



# Convenience

Light weight  
Safety and Clear picture quality



LED GUIDE BEAM



ADDITIONAL FILTER FOR CHILDREN



SAFETY MODE



WIRELESS REMOTE CONTROL

SCREEN DISPLAY



C-TYPE CHARGE

## SPECIFICATION

### X-ray Generator

TYPE	High Frequency Inverter
Power Output	120W ~ 240W
Tube Voltage	40 ~ 80 kV(Increment : 1kV)
Tube Current	3 mA (Fixed)
Exposure Time Range	0.1 sec ~ 1.3 sec (Increment : 0.1 sec)
Focal Spot Size	0.4mm

### Power Adapter

Input	AC 100-240V, 50/60Hz, 1.0A
Output	DC 12V, 2.1A(MAX)

### Battery

TYPE	Lithium-ion Battery
Battery Capacity	2500mAh
Supply Voltage	DC 25.2V