

# Vasoclip<sup>®</sup>

Safe and Reliable Ligation System




## Non-absorbable Polymer Locking Ligation System

- Radiolucent, Inert
- Non-conductive
- Biocompatibility
- Superior ligation security
- No deformation when clamped blood vessel
- No artifacts, No slipping and collapsing
- Minimum tissue damage
- Color-coded cartridge

# Vasoclip®

Safe and Reliable Ligation System



## Endo- Applier



### VasoClip Specification

ART. NO	DESCRIPTION	VESSEL SIZE	COLOR	ENDO-APPLIER ART.NO
✓ VC1601	VasoClip M	2mm - 7mm	Blue	-
✓ VC1602	VasoClip ML	3mm - 10mm	Green	VC1602A
✓ VC1603	VasoClip L	5mm - 13mm	Purple	VC1603A
✓ VC1604	VasoClip XL	7mm - 16mm	Gold	VC1604A

\* 6 clips/cartridge, 14 cartridges/box, E.O Sterile

\* Reference number and article number can be different.



### YUWONMEDITECH CO.,LTD.

Rm3~4, donghwa medical instrument complex, 147-1,  
Donghwagongdan-ro, Munmak-eup, Wonju-si, Gangwon-do, Korea  
Tel. +82-33-743-0081 / Fax. +82-33-744-0081 / E-mail. info@yuwonmeditech.com  
www.yuwonmeditech.com