





Korea
Medical Device Companies
Directory Book **2018**

KOREA Medical Device Companies Directory Book



Korea
Medical Device Companies
Directory Book **2018**

Supported by



2018



This brochure illustrates a project promoted by Korean medical device companies wanting to develop a presence in global market with support from Korean government. Inside you will find how Korean medical device companies are reliable partners for global collaboration.

Biomaterial

BioPlus Co., Ltd.	6
Hankil Tech Medical Co., Ltd	7
HANSBIOMED CORP	8
Osteonic Co., LTD	9

In Vitro Diagnostic Devices(IVD)

B-bio Co., Ltd.	12
ENERGIUM CO., LTD	13
G-Tech International	14
NGeneBio Co., Ltd.	15
Redicare Co.,Ltd.	16
SEJONG BIOMED CO., LTD.	17
Z Biotech Inc.	18

Medical and Dental Supplies

ACE Medical Co., Ltd.	20
ALTOCHEM CO.,LTD	21
Bio Implant Technology Co., Ltd	22
BNR CO LTD	23
BS Medical Co.,Ltd	24
BSL co.,Ltd	25
CHUNGWOO CO., LTD.	26
CU Medical Systems, Inc.	27
Cworld	28
Dongkoo Bio&Pharma	29
FEELTECH BIO CO.,LTD.	30
LAMEDITECH	31
MEDEXEL CO., LTD.	32
N-Finders Co.,Ltd.	33
PINE BM	34
SAESHIN PRECISION CO., LTD	35

Medical Image Processing Unit & Biomedical electronic devices

Apro Korea Inc.	38
Bio Medical Technologies Co.,Ltd.	39
BIONOTE	40
Bistos co., Ltd.	41
CHAMMED CO., LTD.	42
CHARMCARE Co., Ltd.	43
DRGEM Corporation	44
DRTECH Co., Ltd.	45
DURICO C&T, INC.	46
ForU	47
GENORAY CO., LTD.	48
Healcerion Co., Ltd.	49
HNT MEDICAL Co., Ltd.	50
IDO Co., Ltd.	51
iiSM Inc.	52
LISTEM Corporation	53
MEDSTAR Co. Ltd.	54
Mega Medical Co.,Ltd.	55
MicroDigital Co., Ltd.	56
MSLINEENG Co., Ltd.	57
N-BIOTEK	58
PIXXGEN Corporation	59
PROVIX	60
RADIANQBIO CO., LTD	61
Samsung Electronics Co., Ltd.	62
Samsung Medison Co., Ltd.	63
SEOIL PACIFIC CORP.	64
SG HealthCare Co., Ltd.	65
Signers Co, Ltd.igners Co, Ltd.	66
Stemlab, Inc.	67
Zeronitec Co., Ltd.	68

Physiotherapy & Home care Equipment

AMI GLOBAL CO., LTD	70
ATEX CO., LTD.	71
ChoisTechnology Co., Ltd.	72
Color Seven Co., Ltd.	73
CYBERMEDIC CO., LTD	74
DAEYANG MECHANICS	75
DAESUNG MAREF CO., LTD.	76
Hanmed Co., Ltd.	77
Homelec Korea Co., Ltd.	78
K1MED CO.,LTD	79
KTMED INC.	80
MAXSTAR INDUSTRIAL CO., LTD.	81
NEOFECT Co.,Ltd.	82
openM inc	83
Optowell Co.,Ltd.	84
P&S Mechanics Co., Ltd.	85
SEALTECH	86
TAEHA MECHATRONICS Co., Ltd.	87
(주)웰뷰텍	88
Won Indutstry.,Co.,Ltd.	89

Surgical & Medical Supplies

ACMEDI KOREA Co., Ltd	92
AMT Engineering Co. Ltd.	93
DAEJUMEDITECH ENGINEERING CO.,LTD	94
DAIHAN LABTECH CO., LTD.	95
DAIWAH CORP., LTD.	96
DENTIS CO.,LTD.	97
ENDOVISION CO.,LTD.	98
Everaid Co., Ltd.	99

EWHA BIOMEDICS Co., Ltd.	100
Firson Co., Ltd	101
GENOSS	102
IDS Ltd	103
IMEDICOM Co., Ltd	104
LIVSMED Inc.	105
LOWTEM Co., Ltd.	106
MEDICALPARK CO., LTD.	107
MEDONICA	108
MEINNTECH CO., LTD.	109
Nanoomtech Co., Ltd.	110
Renosem Co.,Ltd.	111
SAEUM MEDITEC	112
SAMSUNG MEDICAL ENG. CO., LTD.	113
Samyang Biopharmaceuticals Corporation	114
ScienceMedic Co., Ltd.	115
Sejong Medical	116
SM ENG. CO., LTD.	117
StraTek Co.,Ltd.	118
SUNGKWANG MEDITECH CO.,LTD	119
SurgiCore Co., Ltd.	120
TaeWoong Medical	121
VOTEM CO., LTD	122
Wonbiongen Co., Ltd.	123
WONTECH CO., LTD.	124
YOUNG CHEMICAL CO., LTD.	125
ZERONE CO., LTD.	126

KOREA Medical Device Companies Directory Book

Biomaterial

BioPlus Co., Ltd.

Hankil Tech Medical Co., Ltd

HANSBIOMED CORP

Osteonic Co., LTD



BioPlus Co., Ltd.

Address MyungJi Bldg, 3rd floor, 29 Junjoo-ro 63 gil
Gangnam-gu, Seoul, Rep. of Korea, 06230
Tel +82-2-521-1898
Fax +82-2-521-1897
E-mail info@ubioplus.com
Website www.ubioplus.com

Technology and Product(s)



HyalDew All


- Ingredient: Hyaluronic Acid 20mg/mL
- Lidocaine: 0.3% (optional)
- Packaging: 1cc * 1 syringe/box



Hankil Tech Medical Co., Ltd


Address 452-29 Pureundeulpanro Paltan-myeon,
Hwaseong-si, Gyeonggi-do, Republic of Korea
(#18532)
Tel +82--31+354-6369
Fax +82-31-354-8161
E-mail Overseas01@hankiltech.com
Website www.hankiltech.com

Technology and Product(s)



Locking Plate system

HKT Locking Plate System is the first implanting system as a outcome of 40years experience in manufacturing medical device. It acclaimed high popularity with its superior quality and anatomical plate design for very different fracture surgeries. At current, it is one of the top 3 brands in Korea including all domestic and imported brands. High quality instruments for convenience of surgeons are provided and it definitely provides distinctive difference in product experience.



Multi Cable Grip Plate System

Pre-contoured Anatomical designed Plates
Cable and Locking Screw (optional) available at the same time.
Consisting of stand-alone usage type and connecting usage type, it can be selected according to the surgical scope situation.
Cable Fixable Screw improves the ease of pulling the cable and increases the stability (Cable Sliding Prevention)



■ HANSBIOMED CORP

Address 7, Jeongui-ro 8 gil, Songpa-gu Seoul, PO 04779, Korea

Tel +82-2-466-2266 (ext. 555)

Fax +82-2-463-1554

E-mail hkbaik@hansbiomed.com

Website www.hansbiomed.com

○ Technology and Product(s)

DBM Allograft

- Demineralized Bone Matrix (Gel / Putty / Inserter / Flexible Sheet)
- Demineralized bone Matrix with Cancellous Bone (Gel / Putty type)

DFDBA / FDBA

- Demineralized Cancellous Bone (Sponge block / Filler)
- Demineralized Cortical Bone (Chip / Stick)
- Cancellous Bone (Chip)

Acellular Dermal Matrix

Breast Implant for augmentation and reconstruction

Patented Helical Molding structure lifting thread for facial lifting.



■ Osteonic Co., LTD

Address Suite 1206, Ace Techno 3, 38 Digital-ro 29-gil, Guro-gu, Seoul, Korea 08381


Tel +82-2-6902-8400

Fax +82-2-6902-8401

E-mail jasmine@osteonic.com

Website www.osteonic.com

○ Technology and Product(s)



OPTIMUS CMF SYSTEM



QUANTUM TRAUMA & EXTERMITIES SYSTEM



KINEX SPORTS MEDICINE

In Vitro Diagnostic Devices(IVD)

B-bio Co., Ltd.
ENERGIUM CO., LTD
G-Tech International
NGeneBio Co., Ltd.
Redicare Co.,Ltd.
SEJONG BIOMED CO., LTD.
Z Biotech Inc.



B-bio Co., Ltd.

Address #801, Innovalley C-dong, 253, Pangyo-ro, Bundang-gu, Seongnam-si, Gyeonggi-do, 13486, Korea
Tel +82-31-8033-7577
Fax +82-31-8018-8433
E-mail yoon930@b-bio.com
Website <http://www.b-bio.com>

Technology and Product(s)



- B-bio is only company specialized in test strips with Carbon-GDH, Xylose Free and Less Hematocrit Influence technology.
- B-bio's technology fulfills the new ISO 15197:2013 requirements that ensure more accuracy.
- All production process with machinery and equipment is customizable.

Blood Glucose Test Strip Production Line Total Solution

Consulting → Factory&Process Design → Technology Transfer → Staff Training → Optimization → Support for Certification



ENERGIUM CO., LTD

Address B114, 2, GASAN DIGITAL 1-RO, GEUMCHEON-GU, SEOUL, SOUTH KOREA
Tel +82-2-6445-4421~4423
Fax +82-2-6445-4400
E-mail sales@energium.com
Website <http://energium.com>

Technology and Product(s)




GNT

- Suitable for centralized phlebotomy system
- Mass capacity of tube loading
- Easy loading-Using original tube pallet



HENS

- Individual desk top version for phlebotomy
- Flexible installation
- Easy loading-Using original tube pallet



HENm

- Compact design for confined space
- Easy installation in any environment
- Individual equipment management display
- Convenient tube allocation system



HENc

- Rechargeable battery equipped mobile unit



G-Tech International

Address Yonghyun Industrial Complex, 29,Sandan-ro
68Beon-gil, Uijeongbu-city, Gyeonggi-do, Korea
Tel +82-31-852-4070
Fax +82-31-852-4656
E-mail colee@gtech21.net
Website www.gtech21.net


Technology and Product(s)



Weight and Height Measuring Device(BMI)

Weight and Height Measuring Device automatically measures the height and weight of a person who get on the platform of the product, and provides information on height, weight and BMI.

- Weight and Height auto measurement
- Provide BMI which is world obesity standard
- Height touch bar fold function to avoid being broken
- Display can be rotate 180 degree
- Easy to install and operate
- Interface function(RS-232C, Bluetooth)
- Option : Thermal Printer(GL-310P)



Detecting patient absence kit

Detecting patient absence kit is installed on the patient's bed to detect the patient's presence / absence and inform the manager or nurse. So, Immediate action is possible.

- Installing the weight sensor of patient's bed,
- Detecting whether the patient is present or not in bed,
- Communication options: Bluetooth, RS-232, Ethernet, wireless transmitter/receiver
- Protect the human rights of patient by not monitoring by CCTV
- Being used in general hospitals and nursing facilities



NGeneBio Co., Ltd.

Address #1010 Hanhwa Bizmetro, 1104 Digital-ro, Guro-gu, Seoul, Korea
Tel +82-2-867-9798
Fax +82-2-883-9784
E-mail yj.moon@ngenebio.com
Website http://www.ngenebio.com

Technology and Product(s)



BRCAaccuTest™

- NGS IVD kit to detect BRCA1/2 germline mutation for breast/ovarian cancer predisposition.
- CE-IVD marked and Korean Ministry of Food and Drug Safety (MFDS)-IVD approval with clinical trial data in 2017.
- BRCAaccuTest™ PLUS launched for both of somatic & germline BRCA1/2 testing (CE-IVD)



NGeneAnalySys™

- Clinical NGS data analysis SW for variant calling, annotation, QC and clinical report generation.
- CE-IVD marked and Korean MFDS-IVD approval in 2017, combined with BRCAaccuTest™.
- Also compatible with SOLIDaccuTest™&HEMEaccuTest™ (CE-IVD) providing relevant therapy information in clinical report.



HEMEaccuTest™

- Hematologic malignancy NGS kit and clinical analysis SW to cover 108 genes related to AML, ALL, lymphoma and multiple myeloma etc.
- Tested with clinical samples
- CE-IVD marked in Oct/2018



Redicare Co.,Ltd.

Address #1, Ilsanbong-ro, Youngdo-gu, Busan, South Korea
Tel +82-10-2128-1421
Fax +82-51-414-8811
E-mail neil@rdy-care.com
Website www.rdy-care.com

Technology and Product(s)



The system is intended to be used for the quantitative measurement of glucose level in fresh whole blood capillary samples collected from fingertip. The system is intended for self-testing outside the body (in vitro diagnostic use) as an aid to monitor and control your diabetes. Always consult with your healthcare professional about your diabetes management. Use Uneed-G PRO Glucose Meter only with Uneed-G PRO Test Strip.



SEJONG BIOMED CO., LTD.

Address 11, Shinchon 2-ro, Paju-si, Gyeonggi-do, 10880, Korea
Tel +82-31-945-1571
Fax +82--31-945-1570
E-mail info@sejongbiomed.com
Website www.sejongbiomed.com

Technology and Product(s)



SmartMedi™ Prime & Cardiac markers

The SmartMedi™ system comprise of SmartMedi™ Prime analyzer and cartridge for an assay. It provides reliable result to who needs quick answer for making critical decision, and helps for cost-effective treatment decision at the Point of Care(POC). Our "SEJONG BIOMED's Cross-Flow Washing Technology" on the SmartMedi™ cartridge achieve precise and reproductive data at sensitive level of clinical sample. The test result within 18 minute in the SmartMedi™ Prime analyzer offers lab quality result with confidence.



Realtime HPV HR-S Detection

Realtime HPV HR-S detection simultaneously detects 14 high-risk HPVs in a single reaction. High-risk HPVs is cause of cervical cancer. Specially, HPV16 and HPV18 accounted for more than 70% of the cause that leads to cervical cancer. Realtime HPV HR-S detection has been specifically designed to identify HPV16 and HPV18, and concurrently detects the rest of high risk types (HPV 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68) contributing to cervical cancer. It provides reliable result by internal control which checks the entire process from nucleic acid extraction to PCR.



••• Z Biotech Inc.

Address #472, E9 bldg., 1, Chungdae-ro, Seowon-gu, Cheongju-si, Chungcheongbuk-do, Korea 28644
Tel +82-43-278-0075
Fax +82-43-278-0074
E-mail jhson@zetbio.com
Website www.zetbio.com

• Technology and Product(s)



Z Biotech Inc. tries to help patient detect their disease as easy as possible. With our products, AnyLab, You can check your family or pet's health only with small amount of patient sample(10~200ul of blood, feces, eye mucus) in short time(3~15mins)



For human health, Now we could provide 18 kinds of parameter for detecting human's Inflammation, Tumor, Hormone, Thyroid, Cardiac disease. We could also provide 9 kinds of RDT(Rapid Diagnostic Test) products and 4 kinds of qualitative products for veterinary health.

Medical and Dental Supplies

ACE Medical Co., Ltd.
ALTOCHEM CO.,LTD
ALTOCHEM CO.,LTD
Bio Implant Technology Co., Ltd
BNR CO LTD
BS Medical Co.,Ltd
BSL co.,Ltd
CHUNGWOO CO., LTD.
CU Medical Systems, Inc.
Cworld
Dongkoo Bio&Pharma
FEELTECH BIO CO.,LTD.
LAMEDITECH
MEDEXEL CO., LTD.
N-Finders Co.,Ltd.
PINE BM
SAESHIN PRECISION CO., LTD

ACEMEDICAL



ACE Medical Co., Ltd.

Address ACE MEDICAL CORPORATION, LTD. 3F, ACE Ville B/D 25-1, Bomoon-Dong, 7Ga, Sungbuk-Ku, Seoul 02873, Korea
Tel +82-2-920-0055
Fax +82-2-922-9319
E-mail jhlee3@ace-medical.com
 jwkim@ace-medical.com
Website www.ace-medical.com

Technology and Product(s)



PCA Infusion Pump (Model : AutoMed)

- The AutoMed Pain Control System is a valuable solution to manage pain.
- The AutoMed Pump electronically delivers a constant, controlled infusion of local anesthetic directly into the operative site for a relatively pain-free recovery.



Disposable Silicone Balloon Infuser (Model : AutoFuser)

Continuous and/ or intermittent infusion of medication for general infusion use including antibiotic, chemotherapy and pain management therapies.



Blood/Fluid Warming System

- AutoMer II warming system is designed to warm blood, blood products and liquids.
- Supplies infusion solutions or blood at 30,000ml/hr by pressure infuser.
 - The hot plate is warmed from 37°C to 41°C, Blood and Fluid Warming System makes ideal infusion temperature.

ALTOCHEM CO., LTD.



ALTOCHEM CO.,LTD

Address 6, Hwasan-ro 101Beon-gil, Idong-myeon, Cheoin-gu, Yongin-si, Kyeonggi-do
Tel +82-31-336-9347
Fax +82-31-336-9349
E-mail altoall@altochem.com
Website www.altochem.com

Technology and Product(s)



ALTOCAST and ALTO SPLINT are used to fix for fracture of bone. The products have many advantages such as superb durability, light weight and easy molding, high radiolucency and ventilation comparing with other company products.



We have been exporting casting tape and splint manufacturing machines which are fully automated machines. The production capacity is maximum 300 pcs per hour in case of casting tape machine.



Bio Implant Technology Co., Ltd

Address #2101 U-Tower, 1029 Yeongdeok-dong, Gilheung-gu, Yongin-si, Gyeonggi-do, Korea, 16950
Tel +82-31-704-2882
Fax +82-31-704-1881
E-mail sallylee@bioimtech.com
Website www.bioimtech.com

Technology and Product(s)



Kaolin Hemostatic Dressing Gauze

It is high-concentration kaolin-based topical hemostatic dressing product. Kaolin is a naturally occurring inorganic mineral that promotes the body's natural blood coagulation process, resulting in a rapid hemostatic effect. The kaolin ingredients that promotes the hemostatic process is a natural inorganic mineral that does not contain any harmful substances and is clinically a safe substance without fever/ allergy/ adverse reaction/ hypersensitivity reaction.



Collagen Bone Substitute Series

DualPor Collagen scaffold is porous calcium phosphate of dual pore collagen coated that acts as the extracellular matrix (ECM). It is consists of 60% Hydroxyapatite / 40% β -Tricalcium phosphate. DualPor Collagen is a interconnected macro / micro dual pores. Collagen is a material that accounted for about 22% of the bone tissue and helps to ossification osteoprogenitor cell, as well as the role ECM. It allows the blood vessels growth and movement of bone cells by macro pore. It is favorable to blood vessel growth, move and adherence of bone cells by the structural property with existing macro / micro dual pore, this further improves cell adhesion and proliferation by the collagen component of the surface.



BNR CO LTD

Address (Seoul Office) 2309, IS BIZ Tower, 26, Yangpyeong-ro, 21-gil, Yeongdeungpo-gu, Seoul, South Korea (07207)
Tel +82-2-6968-2461
Fax +82-2-6968-2453
E-mail export@vaxcon.com
Website www.vaxcon.com

Technology and Product(s)



VAXCON SAFE PASS Trocar(Blade trocar)

- Transparent sleeve
- Free size sealing system: It enables users to use 4.5~12mm instruments conveniently, without an additional converter.
- Sleeve separation system: It enables a convenient and useful use as the sleeve and main body can be easily separated.



VAXCON Bladeless trocar

- One touch Release button
- Scope insertion
- Bladeless tip to provide higher safety compared with steel blade.



Safe Pouch Bag

Pouch Hold-down device: It enables to manipulate the system with one hand and control the entrance of the bag suitable for a shape of specimen



Safe Wound Retractor

- Single port interface
- Minimal Incision
- Self fixation



BS Medical Co.,Ltd

Address 320, Hyoja-ro, Wansan-gu, Jeonju-si, Jeollabuk-do, Korea.
Tel +82-63-908-9933
Fax +82-63-908-9932
E-mail info@info-bsm.com
Website www.info-bsm.com

Technology and Product(s)



Automatic Assembly line

- Vacuum Blood Collection Tube/Needle Assembly Line
- Syringe Assembly Line (Auto feeding/ Auto Packing)
- Needle Assembly Line • Pen needle Assembly Line
- Dental Needle Assembly Line
- I.V Catheter Assembly Line • Injection Machine & Mold



YES PRP Kit

The extraction method is simple and unique, avoiding air and/or device mediated contamination. Time saving with simple extraction process. You can watch the movement of Buffy-coat which is floating inside plasma with the naked eye.



Prime PDO threads with needle

- Polydioxanone thread
- Cannula type PDO thread



HyaFilia HA filler

- BDDE Crosslinked HA
- Non-animal hyaluronic acid
- Proven safety and efficacy



BSL co.,Ltd

Address 6-13, Chilsan-ro 237beon-gil, Gimhae-si, Gyeongsangnam-do, Republic of Korea
Tel +82-55-328-9235
Fax +82-55-328-9236
E-mail bsl@biosl.com
Website http://biosl.com


Technology and Product(s)



ACS (Automatic Cell Station) is composed of an exclusive centrifugal separator, incubator and unit for automatic stem cell isolation. It offers a solution for automatic stem cell isolation improving isolation efficiency and minimizing human errors.

An operating program for ACS is installed, so automatic stem cell isolation is possible by simply pressing a button. As an incubator optimized for ACS, it consists of drawer-type constant temperature tank and tilting table.

ACPU is an automatic stem cell isolation unit that uses the centrifugal force of a tube shaped container and elastic force of springs. It automatically operates based on the intensity of centrifugal force.



Adinizer are the filtering discs for fat tissue micronizing and homogenizing (A Adinizer Discs: 2400 micons, 1200microns, 600microns).

Adicartrige is the housing adaptor to contain Adinizer disc. When fat tissues are penetrated through Adinizer discs, fat tissue can be very easily and efficiently micronized and homogenized without any emulsification. 5times of penetrations are good enough for adinizing process.



CHUNGWOO CO., LTD.

Address #614, 2 Gasandigital 1ro, Geumcheon-gu, Seoul, Korea
Tel +82-2-2027-2200
Fax +82-2-2027-2207
E-mail Sale2@mycw.co.kr
Website www.mycw.co.kr

Technology and Product(s)



Apolex Tite : Smart-RF Lipolysis Device

- Safe treatment possible due to the real-time temperature sensor
- Due to the RF energy able to melt the adipocyte and regenerate collagen and elastin in same time.
- Skin Lifting & Tightening • RF Lipolysis
- Double Chin • Six-Pack
- Gynecomastia • Vaginal Tightening



Ellisys : Fractional RF Micro Needling device

- By wound healing effect, able to improve the root cause of the skin trouble.
- By RF energy, helps the collagen regeneration.
- Skin Lifting & Tightening • Acne
- Scar • Pore
- Syringoma



Contlex : High Intensity Focused Ultrasound

- Deliver HIFU energy to SMAS layer which is the root cause of the wrinkle.
- Able to deliver HIFU energy to the fixed depth.
- Wrinkle Reduction • Eyebrow Lifting
- Malar Elasticity • Lifting & Tightening
- Collagen Regeneration



CU Medical Systems, Inc.

Address 130-1, Dongwhagongdan-ro, Munmak-eup, Wonju-si, Gangwon-do, S.Korea
Tel +82-31-421-9700
Fax +82-31-421-9911
E-mail angela@cu911.com
Website www.cu911.com

Technology and Product(s)



i-PAD CU-SP1 [Intelligent public access defibrillator]

- Ambient noise detection [Auto volume adjusting]
- CPR guide with voice prompt
- CPR detection
- Pads and battery status indicator
- Multi events recording




LiFEGAIN CU-HD1 [Multifunction Defibrillator / Monitor]

- Efficient e~cube biphasic technology [BTE Type]
- Manual and AED operation
- Defibrillation using paddles and disposable pads
- ECG Monitoring [3-Lead ECG/5-Lead ECG/10-Lead ECG]
- SpO2 pulse oximetry with alarms [NELLCOR]
- Noninvasive pacing mode
- NIBP [Non-Invasive Blood Pressure]
- EtCO2 [End-Tidal CO2]



PLAZE HPS-30 [Low Temperature Plasma Sterilizer]

- Microwave Plasma Sterilization System
- Compact Size
- Securing advanced safety of sterilizing object
- Precise sterilizing process under ultra low vacuum condition allows deep penetration of sterilizing factors, hence provides higher level of sterilization
- High durability with trustful basic functions



Cworld

Address 20, Neungan-ro231beon-gil,Jori-eup, Paju-si, Gyeonggi-do, Republic of Korea.
Tel +82-31-942-1617
Fax +82-31-941-1617
E-mail cworld@choiworld.com
Website www.choiworld.com

Technology and Product(s)



Leadfree apron that Cworld developed is 20% lighter than traditional lead apron in same protection rate.

We have 5 colors with darkblue, blue, olivegreen, pink, orange and 4 types. Customer maded also.

Proven radiation protection
 Out and innerfabrics are waterproof.
 Inner net prevents apron from adhering to the body due to perspiration.

Quality you can trust,
 Identification pocket on apron,
 Cost effective.
 Good design and color



Dongkoo Bio&Pharma

Address B-14F, 114, Beobwon-ro, Songpa-gu, Seoul, Korea
Tel +82-2-2684-5421 (EX 387)
Fax +82-2-2614-0648
E-mail dkbjyp@hotmail.com
Website www.dongkoo.com

Technology and Product(s)



“SmartX” Strong Point

- Convenience : Kit to Syringe method is applied to SmartX for reducing extraction time and easier operation
- Cost-effective : SmartX is a SVF extraction kit that is cost-effective in initial and operational cost
- Safety : Kit to Syringe method is adopted to SmartX in order to minimize contamination and protect from air-borne pollution
- Efficiency : SmartX extracts greater amount of SVF cells without damage or loss in a short time

“SmartX” can be used for ~

- Breast Augmentation
- Parry-Romberg Disease
- Dark Circle
- Hair Growth
- Diabetic Foot Ulcer
- Peripheral Vascular Disease
- Spinal Cord Injury
- Osteoarthritis
- Rheumatoid arthritis
- Crohn's Disease
- Parkinson's Disease



FEELTECH BIO CO.,LTD.

Address 1079-20, Charyeonggogae-ro, Gwangdeok-myeon, Dongnam-gu, Cheonan-si, Chungcheongnam-do, Korea
Tel +82-31-628-8144
Fax +82-31-707-0990
E-mail felix@dbneedle.com
Website www.feeltech.co.kr

Technology and Product(s)



Insulin Syringes

Feel-Ject

Features

- Latex Free, Non-Toxic, Non-Pyrogenic
- Nearly painless injection with comfort
- Super-secure attachment of needle
- Reduced dead space for less waste
- Bold, clear and large numbers for easy readability



Insulin Pen Needles

Feel-fine

Features

- Latex Free, Non-Toxic, Non-Pyrogenic
- Nearly painless injection with comfort
- Super-secure attachment of needle
- E.O Gas sterilization



I.V. Catheter

MF CATH PRO

Features

- Non-toxic, Non-pyrogenic
- Comfort insertion with extra needle cut
- Flexible and resilient polyurethane tube
- Radiopaque stripes
- Blood-stopping aid system



LAMEDITECH

Address Rm1002, Ace High-End 6th, Beotkkot-ro 234, Geumcheon-gu, Seoul, Korea, 08513
Tel +82-2-852-7980
Fax +82-2-852-7983
E-mail sales@lameditech.com
Website www.lameditech.com

Technology and Product(s)



HandyRay LMT-3000

The HandyRay is an innovative laser lancing device for capillary blood collection based on laser technology, that is, without normal needle. This product leads to improve the disadvantages of the traditional needle lancet including pain, callus and cross infection and also enables safer and more comfortable capillary blood sampling for the various blood tests

The HandyRay can be used with glucose level test, POCT and so on

- Laser medium : Er:YAG
- Wavelength : 2940nm
- Pulse durations : 250us
- Laser class : Class IV Laser
- Laser Power Range : 150~270mJ(±20%)
- Laser level range : 1~5 step
- Repetition : Single shot operation
- Spot size : 350um
- Weight : 170g
- Dimension : 34.7mm(W) x 99.5mm(D) x 131.5mm(H)
- Charging Voltage : Max 3.7Vdc
- Battery : ● Rechargeable Li-ion Battery ● Life cycle: 500 cycles(1 cycle: 70 shoots)
- Resting time : 2 min (After 10 times continuous laser irradiation.)

MedExel



MEDEXEL CO., LTD.

Address #252, Geumgwangosan-ro, Geumgwang-myeon, Anseong-si, Gyeonggi-do, KOREA
Tel +82-31-299-1077
Fax +82-31-677-8087
E-mail eunse@medexel.co.kr
Website www.medexel.co.kr

Technology and Product(s)



Autokeeper(Safety Pen Needle)

The Autokeeper is a safety pen needle. Safety pen needle is an insulin pen needle which automatically locked itself after the injection. The purpose of this product is to guard patients from re-using the needle and to prevent the users from secondary infection which is able to occur in the process of disposal. As MedExel keeps develop itself, it has its own design and technology for the safety pen needle. So it owns domestic patent as well.



TopFine(Insulin Pen Needle)

Insulin disposable pen needle is used in conjunction with pen type injector to deliver insulin into the body. It is for diabetic patients to inject safer and more convenient. As the needle is single use only, it prevents secondary infection of users. Also various lengths and thickness of the needle are able to be customized to patient's body condition.

.FINDERS
N-Finders Co., Ltd.



N-Finders Co.,Ltd.

Address #1211, 55, Digital-ro 33-gil, Guro-gu, Seoul [08376]
Tel +82-2-554-5188
Fax +82-2-562-5189
E-mail nfinders@nfinders.com
Website http://www.nfinders.com/eng/

Technology and Product(s)



N-Cog

Safe and simple treatment Cog thread
 Side effect free(no thread dimpling after treatment by multi-cog)
 Biodegradation by Hydrolysis completely



N-Scaffold

Braided threads for volume bringing a good satisfaction in safety and effectivity.



Derma Spring

Spring-shaped PDO with high elasticity by special methods
 Spring Thread with excellent volume effect



Derma Spring Eye

Specially designed as Blunt type for eye
 Bleeding, swelling, dimpling and other side effects are decreased by Blunt type needle



☐ PINE BM


Address 120-11, Techno 2-ro, Yuseong-gu, Daejeon 34028
Rep of Korea
Tel +82-42-710-7312
Fax +82-42-671-9078
E-mail info@pinesutures.com
Website www.pinesutures.com

🕒 Technology and Product(s)



PINE SUTURES

PDO THREAD, PLLA THREAD, PDO FILLER, INJECTOR for wrinkle removal, skin lightening, skin tightening, face lift are proven to be reliable enough to be used in cardiac surgery. PINE SUTURES manufactures these sutures as the leading medical technology in Korea is leading the market in using these sutures for cosmetics. Vsorb, PS V Line JELLYFILL are well know brand of PINESUTURES.



Items of Exhibition

Health Beauty Products, Medical Beauty Products, Packaging



☐ SAESHIN PRECISION CO., LTD

Address #52 SEICHEON-RO, 1-GIL, DASA-EUP, DALSEONG-GUN, DAEGU 42921, REPUBLIC OF KOREA
Tel +82-53-587-2341
Fax +82-53-580-0999
E-mail traus@saeshin.com
Website www.saeshin.com

🕒 Technology and Product(s)



TRAUSUS10

TRAUSUS10 is piezo & surgical implant dual machine. This ultrasonic equipment is designed for osteoplastic Amputation, including dental implant surgery. The user can get benefits of TWO motors in one unit.



TRAUSSSG10

TRAUSSSG10 is medical device in the field of plastic surgery, neurosurgery, orthopedics, dental, maxillofacial surgery, ENT. This model includes 3 x surgical angles, 3 x saw angles and 1 shaver with multifunctional micro hand-piece and micro LDC motor.

Medical Image Processing Unit & Biomedical electronic devices

Apro Korea Inc.

Bio Medical Technologies Co.,Ltd.

BIONOTE

Bistos co., Ltd.

CHAMMED CO., LTD.

CHARMCARE Co., Ltd.

DRGEM Corporation

DRTECH Co., Ltd.

DURICO C&T, INC.

ForU

GENORAY CO., LTD.

Healcerion Co., Ltd.

HNT MEDICAL Co., Ltd.

IDO Co., Ltd.

iiSM Inc.

LISTEM Corporation

MEDSTAR Co. Ltd.

Mega Medical Co.,Ltd.

MicroDigital Co., Ltd.

MSLINEENG Co., Ltd.

N-BIOTEK

PIXXGEN Corporation

PROVIX

RADIANQBIO CO., LTD

Samsung Electronics Co., Ltd.

Samsung Medison Co., Ltd.

SEOIL PACIFIC CORP.

SG HealthCare Co., Ltd.

Signers Co, Ltd.igners Co, Ltd.

Stemlab, Inc.

Zeronitec Co., Ltd.



Apro Korea Inc.

Address 3 floor 3 sanbonro 86 beongil Gunpo
Gyeonggido 15847 Korea
Tel +82-31-477-7377
Fax +82-31-477-7977
E-mail global@apro-korea.co.kr
Website www.apro-korea.co.kr

Technology and Product(s)



CoATherm AK-F50

Radio Frequency Pain Management System
Minimal Invasive Surgery

- Pain control(Cervical/ Lumbar spinal nerve block)
- IDET(Intradiscal Electrothermal Therapy)
- Muscle Reduction



CoATherm AK-F200

Radio Frequency Thermal Ablation System
Minimal Invasive Surgery

- Liver/Lung/Kidney Cancer Ablation
- Thyroid/Breast Benign Nodule Ablation
- Myoma Mass



APRO-110

Vaginal Irrigation System(Vaginal Douche Apparatus)

- Irrigation and suction can be working at the same time to clean vagina
- Temperature control of speculum and irrigation water
- Two spray guns for spraying solution



Bio Medical Technologies Co.,Ltd.

Address 101 Megu BLDG, 21-10 dongho-ro 15-gil, Jung-gu, Seoul, KOREA
Tel +82-2-2234-0780
Fax +82-2-2234-0748
E-mail bmtech@bm-tech.co.kr
Website www.bm-tech.co.kr


Technology and Product(s)



Pulse Oximetry sensor – single use/ reusable sensors / Extension cables



BME-2002 Adult BIS compatible EEG sensor



BME-2003 Pediatric BIS compatible EEG sensor



Disposable & Reusable NIBP cuff



BIONOTE

Address 34 Samsung1ro-4gil, Hwaseong-si, Gyeonggi-do, 18449, Republic of Korea
Tel +82-31-211-0516
Fax +82-31-8003-0618
E-mail En.bionote.co.kr
Website En.bionote.co.kr

Technology and Product(s)



Vcheck(Immunofluorescent analyzer)

Vcheck is animal immunofluorescent analyzer which could test inflammation markers, pancreatitis markers, hormone markers and semi-quantitative infectious disease marker.



Anigen rapid test kit

Animal Rapid test kit for infectious disease with could be tested within 10 mins with small volume of sample. 54 Companion animal, 25 Livestock products are available.



Bistos co., Ltd.

Address 7th Fl., A Bldg., Woolim Lions Valley 5-cha, 144-3, Sangdaewon-dong, Jungwon-gu Seongnam-si, Gyeonggi-do, Korea
Tel +82-31-750-0340
Fax +82-31-750-0344
E-mail BISTOS@BISTOS.CO.KR
Website www.bistos.co.kr

Technology and Product(s)



Fetal Monitors : Dual pulsed Doppler

- BT-300 : LCD(color) / LED type
- BT-350
 - Rechargeable battery
 - Waterproof probe
 - Central monitoring system (BCM350)



Patient Monitor : 3.2 TFT color LCD display

- BT-700 : SpO2, NIBP, Temperature
- BT-750 :
 - Data saving up to 72 hours
 - Rechargeable NIMH battery



Incubator : 7" Color LCD display

- BT-500
 - Air/Skin Temp. & Humidity Control
 - O2 measurement
 - Monitoring Camera (Optional)



CHAMMED CO., LTD.

Address 13Fl. A-DONG, GUNPO IT VALLEY Bldg. 17,
GOSAN Street 148beon-gil, Gunpo-si, Gyeonggi-
do, Korea

Tel +82-31-8068-1406

Fax +82-31-8068-1406

E-mail mjpark@donga.co.kr


Website www.chammed.co.kr

Technology and Product(s)



ENT Unit XU7

- Power source: AC220V, 50/60Hz
- Power consumption: 1000~1100VA(Main body)
- Air compressor: 250W, 2-way head. Max. air pressure 6Kgf/cm²
- Suction motor: 650~680mmHg negative pressure
- Dimensions: 1123(w) * 806(d) * 908(h)mm ± 10%
- Weight: Approx. 80kg ± 10% by different type(I or II)
- Irrigator temperature: 30~40°C adjustable
- Irrigator operating time: 10~990 seconds
- Irrigation water tank storage: 4.5liter water tank, auto drainage
- Switches: 7 inch LCD touch control panel
- Color: White & dark blue



CHARMCARE Co., Ltd.

Address Rm714, 2, Gasan digital 1-ro, Geumcheon-gu,
Seoul, 08591, Korea

Tel +82-2-862-5052

Fax +82-2-862-5065

E-mail info@charmcare.com / info@h2care.com

Website www.charmcare.com
www.h2care.com

Technology and Product(s)



Wristband-type blood pressure monitor

Certified by KFDA as a medical grade blood pressure monitor.
100 times measuring battery life. Measuring SYS, DIA and pulse rate.
Smallest automatic blood pressure.



Pulse Oximeter

We offer you multi-positioning type, cost-economical type, color handheld type and bedside type pulse oximeter. All products feature are as following:

- In Low perfusion (0.05~20%)
- Nellcor compatible probe
- Available for all ages.
- Adaptor or rechargeable battery operated
- Audible and visible alarm



Patient Monitor & EtCO₂ Monitor

- Customized options with various types
- Compact and light
- 7" wide color high resolution TFT screen
- Storage of 72 hours trend and alarm information
- In case of EtCO₂, compatible with RESPIRONICS® LoFlo™ and CAPNOSTAT5™



DRGEM Corporation

Address 7F, E-B/D Gwangmyeong Techno-Park, 60 Haan-ro, Gwangmyeong-si, Gyeonggi-do, Korea [14322]
Tel +82-2-869-8566
Fax +82-2-869-8567
E-mail sales@drgem.co.kr
Website www.drgem.co.kr

Technology and Product(s)



Digital Radiography X-ray system(Full Auto Ceiling System)

The GXR-Full Auto Ceiling Series is a diagnostic digital radiography system that provides reliable high quality digital radiographic images with a reduced dose and the Full auto positioning function will help you in practical usage



Mobile DR System (TOPAZ)

Integrating high technology and over a decade of experience in conventional and digital radiography systems, designed with maximum comfort for patients and user in mind.



Portable X-ray system(JADE)

JADE is an innovative, highly versatile portable X-ray system suitable for a variety of clinical uses. JADE is one of the lightest portable X-ray systems on the market, allowing it to be used in any imaginable way including bedside, operating rooms, intensive care units and in veterinary fields.



X-ray Generator: 32kW to 82kW

Full Auto positioning
 Control of Synchronization with wall bucky stand and patient table
 Table Bucky Tracking function
 Auto Stitching with wall bucky stand and patient table



DRTECH Co., Ltd.

Address 2FL/6FL SPG Dream Building, 166, Jeongjail-ro, Bundang-gu, Seoungnam-si, Gyeonggi-do, 13558, Korea
Tel +82-31-779-7400
Fax +82-31-779-7400
E-mail drlee@drtech.co.kr
Website http://www.drtech.co.kr/eng/

Technology and Product(s)



EVS 4343 Series

- Slim cassette size(ISO4090) detector with Portable handle (Optional)
- Fast preview time (≤ 2 sec.)
- Lossless AED & Perfect X-ray sync with AWC (Adaptive exposure Window Control)



EVS 2430 Series

- Slim wired and wireless flat panel detector for general radiography, mobile/portable radiography, dental and vet system.
- High resolution imaging with smallest 76µm pixel size
- OLED display and 3-button design provides users ultimate portability.



EVS 3643

- Slim cassette wireless flat panel detector for general radiography, mobile/portable radiography, and vet system.
- Embedded with wireless battery charging system enables seamless 24 hours wireless operation for perfect portability.



RSM 2430C

- Slim cassette type digital Mammo. upgrade solution
- Exceptional image quality with 76µm pixel size(smallest in CSI)
- Image sharpness increased by 30% with DRTECH's image processing algorithm TRUEVIEW® ART



••• DURICO C&T, INC.

Address 33, Oedap 6-Gil, Sangju-Si, Gyeongsangbuk-Do
37240 Republic of Korea
Tel +82-2-525-8405
Fax +82-2-525-7461
E-mail stp@durico.co.kr
Website www.durico.co.kr

• Technology and Product(s)



- Printing Black & White Video Image on Medical Ultrasound Diagnostic System.
- Precision Printing Suitable for Thermal Sensitive Printing Method.
- Fitting Perfectly to Most Sony & Mitsubishi Video Printer.
- Certified to ISO9001 / ISO13485 / ISO14001.
- FDA registered.
- UV Technology is applied in HD & HG grade.



••• ForU

Address 2-301, Business Incubation center, 28, Halladae-gil, Heungeop-myeon, Wonju-si, Gangwon-do, Republic of Korea
Tel +82-33-763-6209
Fax +82-33-765-6209
E-mail foruhelp0731@gmail.com
Website www.for-u.co.kr

• Technology and Product(s)



Purpose

It is a wireless camera(or endoscope) system and this product is intended for the purpose of photographing skin or body cavities.

Product Features

- Wireless System
- Full-HD Real-time video
- Endoscope Changeable type
- Image capture and record the video
- Data auto saved to PC (Separate S/w)
- Automatic system control by lifting up & putting down
- Auto Focus
- Adjust screen sizes suitable for endoscope size
- Auto white balance

Function

- Because it is a wireless product, it is free to move compared to wired products.
- No wireless interruption.
- No wireless collision with other equipment.
- Easy to use with button operation and automatic system.
- The user can see clear pictures and videos in Full HD.
- The user can change the endoscope of the product according to the diagnosis area.



GENORAY CO., LTD.

Address 512, 560, Dunchon-daero, Jungwon-gu, Seongnam-si, Gyeonggi-do, 13230, Korea
Tel +82-31-627-3900
Fax +82-31-740-4120
E-mail jay@genoray.com
Website www.genoray.com

Technology and Product(s)



Digital C-Arm X-ray, OSCAR 15, has incredible 26x26cm flat panel powered by a 15kW output generator with beautiful design and high performance at a very affordable price.

- 26x26cm flat panel detector can be adapted
- 15kW HFG (Rotating monoblock)
- Intelligent ABC (Automatic Brightness Control) system
- Snap shot mode
- SID(1,000mm)&Free space(800mm), 150 Orbital rotation



Digital C-Arm X-ray, OSCAR Classic, has incredible 9" image Intensifier powered by a 5.3kW output generator with beautiful design and high performance at a very affordable price.

- 9" Image Intensifier adapted
- 5.3kW HFG (Stationary monoblock)
- Intelligent ABC (Automatic Brightness Control) system
- SID(1,000mm)&Free space(800mm), 150 Orbital rotation



New paradigm of digital mammography system, DMX-600, shows new indicator to the screening and diagnosis of breast cancer for all women in the world.

- Excellent & High resolution Image Quality
- Exposure area selection of Large F.O.V.
- User friendly Versatile functions
- Smart Software user alike



Healcerion Co., Ltd.

Address 804ho, 38-21, Digital-ro 31-gil, Guro-gu, Seoul
Tel +82-2-6347-6328
Fax +82-70-7610-6826
E-mail hilda@healcerion.com
Website www.healcerion.com

Technology and Product(s)



SONON is designed to provide primary care physicians with affordability, portability and flexibility. SONON fills a previous void for primary care physicians to improve the efficiency of these procedures while also extending their range of diagnosis. No longer does the primary care physicians need to rely on simply physical diagnosis and tools: the expanded range and increased offerings that SONON provides can also lead to direct financial benefits for primary care physicians. Additionally, SONON's ultracompact, wireless design and affordability mean the device requires minimal space and resources to operate. With its direct mobile device connectivity, SONON is an innovative offering that provides accurate and timely diagnostics, confidence and peace of mind for both primary care physicians and patients. SONON is designed for most common, everyday diagnostic and screening applications in primary care. It is appropriate for pain management and scanning applications include musculoskeletal, tissue inflammation, vascular, and general abdominal scanning.

SONON also assists in procedures like mass removals, guided injections and fluid aspirations, which can be performed by the primary care physicians to save time and money for the patient. With SONON, diagnostic and screening time is reduced, and the increased efficiency and clearer information provided to secondary and tertiary care from primary care physicians leads to a shorter overall process from early diagnosis to final diagnosis. Additionally, the overuse of conventional ultrasound devices is alleviated as primary care physicians can grow their diagnostic and imaging capabilities, providing more value for the patient and expanding the primary care physician's offerings and value provided to the patient.



HNT MEDICAL Co., Ltd.

Address 20th floor, 70, Gasandigital 2-ro, Geujeon-Gu, Seoul, Korea (08589)
Tel +82-2-6952-2588
Fax +82-2-3289-0705
E-mail hntsales@hntmedical.co.kr
Website www.hntmedical.co.kr

Technology and Product(s)



ESWT HYDRON-QT

Extracorporeal Shock Wave Therapy(ESWT)
 EH(ElectroHydraulic) Type
 Focus/ Multi Focus (Option)
 • Strong Shock wave
 • Easy to change Electrode



ESWT REGENWAVE

Extracorporeal Shock Wave Therapy(ESWT)
 EM(ElectroMagnetic) Type
 Focus/ Multi Focus (Option)
 • Low noise
 • Proved clinical safety & Low pain



ESWT Carewave

Extracorporeal Shock Wave Therapy(ESWT)
 EM(ElectroMagnetic) Type
 Focus/ Multi Focus (Option)
 • Easy to change EM transducer by an operator
 • Light & Small Volume by a compact design



IDO Co., Ltd.


Address Works A 1010-1012 Garden 5 Bldg., 52 Chungmin-ro, Songpa-gu, Seoul Korea, 05839
Tel +82-2-401-7799
Fax +82-2-400-4711
E-mail jykang@ntllab.com
Website www.ntllab.com

Technology and Product(s)



Dr.Cervicam® WiFi - A Wireless Cervical Camera

It is composed of a digital cervical camera, medical software programs, and specialized evaluators for examination.



Medical S/W Programs



iiSM Inc.

Address 402, 7-16, Naseongnam-ro, Sejong-si, 30129, Republic of Korea
Tel +82-44-863-9852
Fax +82-44-863-9853
E-mail ism4u@iism.co.kr
Website www.iism.co.kr



LISTEM Corporation


Address 94, Donghwagongdan-ro, Munmak-eup, Wonju-si, Gangwon-do, 26365 Republic of Korea
Tel +82-2-860-9596
Fax +82-2-860-9569
E-mail ex1@listem.co.kr
Website www.listem.co.kr

Technology and Product(s)

Technology and Product(s)



Vein Probe is a medical device, which can be used to visualize blood vessels using near-infrared technology. Vein Probe displays the blood vessels without any distortion. Our device can be used for infants, elderly and obese patients regardless of their skin tone. Vein Probe helps to find blood vessels promptly and precisely. Vein Probe can decrease the number of unnecessary repetitive injections. Vein Probe will increase the safety and satisfaction of patients as well as the work efficiency of medical professionals.




CLASCAN is a digital R/F X-ray taking device offering a wide range of functions and features for maximum patient convenience. Also, it is a remote-controlled system that enables the operation to be done from the main and local consoles, thereby enhancing convenience and exposing the operator less to X-ray radiation



SMART DR - Now, you can upgrade your existing analog X-ray system to fast and convenient DR. SMART DR features brilliant automatic positioning of standing mode and table mode. Small footprint of SMART DR is perfect for any hospitals/clinics environments. LISTEM's proven image quality and intuitive imaging software come standard.



REX-R features versatility and high performance which can be applicable to a wide variety of general radiography with remarkable convenience, productivity, and patient comfort. The impressive convenience for patient and technician will be maximized by using of four-way floating or six-way floating top bucky table and combining with various tube stands from floor-mounted type to ceiling-suspended type.



MEDSTAR Co. Ltd.

Address 35-12, Goraedeul-gil, Uiwang-si, Gyeonggi-do, Korea, 16071
Tel +82-32-457-1111
Fax +82-32-457-2670
E-mail president@medstar.co.kr
Website www.medstar.co.kr

Technology and Product(s)



- Compact & Elegant, Modular Integral ENT Workstation.
- Built-In Ear Washing, Automatic Drainage System, Powerful Suction,
- Spray w/ Polizer Olive Type, Controlling Press, Gauge.
- Clear & Definitive Full HD / HD Camera w/ LED, Endoscopes Holder w/ Warmer.
- One-Touch LED Control Panel
- A Wide Range of Surgical Instruments Drawers w/ Warmer, Built-in RF Machine.



- Double / Single Full HD / HD Camera with Two / One LED + Light Cable
- All in One Functioned, into 24 inch LED Monitor, Camera Switched
- Single / Double / Quadruple Images Taken & Displayed
- Live Video Saved & Replayed onto USB & Monitor, respectively
- Stand / Cart Type Chosen



Mega Medical Co.,Ltd.

Address 110-1, Donghwagongdan-ro, Munmak-eup, Wonju-si, Gangwon-do, [zip 26365] Republic of Korea
Tel +82-2-3662-4493
Fax +82-2-3661-0120
E-mail webmaster@megamedical.co.kr
Website www.megamedical.co.kr

Technology and Product(s)



NET-1100

Features

- Best compacted size
- Portable type
- High quality endurance & noiseless suction motor
- Various useful options
- Well-designed instrument cart (Option)

Specifications

- Compressor pump : 0.9 - 2.0kg/cm²
- Vacuum pump : 600 - 680mmHg
- Anti-Fog device : 400W
- Power consumption : 850W
- Voltage : 110/220V
- Frequency : 50/60Hz
- Weight : 80kg (Full set)
- Dimensions (W x D x H) : Full Setting : 1,050 x 590 x 1,250 mm
Main : 500 x 422 x 934 mm / Tray : 500 x 500 x 975 mm



NET-600A

Features

- Various useful options
- Large and sufficient storage cart
- Central Water Supply, Auto-Drainage System (Option)
- Left-handed type is available

Specifications

- Compressor pump : 0.9 - 2.0kg/cm²
- Vacuum pump : 600 - 680mmHg(Max)
- Main Suction bottle : 2,000 - 3,000cc
- Sub Suction bottle : 500 - 1,000cc
- Anti-Fog device : 400W
- Power consumption : 1,000W (Max)
- Voltage : 110/220V
- Frequency : 50/60Hz
- Weight : 130kg
- Dimensions (W x D x H) : Full Setting : 1,295 x 445 x 865 mm
- Chart table size : 435 x 445 x 865 mm



MicroDigital Co., Ltd.

Address WINS-7th floor, 15, Pangyo-ro 228beon-gil, Bundang-gu, Seongnam-si, Gyeonggi-do, 13487, Rep. of KOREA
Tel +82-70-5138-8507
Fax +82-31-702-2225
E-mail hjyoon@md-best.com
Website www.md-best.com

Technology and Product(s)



Nabi
UV/Vis Nano Spectrophotometer



Diamond
Fully Automatic Immuno-Analyzer for ELISA and CLIA



Lumi
Single Tube Luminometer



iSBS
Integrated Sample Banking System



LuBi
Microplate Luminometer



MSLINEENG Co., Ltd.

Address 26, Wonjeok-ro 148beon-gil, Sindun-myeon, Icheon-si, Gyeonggi-do, Rep. of Korea
Tel +82-32-216-5022
Fax +82-32-216-5022
E-mail yjkwon1023@mslineeng.com
Website www.mslineeng.com

Technology and Product(s)

CVP-2
THE WORLD'S FIRST ANTI-SCATTERING RADIATION FILTER FOR C-ARM

Radiation dermatitis developed on the hands of an orthopedic surgeon

ARE YOU SAFE FROM RADIATION EXPOSURE?

Performance Scatter Map Result

Without CVP (mR)	*With CVP (mR)
0.25	0.075

TRANSMITTED RADIATION LEVELS

Organ	Without CVP-2 (mSv)	With CVP-2 (mSv)	Reduction Rate
DOCTORS' EYE	0.208	0.062	60%
DOCTORS' THYROID	0.158	0.047	62%
DOCTORS' LIVER	0.185	0.055	60%

Patient Radiation Dose(Direct Beam)

Without CVP-2 (mSv/min)	With CVP-2 (mSv/min)	Reduction Rate
46.33	12.45	73%

Reference : IOPSCIENCE(Journal of Radiological Protection)
 2017 MSLINEENG. ALL RIGHTS RESERVED.



■ N-BIOTEK

Address 402-803 Techno-Park, Pyeongcheon-ro 655, Wonmi-Gu Bucheon-Si, GyeongGi-Do, KOREA < zip code: 14502>
Tel +82-32-321-2100
Fax +82-32-328-2372
E-mail export@n-biotek.com
Website www.n-biotek.com




■ PIXXGEN Corporation

Address 5F, SMART BAY, 123, Beolmal-ro, Dongan-gu, Anyang-si, Gyeonggi-do, 14056, Korea
Tel +82-70-4846-8888
Fax +82-70-6455-2905
E-mail sales@pixxgen.com
Website www.pixxgen.com




● Technology and Product(s)



We, N-BIOTEK are a main manufacturer of biomedical and lab equipment in the world. Since we have been supporting to establish personal labs by developing compact and customized products under the concept of “Handy Lab”, our products are renowned for their attractive design, unique features, outstanding quality and competitive price.

Since we have expanded our business to stem cell field, we are providing stem cell total solution which provides everything for establishing a stem processing facility. Also, we are providing training programs which instruct stem cell processing technology, stem cell activation cosmetics and isolating consumables as well.

● Technology and Product(s)



PIXXGEN introduces the first Cassette Sized, 17”x17”, UltraLight, Wired & Wireless, flatpanel detector. Our detector provides No Calibration (ACC) & No Connection to the xray system (AED). Its long battery time, the lightest weight, Cassette Size, AED, ACC and high speed data transmission make it as easy to use and install. Also provides simple HDMI type Data link cable & auto-transforms from Wireless to Wired.



PROVIX

Address 601, Bucheon Chunui Techno Park 3, 80 Jomaru-ro 385 beon-gil, wonmi-gu, Bucheon-si, Gyeonggi-do, South Korea
Tel +82-32-661-8135
Fax +82-32-661-8138
E-mail daniel@provix.co.kr
Website www.provix.co.kr

Technology and Product(s)



This Smart 4K model is designed for Various Endoscopy Applications such as Otolaryngology, Arthroscopy, Laparoscopy, Urology, Gynecology etc. As this 4K UHD Device shows the highest resolution in our product line, we are sure that the performance of this device can satisfy and meet any demands & requirements from the customers. This product will soon be integrated with the AI system.



- Android OS CCU control
- 4K UHD resolution 3840 x 2160 pixels
- HDR10 Video Processing
- FHD 1920*1080 60fps recording
- FHD 1920*1080 60fps playback
- Multi-Language OSD Support
- Capture - PIP/ Resizing support
- Moire Reduction function



RADIANTQBIO CO., LTD

Address #1804, HALLA SIGMA VALLEY, 53, Gasan Digital 2-ro, Geumcheon-gu, Seoul, 08588, Korea
Tel +82-2-852-1122
Fax +82-2-6499-3076
E-mail overseas@radianqbio.com
Website http://www.radian.co.kr

Technology and Product(s)



RADIANTQBIO AED (Automated External Defibrillator) helps to recover normal cardiac function upon sharp heart disease, such as sudden cardiac arrest, by providing cardiac function analysis and, if necessary, electroshock. HeartGuardian™ as an AED, provides precise analysis of victims' body conditions and configures exact amount of shocks needed, promising safety and quicker responses.



Easy to use and accessible by all without special training. When HeartGuardian™ is opened, auto-guidance will be provided for every stage of the operation. In case of difficulty in listening to audio-guidance, there are LED indicators as well. From analysis to charging, all services are automatic. All you need to do is to attach the pads at the right spot as indicated in visual instruction and press the button to turn on the AED. It is compatible for both Adult and Child by simply pressing adult/child mode button. HeartGuardian™ provides precise amount of shock energy tailored to patients' body size and condition of two different levels. It offers data extraction with the embedded memory. HeartGuardian™ AED offers «defibrillation black box» by recording operation history. Once recorded, data can be transferred out and analyzed through Bluetooth-enabled functions.



SAMSUNG



Samsung Electronics Co., Ltd.

Address (Maetan dong) 129, Samsung-ro, Yeongtong-gu, Suwon-si, Gyeonggi-do 16677, Korea
Tel +82-2-2194-1000
Fax +82-2-563-952
E-mail healthcare.sec@samsung.com
Website www.samsunghealthcare.com

Technology and Product(s)



GC85A

GC85A is a digital radiography system offering streamlined workflow, diagnostic confidence and high image quality with technologies such as S-Detector™ and S-Align™. The device adopts Low Dose technology for enhanced patient safety. Samsung's image processing engine S-Vue™ 3.02 provides advanced de-noising technology to deliver the same image quality as the predicate device, but with half the radiation exposure.



GM85

GM85 is a mobile digital radiography system which allows easy navigation with its compact design. 555mm narrow width and 349kg light weight allows easy access around tight spaces while the collapsible column gives users clear visibility when moving the system. Adaptive soft-driving control and front-bumper sensor also help make navigation safer. S-Align™ displays the detector's angle to the THU (Tube Head Unit) for precise alignment along with S-Vue™ 3.02 to provide high quality images.

SAMSUNG MEDISON



Samsung Medison Co., Ltd.

Address 3366, Hanseo-ro, Yangdeokwon-ri, Nam-myeon, Hongcheon-gun, Gangwon-do, Republic of Korea
Tel +82-2-2194-1000
Fax +82-2-2194-1028
E-mail mkt.sm@samsungmedison.com
Website www.samsunghealthcare.com

Technology and Product(s)



RS85

RS85 is a radiology-designated ultrasound system that provides enhanced image quality, usability, and convenience for medical and radiology professionals. Built with exquisite image quality and expert tools, such as MV-Flow™ for detecting blood flow in microvascular, S-Shearwave Imaging™ for quantifying the elasticity of tissues and S-Fusion™ that enables simultaneous localization of a lesion using real-time ultrasound in conjunction with other volumetric imaging modalities, empowers professionals to make faster and more confident decisions.




WS80A with Elite

WS80A with Elite is premium ultrasound system for OB/GYN applications, which is designed to provide superior image and exceptional productivity in workflow through highly useful features. It is equipped with a new ultrasound image processing engine called CrystalLive™, which provides enhanced structure expression. Built with innovative features, such as 5D features, CrystalVue™ and IOTA-ADNEX support healthcare professionals in making faster and more accurate decisions for women's health issues.



HM70A with Plus

The HM70A with Plus, built in a compact hardware model delivers mobility experience, helps make ultrasound exams and ultrasound guided procedures more accurate and streamlined with its image performance and efficient, easy-to-use features. The HM70A with Plus' slim and compact design offers versatile portability and increases productivity.



SEOIL PACIFIC CORP.

Address Rm.515, 516, Ace High-End Tower, 5, Digital-ro
26-gil, Guro-gu, Seoul, Korea
Tel +82-2-325-2106
Fax +82-2-325-2107
E-mail sungha@seoilpacific.co.kr
Website www.comfortcough.com

Technology and Product(s)



Comfort Cough is to help the patients who has difficulty with secretion clearance or inability to cough to clear the secretions in the airway. Comfort Cough mechanically stimulates the natural cough by shifting positive pressure to negative pressure rapidly. There are several alternatives to clear the airways, suctioning, bumping the back of the patients, but these alternatives are making pains to the patients while the result is not so effective. But Comfort Cough is using the finely controlled air pressure. It allows no pains neither no contamination. It gives perfect results which are almost same as those of natural cough. Comfort Cough has several supporting features which is the most upgraded version among available in the market. It has an Oscillation function during inhale and exhale phase. This will help to segregate the mucus from the airway wall and will enhance the clearance efficiency. It has a Percussoring function in inhale phase. This gives more powerful impact on the airway and will accelerate the clearance efficiency.



SG HealthCare Co., Ltd.

Address 10 Yangnyeong-Ro, Yangpyeong-Ri, Moga-Myeon,
Icheon-Si, Gyeonggi-Do, Korea 17408
Tel +82-2-449-6908
Fax +82-2-449-6909
E-mail sales@sghealthcare.com
Website http://sghealthcare.co.kr

Technology and Product(s)



Jumong General

For daily image examination, especially when increased patient volume pushes your existing radiography resources to the limit, this system can help. Available with a variety of options, it is capable of handling a full range of applications including chest, abdominal and extremity examinations.

Jumong F

High reliability, ease of use and diagnostic confidence. Enjoy a radiography system designed for reliable, consistent performance and low cost ownership. With a flexible range of positions, Jumong F provides your patients a fast and smooth examination.

Jumong M

Fully featured and designed for chest, supine and other radiography examinations. Whatever you need, you can find the right solution with ease.



Signers Co, Ltd.igners Co, Ltd.

Address A-208, Keumkang-Penterium IT Tower, 215, Galmachi-ro, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea, 13217
Tel +82-70-8808-7478
Fax +82-31-754-7479
E-mail libby@signers.kr, speedgun72@signers.co.kr
Website <http://www.signers.co.kr>

Technology and Product(s)



SUPiA CR System offers such a better clinic environment with high resolution image quality and affordability. It is very light and fast. Our console software 'gram' is simple, fast, and compatible. You can connect with any DICOM printer and paper printer. SUPiA gram offers easy multilingual support and special tools.

- Compact4 & Lightweight : 436 x 636 x 196mm, 21.5kg (47.4lbs)
- High Throughput : Up to 72 IPs/Hour (14 x 17" / 160 μm)
- Extremely Simple Structure Design
- Authentic Rigid Type among Clinic CRs
- Able to bear up to 150kg of weight on cassette
- Single Cassette Feed
- 14 x 17" : 2.1kg, 10 x 12" : 1.3kg
- Single Phase 50 ~ 60Hz, AC 90 ~ 264V



Stemlab, Inc.

Address #402, Suntech City Building, 474, Dunchon-daero, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea
Tel +82-70-7586-9615
Fax +82-70-7586-9696
E-mail sales@sm-eng.net
Website <http://www.stemlab.co.kr>

Technology and Product(s)



SUPRA ELITE (BMD)

SUPRA ELITE is a pencil-beam densitometer using DEXA (Dual energy X-ray Absorptiometry) based on central bone type for a higher degree of accuracy, intuitive user interface and outstanding cost efficiency. In addition, the system allows users to inspect high-risk regions including spinal column, pelvis and wrist, and It provides reliable analytical results for both clinical and research purposes.

SUPRA ELITE uses an auto-calibration algorithm to adjust counter values and maintain its best performance in addition to user-friendly menu driven and easy operation. The result sheet includes BMD data, other useful information and trend chart that may provide important assistance in clinical treatment.



DEXINO (BMD)

DEXINO is a central bone type equipment of the pencil beam type that uses dual-energy X-ray absorptiometry (DEXA), DEXINO enables a more effective analysis of results and forecasting of trends. With DEXINO, scanning particular parts, such as the spine and the femur region, would be so much easier. DEXINO is more accurate bone density measurement, which is equipped with various and effective features.



Zeronitec Co., Ltd.

Address #202-602, 388, Songnae-daero, Bucheon-si, Gyeonggi-do, Republic of Korea.

Tel +82-32-811-8890

Fax +82-32-819-8890

E-mail info@zeronitec.com, zeronitec@gmail.com

Website www.zeronitec.com

Technology and Product(s)



Product portfolio

- PURE-50 (50L) : Floorstanding type
- PURE-120 (120L) : Floorstanding type, Automatic sliding door
- PURE-120D (120L): Double door system

Features

- Economical cost & Reliable quality product (CE & ISO) : We provide our global partners with reliable & high-end quality products at more competitiveness. Our products are in compliance with CE and use high-end quality components for our sterilizer system which give our partners and customers benefits of our product reliability.
- Safe cassette-type sterilization agent : PURE sterilization agents (H₂O₂) are absolutely sealed cassette type and multiple-cycle(7cycles/ cassette) cassette, it's safe to be transported and stored and convenient and safe to be treated by user.
- Graphic User interface : Most of operating menus are displayed both graphically and numerically. User can recognize easily the status of sterilizer such as sterilant cassette, process monitoring, process parameters with less human error.
- Front large process monitoring display
Large process monitoring LED device located in sterilizer door can provide user with more convenient use of our sterilizer system. So, user can know which cycle phase is now in process and can recognize remaining cycle time easily even from a distance..

Physiotherapy & Home care Equipment

AMI GLOBAL CO., LTD

ATEX CO., LTD.

ChoiTechnology Co., Ltd.

Color Seven Co., Ltd.

CYBERMEDIC CO., LTD

DAEYANG MECHANICS

DAESUNG MAREF CO., LTD.

Hanmed Co., Ltd.

Homelec Korea Co., Ltd.

K1MED CO.,LTD

KTMED INC.

MAXSTAR INDUSTRIAL CO., LTD.

NEOFECT Co.,Ltd.

openM inc

Optowell Co.,Ltd.

P&S Mechanics Co., Ltd.

SEALTECH

TAEHA MECHATRONICS Co., Ltd.

(주)웰뷰텍

Won Industry.,Co.,Ltd.



AMI GLOBAL CO., LTD

Address 22, JANGNIMBEONYEONG-RO, 104BEON-GIL, SAHA-GU, BUSAN, 49465, SOUTH KOREA
Tel +82-51-293-6472~3
Fax +82-51-293-6474
E-mail drmedkorea@gmail.com
Website www.drmedbrace.com

Technology and Product(s)



We AMI Global Company, have been a global medical device company specialized and experienced in the field of comfortable orthopedics, rehabilitation, healthcare and sporting supports since June, 1965 with its medical line, Dr.MED and sport line. The operating system of factories is accredited by ISO13485 and its products are certified by FDA(USA) in 2003, CE mark in 2006 and GMO(medical device). Our own brand DR.MED has become various international partners for all those working at medical support/brace all over the world and will bring you improvement of person's life quality, comfort and pleasure. And then you can see more detail about our product line in the website as below.
 (www.demedbrace.com)



ATEX CO., LTD.

Address 144, Jangjisan-ro, Gwangtan-myeon, Paju-si, Gyeonggi-do, 10955, Korea
Tel +82-31-944-0546
Fax +82-31-942-6182
E-mail atexsales@atexmedical.com
Website www.AtexSports.com

Technology and Product(s)




Kinesiology Tape is designed to mimic the human skin in both thickness and elasticity to allow a wide range of mobility and flexibility without restricting your body movements. When the tape is properly applied, it will lift the skin to create increased circulation of blood, lymph, and tissue fluids which results in faster relief of muscle and joint pain. It is becoming one of the must have essential tools in today's pain management and preventive practices.



ChoisTechnology Co., Ltd.


Address #S-1404, Songdo Technopark IT center, 32, Songdogwahak-ro, Yeonsu-gu, Incheon, Korea (21984)
Tel +82-32-346-3409
Fax +82-32-246-3406
E-mail inquiry@choistec.com
Website www.choistec.com

Technology and Product(s)



Smart monitoring thermometer 'Thermosafer(XST200)':

- (via dedicated app) It has real-time temperature monitoring and abnormal body temperature alarm and it can record the temperature measurement result and check history with graph
- More stable and accurate temperature measurement is possible comparing with Non-contact type thermometer (it's error range is $\pm 0.1^{\circ}\text{C}$)
- Superior measurement distance by applying Bluetooth technology (max 15m)
- Using low irritant silicon patch for sensitive skin.



Color Seven Co., Ltd.

Address 4F, Geoam Bldg, 62, Banpo-daero 4-gil, Seocho-gu, Seoul, Korea
Tel +82-2-3462-5618
Fax +82-2-3462-5619
E-mail sales@color7.co.kr
Website www.color7.co.kr/eng

Technology and Product(s)



VPOD - New Light Therapy Patch

VPOD is a unique healthcare light therapy device which is consumable model. It is a skin adhesive Light Therapy Device based on the Low Level Light Therapy technology. This is an effective and safe electric device for women & men's health and also for beauty purpose. The principle is called PAMS technology (Color Seven's own intellectual property). Light from the device relaxes the smooth muscles in the tissue and recovers blood flow.



Woman Stress Free (COLOR DNA WSF)

Woman Stress Free is the first and unique in the world for Dysmenorrhea treatment medical device. It is a skin adhesive Color Light Therapy based on the Low Level Light Therapy technology that proved scientifically the principal of treatment. This is an effective and safe medical device for women who have painful menstruation stress due to PMS (PreMenstrual Syndrome) & dysmenorrhea. The principle of treatment is based on PAMS technology (Color Seven's own intellectual property) to relax the smooth muscles in uterus tissue. Women Stress Free was medically verified through clinical evaluation with menstrual pain patients in two University hospitals and the result has been published in Springer, Vol. 286:947-952, also this evaluation presented "Issues in Gynecology, Obstetrics, Fertility and Pregnancy Research in 2013 edition"



CYBERMEDIC CO., LTD

Address 203/1, Standard factory, 20-23, Yackchon-ro 8-gil, Iksan-si, Jeollabuk-do, Republic of Korea
Tel +82-63-835-5114
Fax +82-63-834-0154
E-mail Daniel.lee@cybermedic.co.kr
Website www.cybermedic.co.kr

Technology and Product(s)



MotoCog is Interactive ADL & Cognitive rehabilitation training system for elderly, hypofunction of ADL, damaged of cognitive function, central nerve injury and dementia patients



An instrument for analyzing and strengthening the anterior and posterior movement of the tongue to overcome dysphagia. Improved dysphagia by increasing strength and coordination of tongue, which is important for food movement also pressure making and holding, when swallowed. An instrument for measuring and strengthening laryngeal muscle resistance for overcoming dysphagia.



A System of exercise for muscle strength, muscular hypertrophy, muscle endurance, pain modulation and muscle coordination training according to measurement of maximum muscle strength using electronic pulley(Apply center of professional sport and musculoskeletal center and physio therapy exercise).



The FES for dysphagia patient who cannot easily swallow food and liquid, has a built-in treatment pattern and stimulate the muscles necessary to swallow. Exclusive pattern helps to exercise a swallowing and to reach a normal swallow.



DAEYANG MECHANICS

Address DaeYang Mechanics 32 Muhak-ro 35 Gil , Dongdaemun-gu, Seoul 130-823 South Korea
Tel +82-2-929-2222
Fax +82-2-928-4411
E-mail info_hubex@naver.com
Website www.hubex.com

Technology and Product(s)



3D Care Exercise System for the Pelvis, Back and Abdomen

A multi rehabilitation system that trains the core area with the extensor and refractive exercise on the hip joint and the lumbar joint while completely relaxing the body muscles. As the upper body maintains zero weight by the gravity offset program, the neck, shoulder and back muscles do not become tense but relaxed. The back and the knees are safely protected so the elderly and patients can easily respond to the musculoskeletal disorders. Hospital, public health, occupational health, women's fitness and silver care



DAESUNG MAREF CO., LTD.

Address 298-24, Gongdan-ro, Gunpo-si, Gyeonggi-do, Korea.
Tel +82-70-8740-6056
Fax +82-31-459-7215
E-mail medical@dsmaref.com
Website www.dsmaref.com

Technology and Product(s)



DTS-3000 (Smart Tourniquet System)

We are providing the highest functional pneumatic tourniquet systems in market. It provides suitable functions for an operating room under the optimized conditions. (Recording, IOP sensing, LCD Touch Screen, Safety Alarm System, Battery mode).



DVT-2600 (Venous Assist System)

DVT-2600 IPC system can be used individually as an alternative medicine for a stroke patient, or a patient who cannot use orally administer drugs, or any patient who cannot use anticoagulants such as heparin or warfarin, to prevent DVT and PE(Pulmonary Embolism)



Hanmed Co., Ltd.

Address 80-121, Golden Root-Ro, Juchon-Myeon, Gimhae-Si Gyeongsangnam-Do, KOREA
Tel +82-55-331-0575
Fax +82-55-331-0547
E-mail kinetrac@naver.com
Website www.kinetrac.com
 www.hanmed.net

Technology and Product(s)



Decompression System Kinetrac DAVINCI

- Traction With Twist
Traction With Twist directly focused on the Axis of Spine
- Traction With Target Lordosis a new-concept Kinetic Decompression technology have created the new theory of disc decompression.
- Traction With Muscle relaxation therapy
Pelvic joint, Lumbar Muscle Flexion/Extension.



3D Gravitation Muscle Training and evaluation System 3D NEWTON

- 3D Newton is a new conception of exercise system for rehabilitation with 3D Gravitation Training an system. It is possible of 3D space exercise both 360° rotation and 0° ~ 60° inclination control.
- Real-time Bidirectional Monitoring training system



Homelec Korea Co., Ltd.

Address A-522, Tera Tower, Songpa-daero 167, Songpa-gu, Seoul, Republic of Korea
Tel +82-2-521-2288
Fax +82-2-2054-3588
E-mail info@homelec.co.kr
Website www.homelec.co.kr

Technology and Product(s)



CM07 : The Best Premium Cordless Massager

- Dual motor : powerful percussion motor with 3,700RPM, gentle micro-vibration motor with 11,000RPM
- Adjustable massage strength with dial button
- Patented air cushion stick for joints and bones
- Slim design : easy to maneuver over the body with an extra long handle for hard to reach area
- Rechargeable with quality Li-ion battery : 180 minutes continuous use
- 780g, ultra-light weight for easy use
- 6 Multiple massage tips



K1MED CO.,LTD

Address Rm.205, Woolim E-Biz Center 2, 12, Digital-ro 33 gil, Guro-gu, Seoul, 08377, Korea.
Tel +82-2-871-0657
Fax +82-2-871-0658
E-mail overseas@k1med.com
Website www.k1med.com

Technology and Product(s)



LUBISONE is new technology to recover lesion in body tissue.

It promotes revascularization and removing calcareous near tendons or joints though physical changing of cell membranes while stimulates musculoskeletal system where has lesion by shock wave.

- Name of product : LUBISONE
- Model Number : K1W-1501
- Rated Voltage : 230V~, 50/60Hz
- Power Consumption : 150VA
- Weight (Main unit only) : 10kg
- Dimension (Main unit only) : 397*370*200mm



KTMED INC.

Address #602, Byucksan Digital Vally 1-cha, 24, Digital-ro
27-gil, Guro-gu, Seoul, Korea
Tel +82-2-2108-7890
Fax +82-2-2108-7894
E-mail sales@ktmed.com
Website www.ktmed.com

Technology and Product(s)



NePlus(NE-SM1)

SUPER MESH NEBULIZER “Clinical Test Completed!”
(Doi: 10.1111/bcpt.12499)

- Most strong mesh(No metal dust & No hole block)
- MMAD 5µm particles sprays more than 55%
- Noiseless operation
- Preventing infection by boiling or autoclave for Mask, Mouth Piece
- No disposable kit(Reusable mask, medication bottle, battery)
- Very portable size with rechargeable Li-Ion Battery
- Residual volume is under 1% by world patented technology



AnyNeb(NE-SM2)

SUPER MESH NEBULIZER “Simple is the Best”

- Most strong mesh(No metal dust & No hole block)
- MMAD 5µm particles sprays more than 55%
- Noiseless operation
- Preventing infection by boiling or autoclave for Mask, Mouth Piece
- No disposable kit(Reusable mask, medication bottle)
- Very portable size with normal Alkaline batteries.



MAXSTAR INDUSTRIAL CO., LTD.

Address 152-12, HWANGGEUM-RO 23 BEON-GIL,
YANGCHON-EUP, GIMPO-SI, GYEONGGI-DO,
KOREA, 10048
Tel +82-31-989-3543
Fax +82-31-989-3548
E-mail overseas@maxstar.co.kr
Website www.maxstar.co.kr

Technology and Product(s)



UAM-12000: BEAUTY USE



UAM-8500 : MEDICAL USE



UAM-8100 : HOME USE



UAM-8300N : HOME USE



UAM-9306-BK: SPORTS USE



UAM-9406-BK: SPORTS USE



NEOFECT Co.,Ltd.

Address 401 West Hall, Dankook University, 152, Jukjeon-ro, Suji-gu, Yongin-si, Gyeonggi-do, Korea, Republic of
Tel +82-31-889-8521
Fax +82-31-8005-4121
E-mail contact.eu@neofect.com
Website www.neofect.com/en

Technology and Product(s)



RAPAEL Smart Glove

The wearable device-based digital solution for hand rehabilitation




RAPAEL Smart Kids

The wearable device-based digital rehabilitation solution for pediatric patients




RAPAEL Smart Board

The digital solution for upper limb rehabilitation



RAPAEL Smart Pegboard


Digital pegboard for functional and cognitive rehabilitation



openM inc

Address 413, Expo-ro, Yuseong-gu, Daejeon, 34051, Republic of Korea
Tel +82-42-823-7301
Fax +82-42 823 7302
E-mail sales@openm.com
Website www.openm.com

Technology and Product(s)



OPENCAST is a next generation, shower-friendly support frame.

Benefits of OPENCAST

Shows with cast on

OPENCAST allows daily washing and showers with cast on. Water can be easily dried off with a towel and a hair dryer.

Odor and Itch free

Perspiration quickly evaporates to prevent itchiness and unwanted odor.

Visual check

Doctors and patients can visually check the condition of the skin through openings with cast on.

Repeated wear and removal

OPENCAST can be repeatedly opened and be reworn as needed.

Stylish design

OPENCAST is made of fine materials and has a stylish design.

Comfortable and Hygienic

OPENCAST is Comfortable and Hygienic.



Optowell Co.,Ltd.

Address Optowell Co.,Ltd.
Tel +82-63-214-1186
Fax +82-63-214-1186
E-mail gabriel@optowell.com
Website www.comalgum.com

Technology and Product(s)



Comalgum is the rhinitis curer applying Low Level Laser Therapy(LLLT) which is not harm to human body. Comalgum is able to cure trouble in nasal cavity that is the root cause of rhinitis and intensified trouble in nasal inside in the meantime because it emit over 2 wave-length of laser from the laser detector. Comalgum makes easy to breath by making the breath hole in the nasal inserted area.

Comalgum is recommended for

Children

Children’s health and fitness may be impaired if they suffer from frequent colds which can lower their immunity to further infections and weaken their constitution. Children’s normal growth may be adversely affected by constant nasal infections and consequent breathing problems. Nasal congestion can inhibit concentration and lead to developmental and concentration problems.

People who taking examinations

Rhinitis, hay fever, runny noses and excessive sneezing are all well-documented causes of memory and concentration problems. The resulting fatigue also impairs exam performance.

Office workers

The office environment, particularly air-conditioned offices is a well-known source of nasal problems and infections. Acute rhinitis and allergic rhinitis frequently develop as a result of colds caught in the office



P&S Mechanics Co., Ltd.

Address Rm. 509, Woori Venture Town II, 70, Seonyu-Ro, Yeongdeungpo-Gu, Seoul, South Korea (zip) 07294
Tel +82-2-2632-3084
Fax +82-2-2632-3085
E-mail walkbot@walkbot.co.kr
Website www.walkbot.co.kr

Technology and Product(s)



Walkbot is robot-assisted gait rehabilitation training system for locomotion-impaired patients who are suffered from stroke, SCI, TBI, multiple sclerosis and cerebral palsy, etc. Products are lined up Walkbot_S (for adult), Walkbot_K (for pediatric), Walkbot_G (for both adult and pediatric, module type) and recently launched Walkbot_Premium (for adult) having the most advanced features of Walkbot series. Based on effectiveness on the training and the excellent quality of Walkbot, it has been installed in USA, Europe, Russia, Asia region and Korea.

- The world-first activated hip/knee/ankle joint drive motors based on kinematic study makes accurate and natural gait pattern close to the human being.
- Individually customized gait pattern maximizes its clinical effects
- More sessions and efficient therapy by automatic leg length adjustment
- Useful training mode (interactive, assisted training mode)
- Auditory(Metronome), visual(Front Camera), tactile(Game/Sport) stimulation for the neuroplasticity encouragement
- Various assessment tools (ROM, Stiffness, Force)
- Real-time biofeedback through augmented reality software
- CE, FDA, MET, CFDA, ISO 13485, ISO 9001 certified
- Clinically proven equipment



SEALTECH

Address 51-13 Myobong-ro 173beon-gil Idong-eup
Cheoin-gu Yongin-si Gyeonggi-do
Tel +82-31-431-7915
Fax +82-505-009-7916
E-mail evoss@esealtech.com
Website www.evoss.co.kr

Technology and Product(s)




1. Brand : EVOSS
2. Category : Anti-bedsore cushion and mattress
3. Raw material : 100% natural rubber
4. Type : Air-cell cushion and mattress
5. Product
 - 5-1. Cushion
 - 5-1-1. Octagonal 8*9 380*420*60mm
 - 5-1-2. Octagonal 7*8 460*410*80mm
 - 5-1-3. Octagonal 7*7 410*410*80mm
 - 5-1-4. Hexagonal 8*9 420*380*60mm
 - 5-1-5. Hexagonal 8*9 450*410*60mm
 - 5-1-6. Hexagonal 8*9 510*490*61mm
 - 5-1-7. Hexagonal 8*9 510*490*87mm
 - 5-2. Mattress : Octagonal 810*1840*80mm
6. Incoterms : EX-WORK, FOB, CIF, DAP, DDP
7. Payment : T/T in advance, LC
8. Price : The price depends on the quantity.
9. Country of origin : Made in Korea
10. Address : 51-13 Myobong-ro 173beon-gil Idong-eup
Cheoin-gu Yongin-si Gyeonggi-do



TAEHA MECHATRONICS Co., Ltd.

Address 421, Bakdal-ro, Manan-gu, Anyangsi, Gyeonggi-do
Tel +82-31-463-7989
Fax +82-31-463-0529
E-mail petergo@taeha.co.kr
Website www.stexfitness.com

Technology and Product(s)



Simple and User Friendly Design

- Easy to use
- 19 cm height of step
- Guide hand rail for safety
- Track-Master protocol with RS232 connection

Main Specification

- 2.5HP Treadmill Duty AC Motor
- Speed: 0.8 Km/n ~ 25 Km/h
- Elevation: 0% ~ 25%
- Maximum User Weight: 180 Kg
- Warranty: 3 years
- Dimension: W 1,010mm * L 2,100mm * H 1,150mm
- Running Space: W540 mm * L1,500mm
- Power: AC110 or AC220V / 50/60Hz
- Product Weight: 170 Kg



주월뷰텍

Address 14, Baekbeom-ro 677 beon-gil, Seo-gu, Incheon, Korea
Tel +82-32-870-0934
Fax +82-32-870-0939
E-mail sales@welbu.com
Website www.Welbu.com

Technology and Product(s)



Seven Liner ZAM ZAM will make you most beautiful body line!
 Seven Liner ZAM ZAM is a limb cyclical pressure device, helping to prevent blood clots in blood vessels with 5 chambers. In order to help increase blood circulation, periodically, high air pressure pumping to inflate the limbs and sleeves prevents swelling and daily fatigue. Especially, Hose Plug Sensor Switch (cuff connection recognition) and Cuff Deflation function (removal of remaining air in the cuff after use) are our own technology patented in Korea.



Won Industry.,Co.,Ltd.

Address 7-12,7-6,Bangadari-gill, Siheung-si,Gyeonggi-Do, KOREA
Tel +82-31-404-3093
Fax +82-31-318-5573
E-mail lead_care@naver.com
Website www.lead-care.com

Technology and Product(s)




Morning Life

- Relaxes, detoxifies, and rejuvenates; but it can also be used as a pre-workout to loosen up muscles and a post-workout to help your muscles recover by ridding them of lactic acid build-up
- Digital or analog type
- Dimension : 280x218x172mm
- Time Range : 0~60min
- Weigh : 4.2Kg




FoxyL

- Relaxes, detoxifies, and rejuvenates; but it can also be used as a pre-workout to loosen up muscles and a post-workout to help your muscles recover by ridding them of lactic acid build-up
- Dimension : 210x200x125mm
- Weigh : 3.5Kg
- Anti-overheating function
- Comfortable touch display function



Thor

- Rated voltage/Power consumption : 8.4V, 1A / 13W
- Mode : Stage 1 to 5 Dynamic Control
- Operation Time : 15 Minutes Operation Auto Timer
- Configuration : Body + Self massage toolX3
- Three types of replaceable safety bars are available in a variety of areas.
- Easy portable wireless massager




Surginet

*Feature and Benefit

- Unwoven strands
- Washable and Reusable
- Simple to use especially for head, breast, shoulder, waist and thigh
- Hold the gauge placed on the affected part comfortably even when you move.
- Free cutting for your purpose

Surgical & Medical Supplies

ACMEDI KOREA Co., Ltd
AMT Engineering Co. Ltd.
DAEJUMEDITECH
ENGINEERING CO.,LTD
DAIHAN LABTECH CO., LTD.
DAIWHA CORP., LTD.
DENTIS CO.,LTD.
ENDOVISION CO.,LTD.
Everaid Co., Ltd.
EWAH BIOMEDICS Co., Ltd.
Firson Co., Ltd
GENOSS
IDS Ltd
IMEDICOM Co., Ltd
LIVSMED Inc.
LOWTEM Co., Ltd.
MEDICALPARK CO., LTD.
MEDONICA
MEINNTECH CO., LTD.
Nanoomtech Co., Ltd.
Renosem Co.,Ltd.
SAEUM MEDITEC
SAMSUNG MEDICAL ENG. CO., LTD.
Samyang Biopharmaceuticals Corporation
ScienceMedic Co., Ltd.
Sejong Medical
SM ENG. CO., LTD.
StraTek Co.,Ltd.
SUNGKWANG MEDITECH CO.,LTD
SurgiCore Co., Ltd.
TaeWoong Medical
VOTEM CO., LTD
Wonbiongen Co., Ltd.
WONTECH CO., LTD.
YOUNG CHEMICAL CO., LTD.
ZERONE CO., LTD.



ACMEDI KOREA Co., Ltd

Address #611, Hyundai 41 Tower, 293, Mokdong dong-ro, Yangcheon-gu, Seoul, Korea

Tel +82-70-5158-3446

Fax +82-2-6971-9764

E-mail sales@acmedikorea.com

Website www.acmedikorea.com


Technology and Product(s)



VegaPlax

Advanced Plasma Technology for Medical Purpose. It is the ideal wireless Plasma device made in Korea for Blepharoplasty with Plasma Arc. It will help patients to improve eye ptosis without surgery. Moreover, it can be used for skin lifting by delivering plasma energy into the skin tissue.

- Output Frequency : 70kHz
- Power Level 3 Steps (0.7W, 1.0W, 2.0W)
- Output Mode : Continuous & Pulse
- Pulse Mode Frequency 10, 20, 40, 100, 200Hz
- Running Time : When Full charged, Max 3 hours
- Weight : Body 130g
- Size : 197 x 40 x 38mm
- Screen : 1.5" Full Color OLED
- This product is guaranteed for 1 year from the date of purchase.



AMT Engineering Co. Ltd.

Address #811, Byoksan Technopia, 560 Dunchon-daero, Jungwon-gu, Seongnam-si, Gyeonggi-do, 13230 KOREA

Tel +82-31-737-8066

Fax +82-31-737-8069

E-mail wikim@amtpro.co.kr

Website www.amtpro.co.kr

Technology and Product(s)



Lucas plus: Q-Switched ND:YAG Laser System

LUCAS plus is the best Tattoo removal & Laser Toning device.



Olive plus: 810nm Diode laser device

OLIVE PLUS is the best hair removal laser device.



A-TONE: "PREMIUM" Multi-Waved(I2PL) + PIP Toning System

A-TONE is multi-functional aesthetic devices for toning, hair removal, pigment and vascular lesions.



BOTOSONIC PLUS: Aesthetic HIFU Device

Botosonic is best Face Lifting & Tightening.



DAEJUMEDITECH ENGINEERING CO.,LTD

Address (02066) 487,Mangu-ro, Jungrang-gu, Seoul, Korea

Tel +82-2-2208-0905

Fax +82-2-438-0905

E-mail infilux@daejumedi.co.kr

Website <http://www.infilux.co.kr>

Technology and Product(s)



Co2 Fractional Laser

It is an ULTRA PULSE CO2 FRACTIONAL LASER , which is a new generational version of the CO2 Fractional laser and its wave length is 10600nm, which was replaced by a stronger laser for severe scars, enlarged pores etc. It is Ultra Pulse CO2 Fractional Laser, which is a new generational version of the CO2 Fractional laser for severe scars, enlarged pores, and various kinds of skin disorders.



Diode Hair removal Laser 600W

- What is hair removal laser?
Our Diode Hair removal Laser 600w has high energy pulsed 808nm diode laser for hair removal. Aroma machine can reduce pain while reducing skin damage by Cooling effect to reduce skin's pain and can get perfect result by quick and simple treatment of hair removal.



HIFU for face lifting and tightening and body

- What is HIFU?
HIPRO is HIFU machine that has various cartridges and HIPRO delivers HIFU energy to SMAS layer which is located between subcutaneous tissues and muscle. And then, make 0.7~1.0mm size of TCP(Thermal Coagulation Points) zone without surgery or any disrupting skin surface to stimulate collagen growth.



DAIHAN LABTECH CO., LTD.

Address 247-16 GAGOK-RI, HWADO-EUP, NAMYANGJU-CITY, KYONGGI-DO, KOREA

Tel +82-31-593-5961

Fax +82-31-593-9563

E-mail Info2@labtech.co.kr

Website www.labtech.co.kr

Technology and Product(s)



Bench Top Autoclave (Class B)

With automatic dual door & LCD monitor display for whole procedure.
Available capacity: 12/25/40 liters



Medium & Large autoclave

With pre & post vacuum system & strong hatch door for user safety
Available capacity: 40/60/80/100/200/300 liters



Pass through door autoclave (Class B)

With full touch screen & round/square chamber selectable
Available capacity: 200/300/500 liters



Plasma Sterilizer

With full touch screen & square chamber, it's suitable for sterilizing medical tools/equipments that are weak at heat/high temperature such as endoscopes, catheters, etc.

DAIWHA
www.daiwha.com



DAIWHA CORP., LTD.

Address Daiwha Bldg, 3, Yeoksam-ro 33 gil, Gangnam-gu, Seoul, South Korea
Tel +82-2-558-1711
Fax +82-2-554-0317
E-mail overseas03@daiwha.com
Website www.daiwha.com

Technology and Product(s)



DI-6000/DS-5000 (Open/Close System)

- Innovative Design
- High Convenience – TFT-LCD, Jog Dial
- Various Safety Functions : IV Set Loading(Infusion), Easy Menu, Auto Syringe Detection(Syringe)
- 1~1500ml/h Range of flow rate
- Smart Docking Station System using multiple pumps(~4 units).
- Battery : Operating 10 hours, Recharging 6 hours
- Weight : Approx. 1.9(±0.05)kg



DE-400 MEDITOM RF Generator with APC

- Modern Design : improved convenience LCD touch screen
- Various Safety Function : D.R.M(Daiwha patient Return plate monitoring system)
- Various Output mode
 - 2 modes of Bipolar Cutting(APC-General Cut)
 - 2 modes of Bipolar Coagulation(APC-Standard, Low)
 - 4 modes of Mono-polar Cutting(Pure,Blend 1/2/3)
 - 3 modes of Mono-Polar Coagulation(Forced, soft Low)
- 2 Monopolar and 2 Bipolar Plug with play instrument



DUB-100 & DB-137/ DU137/127/Scissor

- DUB-100 : 1 Ultrasonic Plug, 1 Bipolar Plug, 1 Monopolar Plug
- DU-137 : Bipolar Vessel Sealing Energy Device
- DU-137 : Ultrasonic Energy Device
- DU-127/Scissor : Ultrasonic Open-Type Energy Device
- At Target Hospitals- early 2019

LUVIS



DENTIS CO.,LTD.

Address Seong Seo Seo Ro 99, DalSeo Gu, DaeGu, Korea
Tel +82-2-869-6020
Fax +82-2-869-6021
E-mail luvis@dentis.co.kr
Website http://luvis.co.kr/eng/index.php

Technology and Product(s)



LUVIS-L200(Surgical light)

- Intensity : 160,000 Lux
- Color Temperature : 3800K / 4300K / 4800K
- CRI : 96



LUVIS-M200(Surgical light)

- Intensity : 120,000 Lux
- Color Temperature : 4300K
- CRI : 95



LUVIS-S200(Minor surgical light)

- Intensity : 100,000 Lux
- Color Temperature : 3800K / 4300K / 4800K
- CRI : 90



LUVIS-E100(Examination light)

- Intensity : 60,000 Lux
- Color Temperature : 3500K / 4300K / 4800K
- CRI : 90

endovision
infinite vision of the medical device



ENDOVISION CO.,LTD.


Address #214, 2nd Floor MJ Technopark, 20, Dalseo-daero 109-gil, Dalseo-gu, Daegu, 42709, Korea
Tel +82-53-751-1955
Fax +82-53-751-1964
E-mail braelyn@endovision.co.kr
Website www.endovision.co.kr

Technology and Product(s)



NEXOSEAL

Hemorrhagic hemostasis product during general surgery, easy use during operation and securing visibility, fast hemostatic effect and excellent prevention of additional hemostasis



DUAL PATTY

Bleeding management during spinal or brain surgery enables quick and accurate management of bleeding during surgery and easy protection of surrounding organs



UBE Tool kit

A set of bi-directional spinal endoscopic instruments for various spine surgery.



Flex-car R/C

RACZ catheter used for neuroplasty (PEN) with excellent pushability and trackability



TOME FORCEP

Rongeur used to remove bones or remove tissue during spinal surgery



FREETOME PUNCH

It is a punch used for removal of bone and tissue during spine surgery. It is convenient to rotate due to shaft rotation

EVERAID



Everaid Co., Ltd.

Address 36, Dulleom-gil, Yangsan-si, Gyeongsangnam-do, Korea
Tel +82-55-363-2122
Fax +82-55-363-4683
E-mail hq@everaid.com
Website www.everaid.com

Technology and Product(s)



AIDFoam is HDI type hydrophilic foam dressing. This HDI foam dressing has a high absorption & retention capacity with its special fluid-lock structure. The foam can absorb fluid about 10 to 14 times of its own weight (EN13726-1), and minimize exudates leaking. The high MVTR top film, about 2000~2500g/24hr (EN13726-2), increase exudates control capacity.

On top of that, AIDFoam comfort, silicone adhesive foam dressing has very high tack adhesive silicone, adhering well on body contours. The silicone adhesive provides long wear time and minimizing pain and damage in changing dressing.



AIDerm is high MVTR PU transparent film dressing. The sterile thin and non-glossy film dressing can perfectly protect skin and wound from outside contaminants such as bacteria, germ and fluid to prevent 2nd infection.



NeoDerm is high tack adhesive hydrocolloid dressing. This high tack adhesive make the dressing adhere well to skin and minimize skin trouble compared to common hydrocolloid dressings. NeoDerm Changes absorbed exudate to gel promptly to create optimal moist wound environment. The high MVTR top film increases exudate control capacity as well.



EWHA BIOMEDICS Co., Ltd.

Address #901, 9F, Sungwoo BLDG, 49, Mapodaero, Mapo-gu, Seoul, South Korea

Tel +82-2-2055-3820

Fax +82-2-2055-3819


E-mail trade@ewhainc.co.kr

Website http://www.anapa.com/en_index.php

Technology and Product(s)



disposable infusion device ANAPA & ANAPLUS
ANAPA is the world's first and unique disposable infusion device using self-produced CO₂ as power source. By its unique design, it is safe of the many shortcomings of mechanical, electronic and elastomeric devices. It offers the best combination of accuracy, safety and convenience in pain palliation.



- 1ml graduation for easy monitoring of infused medication volume.
- Negligible amount of medication remaining in device after infusion.
- Patented end-cap that allows for safe and time saving self-priming.



ANAPA CTx - Chemotherapy is the world's first and unique disposable infusion device using self-produced CO₂ as power source. By its unique design, it is safe of the many shortcomings of mechanical, electronic and elastomeric devices. It offers the best combination of accuracy, safety and convenience in chemotherapy.

- Unique filter to prevent crystallization of chemotherapy's medication.



Firson Co., Ltd

Address 47, Baekseokgongdan 1-ro, Seobuk-gu, Cheonan-si, Chungcheongnam-do, Korea

Tel +82-1833-8200

Fax +82-2-2138-8280

E-mail global@firson.co.kr

Website www.firson.co.kr

Technology and Product(s)



- Use life : Up to 14 days
- Shelf life : 24 months
- Storage information : 15~30 °C
- Open bottle shelf life : 75 days

- Rapid disinfection (5min)
- Proven safety & Efficacy (CE approval)
- High reusability
- Cost-saving

GENOSS
For Patients & Doctors



GENOSS

Address 1st floor, Gyeonggi R&DB Center, 105,
Gwanggyo-ro, Yeongtong-gu, Suwon-si, Gyeonggi-
do, South Korea

Tel +82-10-3348-9391

Fax +82-31-888-5136

E-mail yrchoi@genoss.com

Website www.genoss.com

Technology and Product(s)



GENOSS MONALISA Filler

MONALISA is a high purity hyaluronic acid filler which keeps the excellent volume effect for a long time at the beginning of the procedure, being decomposed naturally over time. Depending on particle size, you can choose appropriate products from soft, mild to hard type.



GENOSS Inflator B30

- Easy fixing and releasing with the button type Ergonomic design for comfortable handling
- 20cc barrel allows rapid and easy deflation



GENOSS Inflator B40

- Easy fixing and releasing with the button type Ergonomic design for comfortable handling
- 30cc barrel allows rapid and easy deflation



GENOSS PTCA Balloon Catheter

- Great pushability and trackability with PTFE coated hypotube
- High crossability with ultra-thin profile
- Specially engineered tri-wrapped balloon for excellent re-wrapping

ids
Innovation for better life



IDS Ltd

Address #105, Byoksan Digital Valley 3 cha, 271, Digital-
ro, Guro-gu, Seoul, Korea

Tel +82-2-2025-1151

Fax +82-2-2025-1154

E-mail info@ids.co.kr

Website www.ids.co.kr

Technology and Product(s)



SmaXel

SmaXel is a scanning type fractional ultra-pulse CO2 laser, in which the surgeon can control freely Scan Area, Density Level and Beam Degree. SmaXel is an advanced ultra-pulse CO2 laser, which increases treatment effects while minimizing the influence to neighboring tissues as a fine laser beam is irradiated.



Q-Switched Nd:YAG laser (Max : 1.8J)

Application : Tattoos, Vascular lesion Removal, Nevus of OTA, Freckle, Nano peeling, Quasi long peeling, Hair removal using long pulse. Etc



Alexandrite

This hair removal system is a light therapy device specifically designed to eliminate undesired hair from all parts of the body.

- Easy operation
- Less painful
- Treatment: Arm, Beard, Breast, Armpit, Groin, Leg



Bu:tera

High-Intensity Focused Ultrasound (HIFU) is a highly precise medical procedure that applies high-intensity focused ultrasound energy to locally heat and destroy diseased or damaged tissue through ablation.



IMEDICOM Co., Ltd

Address #612, 172 LS-ro, Gunpo-si, Gyeonggi-do,
Republic of Korea

Tel +82-31-479-1156

Fax +82-31-479-1157

E-mail leejames@imedicom.co.kr

Website www.imedicom.co.kr

Technology and Product(s)



Balloon Kyphoplasty System

MEDINAUT/MEDINAUT Plus Balloon Kyphoplasty System is intended to be used as conventional bone tamps for the reduction of Vertebral Compression Fractures (VCFs) and/or creation of a void in cancellous bone in the spine during Balloon Kyphoplasty.



Epidural Catheter

For patients who are afraid to take spine surgery against disc disease or spinal stenosis, SPINAUT-E & SPINAUT-V are good alternative surgical procedures for patients and surgeons. At lower spinal column, SPINAUT-E & SPINAUT-V treat pains with lysis and delivering pain medication.



Powered Acetabular Cup Remover

EZX helps you remove an acetabular cup safely, precisely, easily and efficiently. You can extract a cup with less bone, tissue loss and prevent the creation of additional acetabular fracture.



Surgical Power Handpiece

IMEDICOM has two kinds of power tool system, ZAGUAR and Leopard. ZAGUAR powered instruments system is battery operated system and Leopard power tool system is pneumatic system which is operated by compressed air or nitrogen. These two system is intended for Trauma and Joint replacement surgery.



LIVSMED Inc.

Address 304, D-dong, Pangyo-ro 700, Bundang-gu,
Seongnam-si, Gyeonggi-do 13516, Republic of
Korea


Tel +82-70-4450-2717

Fax +82-31-706-3211

E-mail saebomi_lee@livsmed.com

Website www.livsmed.com

Technology and Product(s)



Since the existing general laparoscopic instruments have no joint structure, instruments cannot approach at certain angle such as side of tissue or lesion. In other words, motions are very limited during operation and surgeons need more time to operate.

- It is newly developed the multi-joint and articulating laparoscopic surgical instrument that can overcome the drawbacks of a straight-type device and allow intuitive manipulation. This product is a manual laparoscopic surgical instrument with intuitive multi-joint function. Product name is <ArtiSential>.
- Due to the multi-joint structure of end tool, <ArtiSential> can move 90 degrees to upward, downward, leftward and rightward. The movement of the end tool, the movement of the control part, and the movement of user's wrist and fingers are exactly synchronized. Because of these advantages, intuitive articulating movement can be achieved.
- And <ArtiSential> allows access to a narrow surgical site and allows access to a variety of angles that conventional laparoscopic surgical instruments can not.



LOWTEM Co., Ltd.

Address #A-301, 302, Gunpo IT Valley, 17, Gosan-ro
148beon-gil, Gunpo-si, Gyeonggi-do, Korea
Tel +82-31-427-0763
Fax +82-31-427-0793
E-mail lowtem@lowtem.com
Website www.lowtem.com

Technology and Product(s)



LOWTEM Series plasma sterilizer Crystal/Smart 30 (33Liters) / Crystal 40R (40Liters) / Crystal 50 (50Liters) / Crystal 100 (100Liters) /Crystal 120M(D) (130Liters) /Smart 150(D) (150Liters) LOWTEM plasma system uses low-temperature, hydrogen peroxide and plasma technology to sterilize a wide range of medical instrument effectively and safely.

- Rapid Warmup & Dry System : The world's first patented allows to overcome the humidity related problem, removing residual moisture and chamber's cold point.
- Convenient and Safe Sterilizing Agent : Multi-cycle "Tank type" sterilant is safe for users to store, to deliver, and to use. The automatic RFID system for using sterilant tank is convenient and reliable for users.
- Auto Interlock & Open System : One gentle push makes the chamber door locked, and the chamber door can be opened automatically when the door sensor detects finger movement.
- Plug and Play : Warming up time is very rapid (Maximum 15 min) after main power is on, so that standby power consumption is zero during off-duty or overnight.
- USB History Memory : Sterilization history can be downloaded to your USB memory stick and reprint anytime you want to.
- Login Function : Providing User ID and PW, so that hospital can manage users who access the sterilizer.



MEDICALPARK CO., LTD.

Address #601/623/624/901, Knowledge-Industry Center,
Bundang Suji U-Tower, 767, Sinsu-ro, Suji-gu,
Yongin-si, Gyeonggi-do, Korea
Tel +82-70-4800-5660
Fax +82-70-8266-1108
E-mail info@medicalpark.co.kr
Website http://bexcore.com

Technology and Product(s)



Bexcore is to inspect or remove a tissue without leaving any incision or scar. It is a method to make a 4~5mm incision and aspirate tissue into a needle by placing it near to a lesion and then the tissue cut by an inner blade of the needle is discharged automatically by vacuum suction. By real-time ultrasound, an accurate diagnosis and treatment are possible.



MEIDONICA

Address 37-9, Techno 1-ro, Yuseong-Gu, Daejeon, Korea
Tel +82-2-808-0234
Fax +82-2-808-0235
E-mail info@medonica.com
Website www.medonica.com

Technology and Product(s)



Plasma sterilizer ZEUS-50 & 130S

Three kinds of cycles can be selected depending on the application and the load volume
 Short cycle : 25minutes (± 5) : ZEUS-50 & 130S
 Standard cycle : 35 minutes (± 5) : ZEUS-50 / 45 minutes (± 5) : ZEUS-130S,
 Advanced cycle : 45 minutes (± 5) : ZEUS-50 / 55 minutes (± 5) : ZEUS-130S,
 Installable in a small area
 Operates with economy and efficiency even with a small load : however, a large load can also be handled
 The best low temperature sterilizer to start with a reasonable price



MEINTECH CO., LTD.

Address 5th Floor, A-dong, Digital Empire B/D, 387 Simindaero, Dongan-gu, Anyang-si, Gyeonggi-do, Korea 14057
Tel +82-31-381-7076
Fax +82-31-381-7053
E-mail overseas@meinntech.com
Website <http://www.meinntech.com/eng/main>

Technology and Product(s)



ANYFUSION is the product combining the syringe pump and infusion pump into one unit - Cylinder pump that solved problems of accidental bolus, free flow, inadvertent bolus, siphoning phenomenon etc. that lasted for over 60 years due to the technical limitation by using rotary motion of pistons on the contrary to direct motion; and capable of medicine injection such as contrast medium with continuous accuracy and precise volume control. Moreover, the cylinder cartridge maintains high accuracy within $\pm 2\%$ deviation for long period without replacing the cartridge for 3days, thus reduces unnecessary in use that distress for users; and increases optimization of and convenience of application in use as the device is not affected by gravity nor its location with advantage of own mechanism.

- High accuracy within $\pm 2\%$ of deviation
- Injection capable regardless of location
- Directly connected to the container without frequent replacement of consumables
- Fundamental prevention of free flow, accidental bolus, inadvertent bolus and siphoning phenomenon
- Minimum alarm, removing stress on over/under dose
- Integrated syringe & infusion mode use



nanomtech Co., Ltd.

Address 57, Cheomdan Venture So-ro, Buk-gu, Gwangju, 61003, South Korea
Tel +82-62-955-8588
Fax +82-62-955-8589
E-mail altoall@altochem.com
Website www.nanoomtech.co.kr

Technology and Product(s)



Heart+ResQNT-381 is a compact, portable device that analyzes the heart's rhythm and delivers a defibrillation shock to re-establish the normal heart rhythm. Heart+ResQNT-381 guides the user through each step of the defibrillation process by providing voice/ visual prompts. Heart+ResQNT-381 of Nanoomtech Co., Ltd provides electrode pads that can be used for both adults and children.
 -5 year-warranty for main unit, 2year-warranty for electrode pads.



Renosem Co.,Ltd.

Address Unit 103-806, 397, Seokcheon-ro, Bucheon-si, Gyeonggi-do, Korea
Tel +82-32-222-0180
Fax +82-32-222-0110
E-mail sale@renosem.com
Website www.renosem.com

Technology and Product(s)



RENO PLASMA STERILIZER is the Best Solution for Infection Control.
 Proven Strong Penetrability
 Fast Sterilization Cycle & Gentle Sterilization
 Non Toxic Residue
 User-friendly Interface & Low running Cost




SAEUM MEDITEC

Address 102-1406, 36, Bucheon-ro 198beon-gil, Bucheon-si, Gyeonggi-do, 14557, Republic of Korea

Tel +82-32-661-0804

Fax +82-32-661-0805

E-mail admin@saeum.net

Website WWW.SAEUMMEDITEC.COM

Technology and Product(s)



EPI-Detection® is Medical pressure monitor to detect the epidural space

The EPI-Detection® is Medical pressure monitor and use it to detect the epidural space for accuracy of procedure of "Epidural block"

When the Epidural needle reaches the epidural space, it detects the pressure change of the muscle layer and epidural space. the pressure sensor recognizes it, and then displayed on LED so that the user can identify it.

It has a structural feature that allows the Racz's catheter to be inserted right away.

To put it simply, the device supports the safe treatment of patient by allowing the Epidural needle to precisely reach the epidural space.

This product was developed in collaboration with Severance Hospital, and so far clinical research has been actively carried out. According to the study 100 cases (EPI-Detection ver.1) in 2016, epi- detection was approved its effectiveness and stability. With this result, In July, 2017, we obtained a KOTRAN FDA approval for a medical pressure gauge.




SAMSUNG MEDICAL ENG. CO., LTD.

Address 12, OKGUCHEONSEO-RO 237BEON-GIL, SIHEUNG-SI, GYEONGGI-DO, KOREA

Tel +82-31-499-0679

Fax +82-31-499-0675

E-mail s040431@chol.com

Website www.ssmed.net

Technology and Product(s)



Wall suction

SME Wall suction is collecting a flood, spill, phlegm and flood phlegm from the patient body after surgery, there are 2 regulated vacuum range, 200mmHg and 700mmHg. It contains unbreakable polycarbonate 1200cc collection bottle, over-flow protector, hose and adaptor, Has CE marked.



Flowmeter with humidifier

SME Flowmeter with humidifier connects a gas outlet and supplies an oxygen or air to the patient who is needed it in the ward, ICU, CCU, OR and ER. The flow rate is 0~15LPM. It contains unbreakable polycarbonate bottle, audible alarm and adaptor. Has CE marked.



Gas outlet

SME Gas outlets are installed on the wall, bed head unit and pendant where needs to delivery a gas to the patient. All of gas outlet obtained the CE mark.



Bed head unit

SME bed head unit is a universal type, and shall be contained power socket, medical gas outlet, nurse call as buyer's requirement. It is customized products.



Samyang Biopharmaceuticals Corporation


Address 295 Pangyo-ro, Bundang-gu, Seongnam-si, Gyeonggi-do 13488, Korea
Tel +82-2-2157-9817
Fax +82-2-2157-9062
E-mail www.samyangbiopharm.com
Website sunwoo.kim@samyang.com

Technology and Product(s)




Synthetic Suture Material (Bulk Thread)

- Absorbable Braided Suture : TRISORB, NEOSORB
- Absorbable Mono Filament Suture : MONOSORB, MONOFAST



Surgical Mesh for Hernia Repair

- Proflex Mesh : Partially Absorbable PP+PGCL mesh



Hemostat

- SurgiGuard Oxidized Regenerated Cellulose (ORC) Hemostat



Adhesion Barrier

- InterGuard for Adhesion Barrier



ScienceMedic Co., Ltd.

Address #5F ScienceBuilding 104, 24beon-gil, Jeonpa-ro, Anyang City, Gyeonggi-do, Korea
Tel +82-31-479-1194
Fax +82-31-479-1194
E-mail sgjun@sciencemedic.com
Website www.sciencemedic.com

Technology and Product(s)



Skin Closure

- No scar for Staple or Suture
- Stronger than sutures in protecting incision from distraction forces
- No staple or suture tissue punctures that may allow bacteria into wound site
- No need to visit Hospital or clinic for removal(Remove patient home)
- Faster than sutures
- Fewer wound complications



Battery Powered Handpiece

Battery-typed power Handpiece is portable, so that it is suitable for delicate surgical operating.



Saw Blade

Saw Blades have a sharp edge which are used for cutting bone by vibrating horizontally. There are Linvatec, Stryker, Aesculap, Synthes, Desoutter and etc. and we could manufacture all spec that customers want. It is very competitive for the price and we are supplying over 60% of local market.



Pneumatic Powered Handpiece

Power Handpiece is used for drilling, reaming and pinning a bone and there are drill, reamer, wire, pin and saw attachments for the Multi Handpiece can be assembled & disassembled easily. It is ergonomic design, strong power, durability and also cheap in price.



Sejong Medical

Address 11 Sinchon 2-ro, Paju-si, Gyeonggi-do, 10880 South Korea
Tel +82-31-945-8191
Fax +82-31-945-8190
E-mail info@sejongmedical.com
Website www.sejongmedical.com

Technology and Product(s)



Laparoscopic Trocar



Laparoscopic Specimen Retrieval Bag



Single Port for Minimal Invasive Surgery



Laparoscopic Suction & Irrigation



SM ENG. CO., LTD.

Address SM ENG. CO., LTD.
Tel +82-51-305-8016
Fax +82-51-305-8021
E-mail sales@sm-eng.net
Website www.sm-eng.net

Technology and Product(s)

Suture needles made of SUS 304 stainless steel are more flexible and has high penetration performance. With our high technology, we also guarantee solidity of needle and thread attachment. We offer complete range of absorbable and non-absorbable suture.



StraTek Co.,Ltd.

Address RM705~707, 721~724, Yucheon Factopia, 30, Jeonpa-ro, Manan-gu, Anyang-si, Gyeonggi-do, Rep. of Korea 14086
Tel +82-31-467-3137~8
Fax +82-31-467-3139
E-mail www.stratek.co.kr
Website stratek1998@gmail.com

Technology and Product(s)



Model: SPMI-330
: Combination Electrical Therapy (ICT+SSP)



Model: STS-1000
Muscle Tap Apparatus (for massage)



Model: CPS-2000
: Knee C.P.M



Model: MC-1
: Arm & Leg Trainer active or passive motion for the rehabilitation



Model: 3D-CENTOUR
: 3D Spinal Stabilization Device Computer System



STL-101
: Digital Therapy Laser (for pain relief)



SUNGKWANG MEDITECH CO.,LTD

Address 1014 - 20, Bukhansan-ro, Jangheung-myeon, Yangju-si, Gyeonggi-do, Korea
Tel +82-31-855-0627
Fax +82-31-855-7549
E-mail skmeditech@skmeditech.com
Website www.skmeditech.com

Technology and Product(s)



SurgiCore Co., Ltd.

Address Bldg. A, 43-10, Jingoksandan 5beon-Ro,
Gwangsan-Gu, Gwangju 62464 Korea

Tel +82-62-714-3850

Fax +82-62-714-3861

E-mail www.surgicore.co.kr

Website www.surgicore.co.kr



Technology and Product(s)



LigaForce

[Specification]

Shaft diameter: 10mm, Clip: M/L, 20 implantable clips preloaded.

[Features]

- 340mm shaft length
- The safety interlock
- Anti clip drop system
- Ergonomic grip design



SurgiTractor

[Available sizes]

XX-Small, X-Small, Small, Medium, Large, X-Large, XX-Large

[Features]

- Providing pure circle-Maximum opening with minimum incision
- Easy to fold
- 12mm port equipped high transparency silicone cover is available
- Easy to attach and detach the cover to the retractor



SurgiPort

[Available sizes]

5.5mm, 11.5mm

[Features]

- Breakage free
- Ergonomic design of the obturator head
- Funnel type cannula entry
- Patented screwing structure
- Easy to screw the device into the thoracic space

TaeWoong Medical

Address Bldg.C 4F, Daebang Triplaon Business Tower,
158 Haneulmaeul-ro, Ilsandong-Gu, Goyang-Si,
Gyeonggi-Do, Korea 10355

Tel +82-31-904-6196 (ext.622)


Fax +82-31-904-6157

E-mail contact@stent.net

Website www.taewoongmedical.com




Technology and Product(s)



TaeWoong Niti-S Biliary Stents

- Consists of the implantable metallic stent and the introducer system.
- Both Covered and Uncovered stents are available.
- The D-Type Stent is an exclusive design which maintains full luminal patency in the tortuous biliary anatomy for optimal drainage.



TaeWoong Niti-S Esophageal Stents

- Consists of the implantable metallic stent and the introducer system.
- Different types and wide ranges of stents provide various options to clinicians for different indications.
- The Double-Type Stent, Proximal release delivery system, and the TTS delivery system are our exclusive designs.



TaeWoong Niti-S Duodenal / Colonic Stent

- Consists of the implantable metallic stent and the introducer system.
- Both Covered and Uncovered stents are available.
- Palliative treatment for malignant Pyloric / Duodenal and colorectal obstructions are available.

VO+EM
www.votem.kr



VOTEM CO., LTD

Address 27, Geodudanji 1-gil, Dongnae-myeon,
Chuncheon-Si, Gangwon-Do, Korea 24398
Tel +82-33-910-0701
Fax +82-33-911-0701
E-mail votem@votem.kr
Website www.votem.kr

Technology and Product(s)



VP-700 / VP-1000 / VP-1200

Display : 7"/10.4"/12.1 TFT LCD(800X600), 10 wave trace on/pff, Full ECG 7 wave display
Parameters : ECG, SpO2, NIBP, 4IBP, Resp., 4TEMP, EtCO2, ICO, EKG, Multi-gas, BIS
Interface : RS-232 port, VGA port, ECG output(opt.) port, SD memory card, 12 Languages
Feature : 7days trend data, 6 kinds of virtual screens

VM15, Modular Patient Monitor

Display : 15.1" TFT LCD
4.3" Transport Monitor, Easy S/W upgrade using USB (16GB)

WON No.1 In the World
BIOGEN

Wonbiongen Co., Ltd.

Address 56, 1 Gongdan-ro 5-gil, Gumi-si,
Gyeongsangbuk-do, Korea
Tel +82-54-462-7965
Fax +82-54-471-7963
E-mail wonbiogen@wonbiogen.co.kr
Website http://www.wonbiogen.co.kr

Technology and Product(s)




Therasorb algiplus Hydrophilic/ Adhesive dressing

Alginate ingredient, when mixed(interacted) with blood in the wound, forms alginate-blood gel and exchange ions with exudation in the wound, thereby producing an effect of stopping bleeding.




Therasorb AD

Therasorb AD is coated with silicone adhesive which is a lot of holes on top of polyurethane foam dressing, allowing for self-driven adhesion. Thus secondary dressing is not needed for fixation.



Therasorb Net

Porous film on wound contact layer does not leave residue on wound site during dressing replacement, and does not cause pain during replacement.
- Prous film minimizes expansion, allowing for excellent retention of exudate, leading to optimized humidity.



Therasorb Silver Hydrophilic Dressing

Polyurethane foam containing silver that has antibacterial effect and hemostasis effect of alginate which included in Therasorb AlgiPlus.



WONTECH CO., LTD.

Address 3F, Pangyo Global R&D Center, 22, Daewangpangyo-ro 712beon-gil, Bundang-gu, Seongnam-si, Gyeonggi-do, Republic of Korea / 13488
Tel +82-31-1670-1450
Fax +82-70-7882-8658
E-mail neosoul7@wtlaser.com
Website www.wtlaser.com

Technology and Product(s)



VEINCARE

Varicose veins laser treatment veincare(1470nm)



VEINCARE1940

New generation for treating varicose veins by using 1940nm wavelength.



GI LASER

ESD (Endoscope submucosal dissection) Laser for colorectal cancer and gastric cancer



THY_S

New generation of surgical thyroid treatment.



YOUNG CHEMICAL CO., LTD.

Address 80-93, Golden root-ro, Juchon-myeon, Gimhae-si, Gyeongsangnam-do, Korea
Tel +82-55-365-1361
Fax +82-55-323-9523
E-mail yc@youngchemical.co.kr
Website www.youngchemical.co.kr

Technology and Product(s)





ZERONE CO., LTD.

Address (Shinil IT UTO Bldg., Dangeong-dong) #810, LS-Ro 13, Gunpo-Si, Gyeonggi-Do, Korea

Tel +82-31-453-3512

Fax +82-31-453-5444

E-mail eddie@01zeus.com

Website www.01zeus.com

Technology and Product(s)



ARGON-Z

Argon Plasma Coagulator, AC100V~AC240V, 50/60Hz, Gas flow: 0.1l/min~9.9l/min, Class 1, Type CF.



ZEUS VISION

400W, 50/60Hz, 99 memories, Class 1, Type CF, Endo-Cut(Intra Cut), Bipolar T.U.R.P.



SENS-Z

Specially designed for E.N.T. Available cut and coagulation in single unit. 30W AC100V~AC240V, 50/60Hz, 10 memories, Class 1, Type CF.



ULTRA-Z

Ultrasonic Surgical Unit for Aesthetic lipo-sculpture, 90W AC110V or AC 230V, 50/60Hz, Class 1, Type BF.



INFU-Z 1000

AC100~240V, 50/60Hz, DC12~18V, Class 1, Type CF
 IPX2 standard, F. range: 0.1~999.9ml/h, V. range: 0.1~9999ml,
 Accuracy: ± 5%, Macro/Micro mode, Battery: Li-ion(OT: 10 hours, RT: 8 hours)