



Multifunctional Feminine Medical Equipment

D+.Lady



BUHEUNG MEDICAL Co., Ltd.

Since 1980~

[Dr.Lady] Clinical trial on SCI Journal !

Dr.Lady

Successful results in clinical trial to Urinary incontinence, Sexual function and Premenstrual syndrome, Listed on SCI journal as a therapeutic device.




Clinical information

Title : 『Effect of Low-Frequency intravaginal electrical stimulation on female urinary incontinence, quality of life, and urinary symptoms』

- **Clinical Institution :** The catholic University of Korea Department of Urology, Seoul St.Mary's Hospital, Kandong Sacred Heart Hospital, Catholic University College of Medicine, Uijeongbu St.Mary's Hospital
- **Clinical period :** May, 2017 ~ Feb, 2018
- **Journal :** LUTS (Lower Urinary Tract Symptoms) / Wiley
- **Registered :** Aug, 2019

Results

- 83.87% improvement in urinary incontinence, acute incontinence, unpredictable incontinence
- All patient symptoms improved after 8 weeks of treatment
- Muscle fiber contraction improvement
- Vaginal muscle contraction improvement
- Premenstrual syndrome improvement
- Physical and social limitation improvement
- Personal relationship, emotions, sleep and vitality, role restriction improvement



Dr.Lady, for the happiness of all ladies



- Top level Safety System
- PCT Patents
- Comfortable, easy and Private Treatment
- For Lady's Satisfaction
- 40 years' Electro Therapy Technology Aggregation
- The state of Art Medical device that no one has experienced

Product introduction

MK-000A Dr. Lady is Therapeutic Electro Device for Women who has Incontinence, Sexual Dysfunction, Vagina Atony, Leukorrhea Flow, Uterus Mass, Genital Cold, Uterine Myoma, Menstrual Pain, Hemorrhoids (Piles). The insertion Style Probe has 6 kinds of combinational Program with 3 Functions (Electro Frequency, Heating, and also Fine Vibration). These functions' combination makes treatment efficacy to woman diseases as Incontinence, Sexual Dysfunction, Vagina Atony, Leukorrhea Flow, Uterus Mass, Genital Cold, Uterine Myoma, Menstrual Pain, Hemorrhoids (Piles). By Contrast, outer electrodes' Programs are specialized in Menstrual Pain reduce. It has 2 programs and the frequencies for each program are very different to each other. Uniquely, this Therapeutic Medical Device is wearable. Very small and light (63mm x 88mm x 29mm / 76gram) so you can treat with easy living activities, at any time & any place. Furthermore, it is hidden in wears as inner wear, T-shirt, sweater, neat or skirt. You don't have to stay beside huge device and also you don't have to mind other people when you treat. We, Buheung Medical is very pleased to give you MK-000A Dr. Lady, the state of art medical Device and serving happier life to Ladies in the world.



MK-000A Dr. Lady

Efficacy



Unique Feature

- Very Small & Light / PORTABLE & WEARABLE Therapeutic Medical Device
- Probe for Vagina Insertion & outer electrode pads for Exterior Skin
- Titanium Plate & Poly Carbonate
- Automatic Program System
- UV Ray Sterilization
- Graphic LCD Display
- Alert System with Colour Display
- Touch Panel UI (User Interface)
- Instruction Style UI
- 3 Users Storage Available
- Skin Connection Checking System
- Specialized Age Programs. (Programs for 10–29 Yrs / 30~49Yrs / over 50 Yrs)
- Function memory & Recall
- Rechargeable battery
- Far–infrared Radiation
- Language Setting available


Electronic Frequency treatment Low, Mid, Interferential frequency (1~1560Hz)

Muscle Rehabilitation

- Urinary incontinence, incontinence, Constipation treatment
- Vaginal muscle and pelvic floor muscle strengthening
- Sexual function improvement
- Intestine function improvement
- lower abdomen deflection improvement
- Residual urination reduction
- Increased blood flow, enhanced vascular muscle
- Kegel exercise

Nerve Rehabilitation

- Urinary incontinence treatment
- Sexual dysfunction treatment
- Intestine dysfunction treatment




Treatment from the inside

- Treatment of vaginal and pelvic muscles with direct access
- 360° cubical frequency effect around probe, muscles and nerves are treated to deep areas
- Optimal effects with FES frequencies (Functional Electric – stimulation)
- Very rapid effect on Genital cold

Increase vaginal temperature

Hyperthermia treatment

- Inflammation treatment and Prevention
- Blood circulation improvement
- Cell function improvement
- Uterine myoma treatment and prevention



- Direct heat stimulation
- 360° cubical effect around probe, muscles and nerves are treated to deep areas
- Accurate heat of 45°C, 42°C

Micro vibration

- Strengthen the effect of Electronic frequency treatment
- Blood circulation improvement
- Nerve rehabilitation
- Helps to establish a position in the vagina


Premenstrual syndrome treatment



Menstrual pain reduce with electronic frequency, heating and vibrating

- Vaginal muscle spasm treatment program
- Pad program for non-insertable people

UV Sterilization Effect



99.9% sterilization effect

- UV Reflective material pouch
- Just press – button to use
- Automaticly turned off after sterilization

Main Unit

- Graphic LCD Screen Display
- Touch Panel UI
- Li-Polymer Rechargeable Battery
- 4 Outer Titanium electrodes (for Menstrual Pain Reduce)
- Very Small & light
- Clip Style Wearable Unit



Graphic LCD Display

- Instruction Guide Style
- Language Setting
- Alert System with Colour Display
- Touch Panel UI (User's Interface)
- No Button, Just Touch, Operational Error Prevention System



Urinary Insertion Probe

- Titanium Plate + Polycarbonate
- Semi-Permanent Use
- Electro Frequency + Heat + Fine Vibration



Product Specification

- Model Name: Dr. Lady
- Model No. MK-000A
- Manufacturer: Buheung Medical
- Origin: Republic of Korea
- Power: Li-Polymer Rechargeable Battery



Size & Weight

Main : 63mm x 88mm x 29mm / 76g

Probe : 115mm, 24Ø / 42g

Major Material : Titanium, Poly Carbonate (PC)



Company Introduction

BuheungMedical One of Well known Professional Medical & Health Equipment Manufacturer. BuheungMedical is Founded on the Christmas Eve of 1980, in Seoul, Korea, with motto (Present Health & Happiness to human being) and have developed with Medical & Health Equipment Since Foundation we have been concentrating on Personal Health care and Wellness Product, with its Professionalism, Strict standard of Quality, Variety of Unique technology (which protected by PATENTS). Since 1997, we are intensively studying and developing ELECTRO PULSE THERAPY as Interferential Current therapy, Electro Acupuncture, High-Frequency, Medium, Frequency and Low-Frequency. In these days, since 2006, we are broadening the range of Electro Pulse Therapy's category with other pulse therapy method as Temperature, Vibration, Ray, etc. Now BuheungMedical is one of the well-known-specialists in Medical Industry. Buheung Medical is the oldest established therapeutic apparatus company in Korea, and also Brands and Products are well known in many countries including Germany, Japan, France, and Korea. Now we are trading with over 30 countries in the world including Germany, Japan, USA, Switzerland, France, Italy, Russia, China, Singapore, Vietnam, Malaysia, Thailand, Indonesia, etc.

Awards


- [Korean World-Class Product Award] : The Ministry of Knowledge & Economy, Republic of Korea
- [Seal of Excellence] : KOTRA
- [Hi Seoul Best Brand] : Seoul City, Korea
- [The Prospective Enterprise of Korea-Japan Trade] : The Ministry of Knowledge & Economy, Republic of Korea
- 2019 Ministry of Trade, Industry and Energy Award
- 2019 World Patent Convention Award
- Korea Medical Equipment Industry Association President Award




Certification

ISO 13485 / 9001 / GMP / UL / FDA
KFDA / TUV / DNV / S-Mark

Patents

Over 50 Patents in the World







176, Butkkot-Road, Geumcheon-Gu, Seoul, Korea

Tel. 02)838-5866

Fax. 02)838-5334

E-mail. sales@buheung.co.kr

kh.moon@buheung.co.kr

homepage. www.buheungmedical.com