

Why FreeQ?

User Benefits Exclusive by <FreeQ>



Free from WIRE

FreeQ is free from wires attached to the electro-pads. It is a completely WIRELESS device that allows for you to attach it to any part of your body. Makes attending to troubled spots simple, effective, and time efficient.



Free from ELECTRO-PAD

No more sticky hassle of electro-pads on your body. It eliminates discomforts related to electro-pads like hair being pulled or skin irritation. FreeQ's innovative vacuum cups design replaces uncomfortable pads and makes it trouble free!



Free from Supplement

Frees you from having to purchase any supplementary items like the electropads.



Free from Sterilization (NANO SILVER material applied)

FreeQ implements the state-of-the-art Nano-applied technology with high effective maintenance for proliferation of bacteria, sterilization and bacterial odor.



Free from Discomfort

FreeQ maximizes the effect of therapy without experiencing the discomfort of conventional TENS devices of uncomfortable stimulation. FreeQ frequency (FreeFrequency: unique variation of middle frequency), resolves This is achieved through the deep penetration into the muscle fiber with higher frequency than normal TENS.



Free to do anytime & anywhere

FreeQ's design make things comfortable, convenient, and portable. The size allows you to take FreeQ anywhere and experience the benefit at anytime. It fits into pockets, handbags, brief cases, gym bags, and others, ready to be used hassle free that delivers the desired results.

BUHEUNG MEDICAL Co.,Ltd,Specializing in Personal Health Care and Wellness Products

Manufactured by:
Buheung Medical Co.,Ltd, Specializing in Personal Health Care and Wellness Products.
 With the mission statement of "To Present HEALTH & BEAUTY to Human Beings" the company strives to be the innovative leader in the industry. The company is protected by patent law supported by rigorous developments of highest professionalism in quality, product standard, and technology.
 Patented products of Buheung are currently exported to over 30 countries worldwide.
 The company is located in Seoul, Korea.
 Located in Seoul, Korea, have manufactured Medical & Health devices for 40 years.



C Type package B Type package H1 Type package H2 Type package H3 Type package

BUHEUNG MEDICAL Co., Ltd.
 www.buheungmedical.com
 176, Buttkot-ro, Geumcheon-Gu, Seoul, 153-814, Korea Tel: +82 2 838 5332
 Fax: +82 2 838 5334 E-mail: enyce8181@naver.com sales@buheung.co.kr

Wire Free! Ever First in the World!
 TENS / ICT / EMS / CUPPING AT ONCE

FreeQ!

Worldwide Patented



BUHEUNG MEDICAL Co., Ltd.
 www.buheungmedical.com ISO 13485 / ISO 9001 CE 0434

Wire Free! Ever First in the World!
 TENS / ICT / EMS / CUPPING AT ONCE

FreeQ!

Worldwide Patented

Nano-Silver Technology applied
 FreeQ's vacuum cup and conductive pad are added with antimicrobial colloid composed of silver-nanoparticles stabilized by polymers and exhibit excellent antimicrobial efficacy to a wide spectrum of microorganisms.



C Type

Multi-Functional Electrical Nerve Stimulator

- FUNCTION 1** ACUTE PAIN / CHRONIC PAIN CONTROL
- FUNCTION 2** BODY TONING % MUSCLE TRAINING
- FUNCTION 3** ELECTRO-ACUPUNCTURE
- FUNCTION 4** SLIMMING / BEAUTY & DIET
- FUNCTION 5** BUHANG (VACUUM CUPPING) THERAPY



B-6 Type



D Type



H1 Type



H2 Type




H3 Type

Ever First in the World! World Patented Product! FreeQ

FreeQ is worldwide patented item designed to promote multiple benefits through the comprehensive medical function design that includes: Middle/low frequency stimulation with ICT(Interferential Current Therapy), TENS(Transcutaneous Electrical Nerve Stimulation), EMS(Electric Muscle Stimulation), Electro Acupuncture and Vacuum Cupping Therapy. High effectiveness consists of application of middle frequency for deeper penetrating power into the body than the standard TENS devices. FreeQ allows people with minor acute or chronic pain to resume with daily life in relief and improve muscle condition and strength through the body toning and muscle training function.


Deep Penetration of FreeFrequency™

FreeFrequency consisting mostly of the middle frequency has more deep penetrating power into the body than TENS Therapy.



More Effective ICT(Interferential Current Therapy) Treatment

Using interferential current generated by two different kinds of middle frequency electric current, it treats discomforts related to muscle, skeleton and nervous system.




Mono polarity of electro-pad has only linear effect on the body surface.

A single vacuum cup works as a bipolar effect as two pads of TENS and with combination of adjacent vacuum cup. It has multiple (2ⁿ) effects on the body.

Buhang (Vacuum Cupping Therapy) Effects

A single vacuum cup works in a bipolar effect equaling to two pads of TENS. When used with combination of vacuum cups, it promotes multiple (2ⁿ) effects on the body.





Please note that the bruise mark can remain, depending on the body condition, longer than 2 hours, sometimes a couple of days.

FreeFrequency™ Therapy


Pain Control

FreeFrequency allows people with minor chronic pain to resume with daily life in relief. Commonly treated areas are:
 • Migraine • Acute Headaches • Chronic shoulder pain • Chest pain • Rheumatoid Arthritis • Chronic Neck and Back Pain • Chronic low back pain due to giving birth • Knee pain




Body Toning & Muscle Training

It is generally known that approximate 2500Hz of electric frequency helps with muscle training and rehabilitation. This application allows for deeper muscle penetration and stronger contraction of muscle fibers. FreeQ-Tone therapy is utilized by professional athletes, physical therapists, and sports doctors due to the excellent level of effectiveness. FreeQ-Tone's rhythmic pulses that contract and relax your muscles is similar to that of actual exercise which results to the same chemical change within the muscle.



Electro-Acupuncture

Application of 1-30Hz stimulation of FreeQ can have an acupuncture effect on the troubled parts of the body. Using lower voltage than TENS gives shaper but lighter stimulation.



Beauty & Diet

The chemical reaction resulted from muscle contraction mobilizes (burns) fat and glycogen in the muscle, which enables you to have better contour.

