

**DS MAREF**  
DAESUNG MAREF CO.,LTD.

# CTC-7

Cryo.Thermo.Compression system



# Specifications

- + Power source : AC 100-240V, 50/60Hz
- + Power consumption : 270VA
- + Mode : Cold, Heat, Contrast
- + Temperature : Cold : 6- 5°C ±3°C / 42.8-19°C ±37.4°F  
Heat : 35-45°C ±3°C / 95-113°C ±37.4°F
- + Time Range : Cold 5 - 60 min / Preset mode  
Heat 5 - 30 min / Preset mode
- + Dimension : 10.7(W) x 7(D) x 9.4(H)ft  
27.2(W) x 18(D) x 23.9(H)cm
- + Weight : 3.7kg / 8.2lb (Empty water container)  
3.9kg / 8.6lb (Full water container)
- + Type of protection IP : IPX0
- + Recommended Liquid : Distilled water
- + Class / type of protection : Class I/BF-type
- + SAFTY : IEC 60601-1
- + EMC : IEC 60601-1-2

# Iceless & Continuous Therapy

Pain Relief

Restraining Inflammation

Reducing Swelling



Early Discharging

Saving Cost



# The CTC-7 Wraps

are designed to maximize the functionality of the CTC-7 system by providing cold or heat water with air compression.

In particular, the compression pad covers the inner water pad so that cold/heat water can flow evenly from side to side over the wound, and it adheres more closely the water pad on the skin through periodic air compression to make cold/heat treatment more effective. In addition, automatic compression and interval time are applied intermittently.

Appropriate straps are attached to the wrap so that it can be adjusted and secured to fit a variety of sizes.

CTC wrap has its own drain valve, so it is easy to keep it dry after draining it when storing it after use.



# Main functions

Features

Iceless Operation

Continuous Use Time

Superior Wraps Design

Protocols

Cold + Compression

Heat + Compression

Contrast + Compression

Convenience

Auto Water Check & Alarm

Steam Cleaning

Auto Drain

# One Person Use Wraps



CTCAFD0



CTCUMD0·CTCULD0



CTCKFD0



CTCEF00



CTCSMD0



CTCFD0

Order No.	One Person Use Wraps List	Size
CTCKFD0	KNEE WRAP	UNIVERSAL
CTCUMD0	UNIVERSAL WRAP	M
CTCULD0	UNIVERSAL WRAP	L
CTCBFD0	BACK WRAP	UNIVERSAL
CTCEF00	ELBOW WRAP	UNIVERSAL
CTCSMD0	SHOULDER WRAP	M
CTCAFD0	ANKLE WRAP	UNIVERSAL

# Resuable Wraps



CTCAFD0



CTCUM00·CTCUM00



CTCKFD0



CTCEF00



CTCSM00



CTCFD00

Order No.	Resuable Wraps List	Size
CTCKFD0	KNEE WRAP	UNIVERSAL
CTCUM00	UNIVERSAL WRAP	M
CTCUL00	UNIVERSAL WRAP	L
CTCBFD0	BACK WRAP	UNIVERSAL
CTCEF00	ELBOW WRAP	UNIVERSAL
CTCSM00	SHOULDER WRAP	M
CTCAFD0	ANKLE WRAP	UNIVERSAL