

# EZ-BIOPSY

EZ-BIOPSY Semi-Auto / soft tissue biopsy needles



## Individual components

- Two notch sizes - 10 or 20mm
- pliable finger grip (Rubber)
- The compact, ergonomically-contoured body
- Four different grip color codes(14G : Green, 16G : White, 18G : Orange, 20G : Yellow)
- Easy to loading

## Size

Gauge	Length 7cm	Length 10cm	Length 13cm	Length 15cm	Length 20cm
14	HWA14-07	HWA14-10	HWA14-13	HWA14-15	HWA14-20
16	HWA16-07	HWA16-10	HWA16-13	HWA16-15	HWA16-20
18	HWA18-07	HWA18-10	HWA18-13	HWA18-15	HWA18-20
20	HWA20-07	HWA20-10	HWA20-13	HWA20-15	HWA20-20



## Soft tissue biopsy needle

- EZ-BIOPSY
- EZ- AUTO BIOPSY

CE  
2195




# EZ- AUTO BIOPSY

EZ-BIOPSY AUTO / soft tissue biopsy needles

Fixation Device for Hernia Repair

# EZ-COAXIAL

Coaxial Needle / Universal coaxial introducer needles



## Individual components

No	Features	
1	Upper case part Needle	
2	Stylet	
3	Notch	
4	1cm, 2cm Lever	1cm, 2cm adjustable
5	Loading Bar	Easy loading
6	Lock, A, S	Lock, A : Fully Automatic, S : Semi Automatic
7	Push Button 1	Shooting Button 1
8	Push Button 2	Shooting Button 2

## Size

Gauge	Length 7cm	Length 10cm	Length 13cm	Length 15cm	Length 20cm
14	HWA14-07	HWA14-10	HWA14-13	HWA14-15	HWA14-20
16	HWA16-07	HWA16-10	HWA16-13	HWA16-15	HWA16-20
18	HWA18-07	HWA18-10	HWA18-13	HWA18-15	HWA18-20
20	HWA20-07	HWA20-10	HWA20-13	HWA20-15	HWA20-20

## Individual components

- Clear luer lock hub
- Centimeter markings
- Light weight and extra thin tubing

## Size

Gauge	Length 7cm	Length 10cm	Length 13cm	Length 15cm	Length 20cm
14	HWA14-07	HWA14-10	HWA14-13	HWA14-15	HWA14-20
16	HWA16-07	HWA16-10	HWA16-13	HWA16-15	HWA16-20
18	HWA18-07	HWA18-10	HWA18-13	HWA18-15	HWA18-20
20	HWA20-07	HWA20-10	HWA20-13	HWA20-15	HWA20-20