


# Irrigation Pump

- HW-1000
- HW-2000



# Irrigation Pump

HW-1000 / HW-2000



HW-1000



HW-2000

## Individual components

No	Features
1	Control Switch      Pressurize / Depressurize
2	Rapid Deflation Switch      On / Off
3	Alarm
4	Setting Pressure Display
5	Real Temperature Display
6	Setting Temperature Display
7	Pressure Regulator
8	Temperature Regulator
9	Latch
10	Cover
11	Air Source Connector
12	Bag Hanger
13	Power      On / Off

## Specification

Classification	
Power	AC 220V, 60Hz
Pressure Range	100mmHg ~ 400mmHg
Pressure Increasing	5mmHg
Temp. Range	0℃ ~ 50℃
Temp. Increasing	1℃
Alarm	Over Temperature (Over 20mmHg than Setting pressure) Under Temperature (Under 20mmHg than Setting pressure) Over Temperature (Over 5℃ than Setting temperature) Under Temperature (Under 5℃ than Setting temperature)
Weight	2kg (Single), 5kg (Dual)
Dimensions	170(w)x142.05(L)x415.8(H)