


# Ligation Automatic Clip Applier System for Laparoscopic

LigaStar Polymer Auto Endo Clip



  
**Ligation Automatic Clip**  
Applier System for Laparoscopic  




# Ligation Automatic Clip


## Applier System for Laparoscopic

### LigaStar Polymer Auto Clip 5mm(M/L), 10mm(L)



### LigaStar Polymer Auto Clip Features

- Pneumatic firing system clip included in the application jaw for improved stability
- Available in the clip sizes (M/L, L)
- Pneumatic reusable automatic Handle
- 6Clips, 9Clips, 12Clips, 15Clips, 18Clips single use clip Shaft (5mm, 10mm)
- Minimize clip dropout rate by applying jaw drop-out prevention structure
- Square pyramid shape inside the clip minimizes slippage both horizontally and vertically
- 35cm Reusable shaft, 5mm
- A Square pyramid-shaped inner clip profile provides an enlarged contact surface



Square pyramid shape

Jaws for enhanced stability



#### ► Order information



### LigaStar Polymer Auto Clip 5mm(M/L), 10mm(L)

Product code	Art No	Description	Packaging
Pneumatic reusable Handle	LGS- H	Reusable Handle	1pcs/box
6 Clips Shaft	LGS-PS-06	5mm, 10mm, M/L, L, 6clips, CO <sub>2</sub> Cylinder, 33cm	10pcs/box, 6boxes/carton
9 Clips Shaft	LGS-PS-09	5mm, 10mm, M/L, L, 9clips, CO <sub>2</sub> Cylinder, 33cm	10pcs/box, 6boxes/carton
12 Clips Shaft	LGS-PS-12	5mm, 10mm, M/L, L, 12clips, CO <sub>2</sub> Cylinder, 33cm	10pcs/box, 6boxes/carton
15 Clips Shaft	LGS-PS-15	5mm, 10mm, M/L, L, 15clips, CO <sub>2</sub> Cylinder, 33cm	10pcs/box, 6boxes/carton
18 Clips Shaft	LGS-PS-18	5mm, 10mm, M/L, L, 18clips, CO <sub>2</sub> Cylinder, 33cm	10pcs/box, 6boxes/carton

# Ligation Automatic Clip

## Applier System for Laparoscopic

LigaStar Titanium Auto Endo Clip



  
**Ligation Automatic Clip**  
Applier System for Laparoscopic





# Ligation Automatic Clip Applier System for Laparoscopic

LigaStar Titanium Auto Clip 5mm(S/M), 10mm(M/L)




► Order information




Product code	Art No	Description	Packaging
Pneumatic reusable Handle	LGS- H	Reusable Handle	1pcs/box
6 Clips Shaft	LGS-TS-06	5mm, 10mm, S/M, M/L, 6clips, CO <sub>2</sub> Cylinder, 33cm	10pcs/box, 30boxes/carton
9 Clips Shaft	LGS-TS-09	5mm, 10mm, S/M, M/L, 9clips, CO <sub>2</sub> Cylinder, 33cm	10pcs/box, 30boxes/carton
12 Clips Shaft	LGS-TS-12	5mm, 10mm, S/M, M/L, 12clips, CO <sub>2</sub> Cylinder, 33cm	10pcs/box, 30boxes/carton

## LigaStar Titanium Auto Clip Features

- Pneumatic firing system clip included in the application jaw for improved stability
- Available in the clip sizes (S/M, M/L)
- Pneumatic reusable automatic Handle
- 6Clips, 9Clips, 12Clips single use clip Cartridge (5mm, 10mm)
- Minimize clip dropout rate by applying jaw drop-out prevention structure
- Comb-shape inside the clip minimizes slippage both horizontally and vertically
- 35cm Reusable shaft, 5mm, 10mm
- Comb-shaped inner clip profile provides an enlarged contact surface



Reusable shaft



Disposable Cartridge  
6Clips, 9Clips, 12Clips contained in a Cartridge

