

### System components



Hand Controlled - Pistol Grip Handle



Hand Controlled - Pencil Grip Handle

#### 34cm Electrodes



Hook Electrode  
Hollow Tip



Curved Dissector Electrode  
Hollow Tip



Spatula Electrode  
Hollow Tip



Only Suction Irrigation  
Hollow Tip



Right Angle Electrode  
Hollow Tip

### Order information

Product Code	Description	Shaft Diameter (mm)	Shaft Length (cm)
<b>Handles</b>			
HK2-P	Hand Controlled, Pencil Grip Handle		
SH-P	Hand Controlled, pistol Grip Handle		
<b>Electrodes</b>			
SIE 01	Hook Electrode, Hollow Tip	5	34
SIE 02	Spatula Electrode, Hollow Tip	5	34
SIE 03	Right Angle Electrode, Hollow Tip	5	34
SIE 04	Curved Dissector Electrode, Hollow Tip	5	34
SIE 05	Only suction irrigation Tip	5	34

## Suction & Irrigation Electrode System

- HK2 S/I Electrode(Pencil grip handle)
- SH S/I Electrode(Pistol grip handle)

CE  
2195



MEDICAL  
EQUIPMENT  
MANUFACTURING



# Suction & Irrigation Electrode System / Disposable


Multi - Functional Electrosurgery System facilitates tissue dissection, coagulation, irrigation and fluid evacuation.

## HK2

Pencil Grip Handle


## SH

Pistol Grip Handle



SH – Pistol Grip Handle

HK2 – Pencil Grip Handle



### Features

- Performs dissection, coagulation, irrigation, and fluid evacuation
- Easy to activate suction and irrigation buttons
- Two handles and five electrodes /only S/I shaft
- 5mm shaft diameter
- Facilitate one-handed use
- Easy to push and pull shaft

### Individual components

No	Features	Benefits
1	Cut	Electrosurgery hand control
2	Coa	
3	Irrigation	Irrigation control , Easy to activate irrigation button
4	Suction	Suction control , Easy to activate suction button
5	Cable	ESU cable
6	Tube	Fluid line fittings
7	Probe button	Smooth movement button of the shaft probes prevents fatigue
8	Housing	Shaft connector
9	Shaft	Shaft
10	Tip	Electrode 4 options ( Hook, Curved Dissector, Spatula, Right Angle)

### Individual components

No	Features	Benefits
1	Coa	Electrosurgery hand control
2	Cut	
3	Irrigation	Irrigation control , Easy to activate irrigation button
4	Suction	Suction control , Easy to activate suction button
5	Cable	ESU cable
6	Tube	Fluid line fittings
7	Probe button	Smooth movement button of the shaft probes prevents fatigue
8	Housing	Shaft connector
9	Shaft	Shaft
10	Tip	Electrode 4 options ( Hook, Curved Dissector, Spatula, Right Angle)
11	Connector	Shaft connector botton