


Vasoclip[®]

Safe and Reliable Ligation System




Locking polymer clip


Usability



Vasoclip dedicated Applier

- The invention of special jaw calibration instrumentation technology.
- Prevent Slip-Down when loading clips.
 - ☞ Stable and safe surgery available
- High compatibility
 - ☞ Optimize / Maximize clip utilization






Stability



Protrusion Hinge

- Increased resistance during clip engagement due to contact deformation of hinge protrusion.
- Equalizes the distribution of force in clip ligation by inducing directionality.



☞ Strengthen tensile to prevent blood leaks and clips from leaving

Safety



Anchor-shaped Teeth

- Increase lateral holding force.
- Prevention of separation and maintenance of fixing power.
- Prevention slipping and dislodging from the vessel by compressing the vessel.

☞ Minimize damage to blood vessels